

NovaPACS

Diagnostic Viewer User Manual

May 2025 | Version 9.1.6 | Part No. L-PL-IUP.022 | Rev. 1.0

U.S. Pat. No. 9978027 | U.S. Pat. No. 9741166 | U.S. Pat. No. 10229500

© Copyright 2025 Novarad Corporation. All rights reserved.

This publication may not be copied, reproduced, or converted to any other electronic or machine-readable format without prior written consent of Novarad[®] Corporation. NovaPACS, Nova RIS[®], and PACS Enterprise Imaging[®] are registered trademarks of Novarad Corporation.

Windows 11 is a registered trademark of Microsoft[®] Corporation. All other products or name brands are the trademarks of their respective owners.

Novarad Corporation USA

3152 North University Avenue | Suite 200 | Provo, Utah 84604

Support: 801-221-5895 | Fax: 801-642-1030 | Toll Free: 877-668-2723

Intended Use: 

NovaPACS is intended for the viewing, archiving, analysis, annotation, registration, distribution, editing, fusion, and processing of digital medical images and data acquired from diagnostic imaging devices and all DICOM devices, including mammography.

NovaPACS is intended for use by trained healthcare professionals, including radiologists, physicians, technologists, clinicians, and nurses. NovaPACS allows the end user to display, manipulate, archive, and evaluate images.

Mobile devices are not intended to replace a full workstation and should be used only when there is no access to a workstation. They are not to be used for mammography or fMRI. Mobile devices are used for diagnosis of medical images from different modalities including CT, MR, US, CR/DX, NM, PT, and XA. For a list of compatible mobile platforms, see NovaPACS Diagnostic Viewer User Manual.

While NovaPACS's full workstation provides tools to assist the healthcare professional determine diagnostic viability, it is the user's responsibility to ensure quality, display contrast, ambient light conditions, and to confirm image compression ratios are consistent with the generally accepted standards of the clinical application.

NovaPACS is intended for providing analysis and visualization of functional MRI data of the human brain, presenting derived properties and parameters in a clinically useful context.

Warning:

Digital mammography images must be acquired from an FDA-cleared or approved Full Field Digital Mammography (FFDM) modality for primary interpretation and diagnosis.

Serious Incidents:

Any serious incident that has occurred in relation to NovaPACS medical device should be reported to Novarad and the competent authority of the Member State in which the user and/or patient is established.

Safety Assurance:

Nova RIS is safe for our consumers and the environment. Nova RIS complies with the General Product Safety Directive (2001/95/EC). This holds Novarad responsible to provide our consumers with quality products that are safe for use.

Installation Instructions:

Installation is the responsibility of Novarad Technical Services. The user is directed by Novarad Technical Services if their assistance is required.

NovaPACS is only intended for use with hardware that meets or exceeds minimum performance requirements, satisfies current safety regulations, and is CE marked as required by national or regional regulations.

Hardware Specifications for NovaPACS

Below are the minimum server and workstation requirements for using the NovaPACS Diagnostic Viewer. To use the NovaPACS Diagnostic Viewer systems at their highest performance, it is best to meet or exceed these recommendations. These system recommendations are subject to change.

a) Minimum Server Specifications:

REQUIRED PERFORMANCE BENCHMARKS (The following minimum benchmarks must be reached using the PassMark Performance Test tool. A free version can be downloaded from <http://www.passmark.com/products/pt.htm>)

CPU Mark: 11,000

Memory Mark: 2,100

Disk Mark: 28,000

Passmark Rating: 2,500

The hardware used to meet the above benchmarks is based on the following Specifications:

Processor: Intel Xeon Silver 4110 2.1G, 8C/16T, 9.6GT/s , 11M Cache, Turbo, HT (85W) DDR4- 2400, Bus Speed 8.00 GT/s DM13, TDP 80W

Operating System: Windows Server 2022

Memory: 32 GB DDR4 DR x4 ECC x2 2666MT (In pairs for Dual Channel configuration)

Drive Controller: 2GB NV Cache

Hard Drives: 1TB M.2 SSDs in Dell Boss Card RAID-1; partitioned for OS (C:), Processing drive (D:) 3x4TB 7200 RPM SATA (RAID-5) for Archive*

CD/DVD: None

Network Card: 1 GB Ethernet NICs

Power Supply: Dual, Hot-plug, and Redundant Power Supply (1+1), 750W

Storage: 200Mbps sustained read/write throughput to the disk systems, at minimum

*Drive size used in production may vary based on archive

b) Minimum Workstation Specifications

Standard (Tech, Clinical, Web Viewer, Telerad, Gateway) Workstation—

Operating System: Windows 11 (latest version)*

Processor: Intel® Core™ i3-12300

Hard Drive: 1 TB M.2 SSD

Ethernet: 1 Gigabit

Graphics Card: on-board video, DP

Memory: 16 GB RAM

Screen Resolution: 1600x1200 pixels or higher

Optical Drive: 8x DVD+/-RW Optical Disk Drive

Enterprise Workstation (no Tomo/Mammo)—

Operating System: Windows 11 (latest version)

Processor: Intel® Xeon W-2223

Hard Drive: 1 TB M.2 SSD

Ethernet: 1 Gigabit

Graphics Card: NVIDIA RTX-A4000, 16 GB, 4 DP

Memory: 32 GB RAM

Screen Resolution: 3 megapixels or higher

Optical Drive: 8x DVD+/-RW Optical Disk Drive

Tomo/Mammo / 3D Workstation / High-End Tech Review Workstation—

Operating System: Windows 11 (latest version)

Processor: Intel® Xeon W-2223

Hard Drive: 1 TB M.2 SSD

Ethernet: 1 Gigabit

Graphics Card: NVIDIA RTX-A4000, 16 GB, 4 DP

Memory: 32 GB RAM

Screen Resolution: 5 megapixels or higher

Optical Drive: 8x DVD+/-RW Optical Disk Drive

c) Accessory Requirements for Use

- *A modern supported Microsoft Operating System*
- *A modern supported internet browser*
- *Server hardware that meets minimum specifications*
- *Workstation hardware that meets minimum specifications*
- *Current Virus Protection*
- *A tablet with a modern supported version of iOS Operating System for the iOS Tablet Viewer*
- *At least a five-button mouse*
- *A supported dictation microphone*
- *A modern USB keyboard*
- *At least a 3MP Monitor for Radiology and a 5MP monitor for Mammography*

Maintenance/Upgrade Procedures:

Maintenance, upgrades, and warranty replacements are the responsibility of Novarad Technical Services and the Original Equipment Manufacturer (OEM). Novarad Technical Services will assist the user if needed.

Hardware Disposal:

Contact Technical Support regarding any hardware warranty replacement or disposal issues before taking action. Disposal of hardware should be authorized by technical support and done in accordance with the procedures specified in the original equipment manufacturer's manual.

All patient data is the responsibility of the acquiring site. When the NovaPACS system is decommissioned, the control, management, and exporting of all imaging data is 100% the responsibility of the customer. Novarad provides access to the NovaPACS software until the migration of data out of NovaPACS and into the new PACS system is complete; a table showing all the Study UIDs can also be provided as a reference, if desired. No other Novarad assistance or involvement occurs. When the migration is complete, any hardware owned by Novarad is required to be wiped by the customer or Novarad Support and returned to Novarad headquarters. In the case where NovaPACS was installed and images are stored on customer owned VMs, the Novarad software is simply decommissioned by Novarad Support after the images are migrated by the customer, with no other Novarad involvement.

Degree of Accuracy for measuring functions:

Measurements done on DICOM images that include pixel-scaling data are precise within two-pixel lengths. NovaPACS does not return a measurement with more precision than is appropriate for the DICOM data and resolution of the images.




Due to the imprecise nature of a computer monitor, a mouse, and a person's hand, calibrating measurements by hand can result in inaccurate measurements. Measurements taken after calibration should be treated as rough estimates. Calibration is only valid for the current image unless the Apply to Series option is selected during calibration.




Environment:

NovaPACS's intended operational environment is located in radiology departments in hospitals and clinics as well as imaging centers. Mobile NovaPACS may be used in any environment where mobile devices may be taken.

LIST OF SYMBOLS

The following symbols appear in this manual:

Symbol	Meaning
	Indicates important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Denotes the manufacturer.
	Denotes the manufacturer's European Community representative.

	<i>Medical Device</i>
	<i>Denotes the Importer.</i>
	<i>Denotes the Distributor.</i>



Manufacturing:

Novarad Corporation
3152 North University Ave.,
Suite 200
Provo, Utah 84604
USA
Tel.: +1.801.642.1001
2012



Authorized Representative:



MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM The Hague
The Netherlands



UK Responsible Person:

MIS Healthcare
Unit 26-27
Capitol Way
Capitol Park
London NW9 0EQ
United Kingdom
Tel: +44 208 205 9500



Philippines Distributor/Importer:



Enterprise Imaging Technologies Inc.
Unit 1208 12/F Centuria Medical Makati
Kalayaan Ave., corner Salamanca St.
Makati City, Metro Manila
Philippines
Registration #: 2021040010266-06
Tel: +632 85304 7921

Australian Sponsor:

Emergo Australia
Level 20 Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000 Australia



Europe Distributor:

Information Storage S.L.
Enric Valor, 3 – Edificio ARAV
46100 Burjassot
Valencia, Spain

Contact Information:

Technical Support & Regional
Distributors
www.novarad.net/contact-us
info@novarad.net
support@novarad.net

LIST OF SYMBOLS.....	iv	5-16—View Study on Other Workstation	50
CHAPTER 1—GETTING STARTED	1	5-17—Launch Orthopedic Tools.....	50
1-1—Introducing the Study Browser and Image Viewer	1	5-18—Teaching File.....	50
1-2—NovaPACS Diagnostic Viewer System	2	5-19—Study of Interest	53
1-3—Accessing the Diagnostic Viewer	2	5-20—Collaboration	54
1-4—Logging into the Diagnostic Viewer	4	CHAPTER 6—USING THE STUDY BROWSER VIEW MENU	58
CHAPTER 2—USING THE STUDY BROWSER.....	7	6-1—Worklist Tab	58
2-1—Filtered Studies List.....	7	6-2—Manage Hanging Protocols.....	58
2-2—Patient History List	12	6-3—Patient History	59
2-3—Series List.....	14	6-4—Series 59	
2-4—Preliminary Impression Lists.....	15	6-5—Refresh Lists.....	59
2-5—Images Tab	18	6-6—Reset Filters	60
2-6—Report Tab	19	6-7—Docking Tabs.....	61
2-7—Study Browser Tab Functionality	19	6-8—Radiologist Productivity Display.....	61
2-8—Hanging Protocol Layouts.....	20	CHAPTER 7—USING THE STUDY BROWSER SETTINGS MENU	63
2-9—The Search Worklist.....	22	7-1—Options	63
2-10—The Assigned Worklist	28	7-2—W/L Presets	80
2-11—Advanced Worklists	30	7-3—Manage Hanging Protocols.....	81
2-12—Customizing Study Browser Lists.....	31	7-4—Configure Monitors.....	87
2-13—Using the Study Browser Menus/Custom Toolbar.....	32	7-5—Preferences.....	88
2-14—Activating the Study Browser	33	7-6—Customize Toolbar	91
2-15—Customizing the Study Browser Color and Font Size.....	34	CHAPTER 8—USING THE STUDY BROWSER HELP MENU....	94
CHAPTER 3—USING THE STUDY BROWSER FILE MENU	36	8-1—Instructions.....	94
3-1—Open from Directory	36	8-2—Current User Information	94
3-2—Open from File	36	8-3—Contact Support.....	94
3-3—Exit the Diagnostic Viewer	36	8-4—Submit Feedback	95
CHAPTER 4—USING THE WORKLIST MENU	37	8-5—User Manual	95
4-1—Introduction to the Assigned Worklist	37	8-6—What’s New	95
4-2—Add to Worklist	37	8-7—About 95	
4-3—Open Worklist Item	38	8-8 Feature Highlight Dialog.....	96
4-4—Next Worklist Item.....	39	CHAPTER 9—USING THE IMAGE VIEWER	97
4-5—Remove Worklist Item.....	39	9-1—Introduction to the Image Viewer.....	97
4-6—Clear Worklist.....	39	9-2—Hanging Protocol Drop-Down Menu	98
CHAPTER 5—USING THE STUDY BROWSER STUDY MENU..	40	9-3—Customizing Image Viewer Toolbars	100
5-1—Select All	40	9-4—Moving Toolbars	101
5-2—Open 40		9-5—Locking Toolbars.....	101
5-3—Close Workspace	41	9-6—Accessing Menu Options.....	101
5-4—Close All Workspaces.....	41	9-7—Image Viewer Icon Buttons	101
5-5—Route 41		9-8—Dragging and Dropping Images	102
5-6—Documents	42	9-9—Copy to Clipboard.....	104
5-7—Exam Report.....	42	CHAPTER 10—USING THE IMAGE VIEWER STUDY MENU	106
5-8—Manage Report Workflow	43	10-1—Patient History.....	106
5-9—RIS Information	43	10-2—Select All Images in a Study	107
5-10—Match Study to Order	44	10-3—Close Study	108
5-11—Retrieve Images	45	10-4—Close Study and Show Study Browser	108
5-12—Patient CD/DVD/USB	46	10-5—Close Workspace	108
5-13—Assign Password	49	10-6—Close All Workspaces	108
5-14—Scan Complete	49	10-7—Close All Workspaces and Show Study Browser	108
5-15—Launch Mammo Viewer.....	49		

10-8—Mark Reviewed, Close Workspace, and Show Study Browser	108
10-9—Manage Hanging Protocols	109
10-10—Add Split	114
10-11—Choose Series	114
10-12—Choose Series Box Docking	115
10-13—Auto Sync	116
10-14—Cross Localize	116
10-15—Cross Localize	117
10-16—Set Sync	117
10-17—Aggregate	117
10-18—Print	118
10-19—Route	119
10-20—Documents	120
10-21—Exam Report	120
10-22—Manage Reports Workflow	120
10-23—RIS Information	121
10-24—Dictate	121
10-25—Launch Mammo Viewer	124
10-26—View Study on Other Workstation	124
10-27—Launch Orthopedic Tools	125
10-28—Assign Password	125
10-29—Notes	125
10-30—Collaboration	125
10-31—Teaching File	127
10-32—Research	129
10-33—Upload to GlobalRad	129
10-34—Save View Settings	129
10-35—Show All Series	129
10-36—Series Layout	129
10-37—Hanging Protocol Layouts	129
10-39—Minimize	130
10-40—Exit	130

CHAPTER 11—USING THE IMAGE VIEWER SERIES MENU 131

11-1—Select All Images	131
11-2—Close	131
11-3—Views	131
11-4—Split Horizontal	134
11-5—Split Vertical	134
11-6—Choose Image Layout	134
11-7—Calcium Scoring	134
11-8—Sort Images	135
11-9—Colorize	135
11-10—Create New View	136
11-11—Subtract	136
11-12—Advanced Analysis	137
11-13—Remove Images	138
11-14—Ejection Fraction	138
11-15—Percent Stenosis	139

CHAPTER 12—USING THE IMAGE VIEWER IMAGE MENU 141

12-1—Zoom	141
-----------	-----

12-2—Zoom Map	141
12-3—Reload	142
12-4—Invert	142
12-5—Subtract Clip Frames	142
12-6—Calibrate Measurement	142
12-7—Details	143
12-8—Flip/Rotate	144
12-9—Annotate	144
Accessing the Annotate Options	144
12-10—Navigate Through Annotated Images	151
12-11—Move Images	151
12-12—Export to File	152
12-13—Magnification	152
12-14—Undo All Selections	152
12-15—Key Image Series	152
12-16—Critical Findings	153
12-17—Equalize Histogram	154
12-18—Set Win/Level	155
12-19—Smooth	155
12-20—Sharpen	155
12-21—Local Contrast Enhancement	155

CHAPTER 13—MISCELLANEOUS INFORMATION 155

13-1—Navigation	156
13-2—The Radius Menu	156
13-3—Additional Software Products	159

Index 165

Document revision history 169

CHAPTER 1—GETTING STARTED


NovaPACS Diagnostic Viewer is a world-class product and is one of the original Picture Archiving and Communication Systems (PACS). The NovaPACS Diagnostic Viewer includes many innovations created by Novarad. NovaPACS is known in the industry as a product that leads to radiology efficiency through a mouse-driven viewer, customizable interface, and easy access to advanced features. NovaPACS Diagnostic Viewer also includes web-based functionality, dual monitors, viewers, technology tools, an administrative console, and on-site seven-year archives.

Novarad understands that workflow and efficiency are highly valued, and all products are built around these concepts. NovaPACS provides radiology departments with the following:

- The best web-based technology utilizes high-level streaming, imaging, algorithms, and proprietary methods for distributing large data sets quickly.
- Images and reports that are available anytime, anywhere, all accessible via the latest web browser version.
- A full Diagnostic Viewer that is accessible from any Microsoft Windows computer and location, and is feature-rich, including advanced hanging protocols, MIP/MPR, and voice annotation.
- A user-friendly, easy-to-learn Diagnostic Viewer that has customizable buttons, toolbars, menus, hotkeys, and mouse settings. These settings are saved as a user-specific profile that provides continuity when changing workstations.
- An Image Viewer that is mouse-driven and allows the radiologist to focus on the images on the screen rather than on keys or buttons.
- A Study Browser that displays the entire patient's image history, the number of series for a study, and the number of images in a series, along with many other features.

1-1—INTRODUCING THE STUDY BROWSER AND IMAGE VIEWER

The **Study Browser** is essentially the control panel for the Diagnostic Viewer. The Study Browser displays information in branches or lists and shows details such as the entire patient's image history, the number of series for a study, and the number of images in a series. In addition to these details, the viewable columns can be set to show studies with existing reports and the name of the reviewing radiologist.



The **Image Viewer** displays the selected images from the Study Browser. Radiologists can design their workspace with the Image Viewer's flexible display options, allowing them to move efficiently through their caseload.



1-2—NOVAPACS DIAGNOSTIC VIEWER SYSTEM

⚠ **Note:** NovaPACS is not intended for use on smartphones. Any unsupported mobile platforms display a persistent on-screen message of “Not intended for diagnostic use.”

The following light level conditions are supported for diagnostic reading:

- A typical radiology reading room (7.3 lux)
- An office (250 lux)
- A bright room (850 lux)

⚠ **Note:** Lighting environments outside the range of 7.3 to 850 lux are not recommended for diagnostic reading.

⚠ **Note:** The brightness setting on the mobile device should always be at 100% for maximum detail. Adjustment outside of the maximum brightness is not recommended.

1-3—ACCESSING THE DIAGNOSTIC VIEWER

The Diagnostic Viewer functions essentially the same, no matter how it is accessed. There are two ways users can access the **Diagnostic Viewer**:


1. From a computer or workstation—referred to as accessing the **Managed Viewer**.
2. From a web browser—referred to as accessing the **Novaweb Launch** page. This page is where users can launch all Novarad products, including NovaPACS.

Note: Users are unable to route studies to a machine using the Web Viewer and studies are only cached temporarily if studies are sent to the machine. When using the Web Viewer from the Novaweb page, users must write access to the machine in order to retrieve images, share patient studies, and to save user preferences locally. All other functionality is identical to using a Managed Viewer.

THE MANAGED VIEWER

The Managed Viewer is often faster and allows users to route images directly to the workstation for faster opening and display. The Managed Viewer gives users full access to all product functions. Many facilities have a combination of both Managed Viewers and Web Viewers. Novarad suggests radiologists and cardiologists use a Managed Viewer for best performance.

LAUNCHING THE MANAGED VIEWER



Novarad's **Image Viewer** icon shows on the computer's desktop after the Managed Viewer has been installed. Double-click the icon to launch the Diagnostic Viewer.


The login dialog box appears.

Note: Administrators can configure the Diagnostic Viewer to use a Windows login as the login for the Viewer. If using this configuration, users will not see a login prompt. Also, if configured, users can select the **Remember Me** checkbox to be automatically logged in for each subsequent visit. It is recommended that only those with a dedicated machine use this option.



THE WEB VIEWER

The Web Viewer is Novarad's thin-client version that downloads Image Viewer files temporarily and is accessed through a secure server using the Internet. Users must download and use Microsoft's .NET Framework 4 Client Profile to use the Web Viewer from the Novaweb page. This Viewer gives users full access to all product functions and can be accessed anytime from any computer through a website and a secure server.



- 1. Diagnostic Viewer**—Novarad's Diagnostic Viewer accessed via the web provides access to DICOM PACS with full product functionality. Radiologists and Radiology departments can use this version to work from an office, clinic, and/or hospital.
- 2. HTML5 Viewer**—The HTML5 Viewer is a web-based viewer.
- 3. Single Study Access**—Users can assign a single-study password in the NovaPACS Diagnostic Viewer.

LAUNCHING THE WEB VIEWER FROM THE NOVAREB PAGE

To launch the Web Viewer from the Novareb page, open the web browser and enter the site's PACS server address in the address bar. Press **Enter**.

Click the **Diagnostic Viewer** link and enter the credentials on the login page to access the Diagnostic Viewer.

- **MPR/3D Plug-in**—System administrators can download this one-time plug-in when using the **Web Viewer**.
- **Cardio Reporting Plug-in**—Users should download and run this one-time plug-in before using the Web Viewer with EncaptureMD.



Important: The Novarad Web Viewers now require multi-factor authentication. When users log in, they will be prompted to register for multi-factor authentication. Users will be sent a one-time passcode and instructions on how to use the MobileRad App to complete multi-factor authentication.



Warning: Novarad's Web Viewers are not intended for diagnostic purposes.

LOGGING INTO THE DIAGNOSTIC VIEWER

Click the **Diagnostic Viewer** option. Enter the username and password and click **Login**.

LOGGING IN AUTOMATICALLY

Users can select the **Remember Me** checkbox in the login dialog box. This saves the username and password for future logins and the Diagnostic Viewer launches automatically.

Note: Administrators can enable/disable the ability to automatically log in on the Machines tab of the NovaPACS Admin Console. If the option to automatically login is enabled, the Remember Me checkbox is available on the login dialog box. If disabled, no Remember Me checkbox displays. Use the Remember Me checkbox with caution, as any user accessing the NovaPACS Diagnostic Viewer is automatically logged in as this user until deselected on the Data Sources page. It is recommended that only those with a dedicated machine use this option.



Note: Administrators should contact Novarad Support to have their users' passwords auto expire after a set number of days.

LOGGING INTO A DATA SOURCE


If the site uses other servers to access images, users receive a prompt to log into the other data sources when the Diagnostic Viewer tries to connect to it. Some sites need to log in multiple times—which logs them into their multiple data sources.

Note: If users click **Cancel** for any data source login prompt, they can access and log into that data source at any time from within the Diagnostic Viewer. Clicking **Cancel** returns users to the Novareb page.

This prompt also appears when selecting an available data source within the Diagnostic Viewer that has not been authenticated.



Note: Contact a system administrator to log in across all data sources using the same username and password for each additional server.



Users can also log into data sources from inside the NovaPACS Diagnostic Viewer. If the user has not logged into an available data source at the time of login, they can click the name of the data source listed at the bottom of the Study Browser window to access the login dialog box for that data source. This enables users to quickly log into data sources at any time from the NovaPACS Diagnostic Viewer.


MULTI-FACTOR AUTHENTICATION

NovaPACS includes multi-factor authentication to help safeguard patient information. NovaPACS multi-factor authentication is integrated with MobileRad. Upon client authentication, users will see a dialog informing them that they must enroll in multi-factor authentication. The dialog provides instructions and a QR code that will link to the MobileRad app. The QR code contains the server name and username which can be used to log in to MobileRad. The server name is displayed above the **Continue** and **Skip for Now** buttons. Clicking **Continue** will initiate two-factor authentication with MobileRad.

After a successful client authentication, users will see the “A multi-factor authentication approval request has been sent” message. If users did not receive a request, click the **Didn't receive a request?** link to obtain a one-time passcode. The one-time passcode prompt will provide instructions and a QR code that will link to MobileRad.

NovaPACS multi-factor authentication is also integrated with Duo Security. When users login, if they are enrolled for multi-factor authentication, they will be sent an approval request to their registered device. Once the request has been approved, users will then have to scan a QR code to receive a one-time passcode to enter.

If users are having trouble with the multi-factor authentication, there is a **Still need help?** link located under the passcode field. This link will walk users through how to find their one-time passcode. Once users have correctly entered their one-time passcode, they will be logged in.




If users select the **Remember Me** button on the login screen, the system will remember the device so users will not have to enter two-factor authentication every time they log in.


Users can see who has been exempted from MFA through the Admin Console's Users Tab. Users will see an MFA exempt column. If a user has an 'X' by their name, they are exempt.

DISABLED DUPLICATE DATA SOURCES



Users can be prevented from logging into duplicate data sources simultaneously—which could lead to studies inadvertently being displayed multiple times. Instead, if the user clicks on a disabled duplicate data source, they receive a prompt asking to switch their current data source to the disabled duplicate; hence, they only use one instance of that data source at any given time.




How to find and use duplicate data sources:



Locate the Duplicates list (if available) in the bottom right-hand corner of the Study Browser window.



Click **Yes** to assign the duplicate as the preferred data source or click **No** to cancel the preference change. Users may be prompted to log in, depending on the configured login options.





The new data source is displayed along the bottom of the Study Browser, and the previous data source is now listed in the Duplicates list in the bottom right-hand corner of the Study Browser.

Note: The site may not have any disabled duplicate data sources if there are no remote connections set up for NovaPACS. Users should always save (export) a copy of user preferences before switching data sources to a disabled duplicate as a failsafe. The NovaPACS Diagnostic Viewer automatically determines if a data source is a duplicate of another and needs to be disabled.

CHAPTER 2—USING THE STUDY BROWSER

The **Study Browser** functions as the control panel for the Diagnostic Viewer. The Study Browser displays information including study information, patient information, patient history, the number of series in a study, the number of images in a series, and applicable reports for studies. In addition to these details, the viewable columns can be set to show information such as studies with existing reports and the name of the reviewing radiologist.

The Study Browser is where users can search for and open their studies into the Image Viewer. Users can customize search filter options on the various worklist tabs in the Study Browser to quickly locate their studies.



Users also have access to multiple worklist tabs including the Search Tab, the Assigned Worklist tab, Advanced Worklist tabs, the Teaching Files tab, and the Preliminary Impressions tab. These worklist tabs offer users additional ways to search for studies and save frequently used search criteria. Users can streamline their workflow and reduce the time spent searching for studies using the various worklist tabs in the Study Browser.


2-1—FILTERED STUDIES LIST

The **Filtered Studies** list is the top list in the Study Browser window by default and displays the list of studies that match the information in the search filters.


Note: Lists can be dragged and dropped to new locations in the Study Browser.

The Filtered Studies list is on all tabs in the Study Browser.

The Filtered Studies toolbar lists items, and right-click menus vary depending on which tab is being used.



DISPLAY



Users can set the Filtered Studies list to display all studies for each of their patients or they can set this list to only show the most recent study for their patients.

If users set the Filtered Studies list to show only the most recent study for their patients, all prior studies are located only on the Patient History list in the Study Browser. The study displayed for the patient in the Filtered Studies list is the most recent study that matches the search filter criteria.

If the Filtered Studies list is set to display all studies for patients, the list shows all studies for patients that match the search filter criteria.


The Filtered Studies list displays information in columns, including (but not limited to) the following:

- **Local**—displays whether studies are cached locally.
- **Open**—displays each study that is open.
- **Last Name**—displays the patient's last name.
- **First Name**—displays the patient's first name.
- **Birth Date**—displays the patient's date of birth.
- **Modality**—displays the modality for the study.
- **Study Description**—displays the study description that was ordered or entered at the modality for the study.

Note: The Study Description column displays the Procedure Name field from Nova RIS, the Study Description from an applicable HL7 message, or the Study Description manually entered at the modality. Novarad Support can configure the system to automatically update the PACS Study Description field when edits/changes are made in Nova RIS. For more information about this configuration, contact Novarad Support.

- **Patient ID**—displays the patient's identification number.
- **Status**—indicates if the study is New, Reviewed, Checked Out, Verified, Imported, or Locked for Edit
- **Accession #**—displays the accession number (or uniquely identifying number from Nova RIS or a HIS/EHR system) for the study.
- **Procedure Status**—displays the Nova RIS status of the procedure.
- **Facility**—displays the facility where the study was performed.
- **Study Date**—displays the date of the study.
- **Report**—indicates if a report has been saved if created in Novarad products.
- **Preliminary Impression Status**—indicates if the study has an Agree, Disagree, Indeterminate, or New preliminary impression status.
- **Report Created By**— This column displays who created the report.

- **Triage**—displays triage values so users can quickly identify studies based on their triage values. For example, a CVA icon will appear for stroke cases.
- **QA**—displays which studies do not have orders that match the study. This will be displayed with a red QA flag.
- **Technologist**—displays the Technologist who performed the study.
- **Validated By**— displays the Technologist or Admin who validated the study.
- **Procedure Codes**— This displays the procedure code for the study. This helps users easily identify similar studies. For instance, if users are reviewing an extensive patient history for a specific procedure, they can search for the procedure code that corresponds to the exact procedure they are looking for. The Procedure Code is available on all study lists.
- **Patient Location**— shows where in the hospital the patient is currently located.




Note: Customize the columns that appear in the list by right-clicking the column heading bar and selecting/deselecting the heading labels in the drop-down list. Rearrange the order of the columns by clicking and dragging the column headings to the preferred locations. Auto-fit the column to its contents by double-clicking the right divider of the column heading.


INCLUDE / EXCLUDE FILTER

The **Include / Exclude** control is used with the Study filters as shown below. The control is found above each respective filter.

- If the control is set to **Include**, the control remains visible while setting the filter. When work in the filter is complete, the control is hidden.



- If the control is set to **Exclude**, the control remains visible while setting the filter, and the filter value(s) displays with a strikethrough.



- If the filter is set to **(All)**, the control does not display. If a filter is set to **Exclude** and the filter value is changed to **(All)**, the control reverts to **Include**.
- The filter controls default to **Include**.


THE INCLUDE / EXCLUDE STATE OF EACH FILTER IS SAVED WITH THAT FILTER AND PERSISTS WITH THE FILTER. ALWAYS SHOW NON-REVIEWED STAT STUDIES

Last Name	First Name	Study Date	Modality	Study Description
Patient	Anonymous	1/27/2014 10:2	CR	THORACIC SPINE
Patient	Sample CR	5/26/2010 1:22:	CR	FOREARM LT

Novarad Support can configure the NovaPACS products to always show non-reviewed STAT studies regardless of search filter settings. This ensures that users see

STAT studies when they come through to give them the proper prioritization.

The Always Show Non-Reviewed STAT Studies feature applies to all Novarad products such as the NovaPACS Diagnostic Viewer, the NovaPACS Admin Console, and the NovaPACS HTML5 Viewer.




Users see STAT studies regardless of their search filters. This prevents users from inadvertently hiding STAT studies while searching in NovaPACS.

In addition to seeing all non-reviewed STAT studies, users may also see STAT studies of all statuses if they meet the specified search criteria. For example, if users are searching for all CT studies, they see all CT studies including reviewed STAT studies as well as all non-reviewed STAT studies from all modalities. Users can sort non-STAT studies by clicking the Status column, but the STAT studies display at the top regardless.

Study Date	Modality	Study Description	Status	Priority
1/27/2014 10:28 AM	CR	THORACIC SPINE	Checked Out	STAT
5/23/2010 1:08 PM	CT	LUMBAR SOINE	Reviewed	STAT
6/29/2010 10:31 AM	CT	CARDIAC	Reviewed	STAT
5/23/2010 12:22 PM	CT	CT ANGIO CHEST	New	STAT

SEARCH TAB FILTERED STUDIES LIST


TOOLBAR



- **Select All**—Selects all the studies in the list
- **Add to Worklist**—Adds the selected study to the current user's Assigned Worklist
- **Open**—Opens the selected study in the Image Viewer
- **Open + Old**—Opens the selected study and the most recent related study in the Image Viewer
- **Open in Current Workspace**—Opens the selected study in the current open workspace in the Image Viewer
 - **Assign**—Allows users to assign studies to their own Assigned Worklist or another user's Assigned Worklist
 - **Studies per Data Source**—Allows users to select the number of studies displayed from their data source(s)


Note: Selecting a large number of studies could cause the product to perform more slowly.

ADVANCED WORKLIST FILTERED STUDIES TOOLBAR



See the descriptions above.

Note: Toolbar buttons may be disabled if a study is not currently selected.



SMART WORKLIST

The Smart Worklist feature allows users to view a study's due date based on contractual agreements and study protocols. The Smart Worklist displays a due date and time in the NovaPACS Study Browser to communicate the sense of urgency for applicable studies. Once a new study is routed from the modality, the image processor calculates a due date. When a new study routes, the turnaround time populates based on the Apply Contract rule(s), as designated in the Admin Console.

After configuration in the Admin Console, users can view the Due column by right clicking the column heading. Once the drop-down list appears, select **Study**, then click **Due** to view due dates, as applicable. To sort the studies by due date, click the column header.

Note: The Due column is hidden by default and must be selected to display correctly.

Users continue to see STAT studies regardless of the search filters and Due column. This prevents users from inadvertently hiding STAT studies. The Study Browser now calls attention to studies that are due soon, or are overdue, by highlighting the entire row golden yellow, to differentiate between STAT studies, which appear red. If the user has selected the STAT Studies on Top site setting, STAT studies appear first, then due/overdue studies appear directly under the STAT studies.

Local	Open	Modality	Study Date	Study Description	Report	Facility
		MR	1/27/2014 10:28:11 AM	MRI TSPINE W/WO		NR-TRAINING
		MR	1/27/2014 10:28:11 AM	THORACIC SPINE		NR-TRAINING

Additionally, when users select Open Next Worklist Item, NovaPACS automatically prioritizes due and overdue studies after STAT studies and before other studies.

2-2—PATIENT HISTORY LIST

The Patient History list displays all completed studies for the selected patient including all time frames and modalities.


Local	Open	Modality	Study Date	Study Description	Report	Facility
		CR	1/27/2014 10:28:11 AM	MRI TSPINE W/WO		NR-TRAINING
		CR	1/27/2014 10:28:11 AM	THORACIC SPINE		NR-TRAINING

This list appears on the left side of the Study Browser window, just below the Filtered Studies list, by default. The Patient History list is on all tabs in the Study Browser (i.e., the Search tab, Assigned Worklist tab, and Advanced Worklist tabs).

Note: Lists can be dragged and dropped to new locations in the Study Browser. Lists can also be toggled on/off at any time.

The Patient History list displays information in columns, including (but not limited to) the following:

- **Local**—displays whether studies are cached locally
- **Open**—displays each study that is open
- **Last Name**—displays the patient's last name



- **First Name**—displays the patient's first name
- **Birth Date**—displays the patient's date of birth
- **Modality**—displays the modality for the study
- **Study Description**—displays the study description that was ordered or entered at the modality for the study

Note: The Study Description field displays the Procedure Name field from Nova RIS, the Study Description from an applicable HL7 message, or the Study Description manually entered at the modality. Novarad Support can configure the system to automatically update the PACS study description field when edits/changes are made in Nova RIS. For more information about this configuration, contact Novarad Support.

- **Patient ID**—displays the patient identification number
- **Status**—indicates if the study is New, Reviewed, Checked Out, Verified, Imported, or Locked for Edit
- **Accession #**—displays the accession number (or uniquely identifying number from Nova RIS or a HIS/EHR system) for the study
- **Procedure Status**—displays the Nova RIS status of the procedure
- **Procedure Code**—displays any Procedure Codes

Note: If Procedure Codes and relationships have not been added on the Procedure Codes tab of the NovaPACS Admin Console, there may not be any information in the Procedure Codes display column of the Patient History List.


- **Facility**—displays the facility where the study was performed
- **Study Date**—displays the date of the study
- **Report**—indicates if a report has been saved—if created in Novarad products
- **Preliminary Impression Status**—indicates if the study has an Agree, Disagree, Indeterminate, or New preliminary impression status

Note: Customize the columns that appear in the list by right clicking the column heading bar and selecting/deselecting the heading labels in the drop-down list. Rearrange the order of the columns by clicking and dragging the column headings to the preferred locations. Auto-fit the column to its contents by double-clicking the right divider of the column heading.

- **Select All**—Selects all the studies in the list
- **Add to Worklist**—Adds the selected study to the current user's Assigned Worklist
- **Open**—Opens the selected study in the Image Viewer

- **Open + Old**—Opens the selected study and the most recent, related study in the Image Viewer
- **Open in Current Workspace**—Opens the selected study in the current open workspace in the Image Viewer

2-3—SERIES LIST




The **Series** list displays all the series that are part of the selected study. The Series list appears below the Filtered Studies list and to the right of the Patient History list in the Study Browser Window, by default. The Series list is on all tabs in the Study Browser (i.e., the Search tab, Assigned Worklist tab, and Advanced Worklist tabs).

When NovaPACS receives a finalized, PDF DICOM report from Nova RIS, it displays as a series in the study titled Novarad Diagnostic Report.

Note: The Novarad Image Processor creates a new placeholder study if one

doesn't exist or if there is a diagnostic report that doesn't match an existing NovaPACS study. Admins can then go back later and fix the studies once they find out why there wasn't a matching NovaPACS study for the report.

The series details such as series description and image count display in the Series list. The Series list displays information in columns, including (but not limited to) the following:



- **Critical Findings**—displays an indicator if a radiologist has created a critical finding for the selected study/series
- **Description**—displays the series description, generated by the modality
- **Image Count**—displays the number of images in each series for the selected study
- **Modality**—displays the modality for the study
- **Series ID**—displays the series identification that provides a unique identification for each series within the study (generated by the modality)
- **Series UID**—displays the unique identifier for each series for the selected study (generated by the modality)
- **Study UID**—displays the unique identifier for the selected study (generated by the modality)


- **Select All**—Selects all the studies in the list
- **Open**—Opens the selected study in the Image Viewer

2-4—PRELIMINARY IMPRESSION LISTS

The NovaPACS Diagnostic Viewer now supports a primary/secondary read workflow where users can review images and record preliminary impressions and where radiologists can view and agree/disagree with the findings.

RECORDING PRELIMINARY IMPRESSIONS


1. First, users need to enable the Preliminary Impression list. Click **View** in the Study Browser and select **Preliminary Impressions Tab**.



Note: By default, this option is disabled in the Study Browser.

2. Once selected, the Preliminary Impression list appears.
3. To add a preliminary impression, select a study in the Study Browser, and click the **Add icon** in the Preliminary Impression list.
4. A Preliminary Findings text box displays. Click in this text box and type the preliminary impression.

5. Click **Save** to record it. Click **Cancel** to cancel out of creating a preliminary impression.




6. The preliminary impression now displays on the Preliminary Impression list under the Preliminary Findings heading. This has a time and date stamp for the user that created the finding. Users cannot edit a recorded preliminary finding.


Note: Once the preliminary finding is saved and recorded, the study has a New Preliminary Impression Status—so users can search for studies that have preliminary impressions. This is discussed in greater detail in later sections.

LOCATING STUDIES WITH PRELIMINARY IMPRESSIONS

There is a Preliminary Impressions worklist tab that is enabled by clicking **View** → **Preliminary Impressions Tab**.



There is also a **Preliminary Impression Status** search filter that can be added to the **Search tab**, **Advanced Worklist tabs**, and the **Preliminary Impression tab**.



Users can also add a **Preliminary Impression Status** display column in the **Study Browser**.


Studies that have a **Preliminary Impression** waiting for approval have a **New Preliminary Impression** status.

REVIEWING PRELIMINARY IMPRESSIONS


Modality	Status	Preliminary Impression Status
MG	New	
CR	Reviewed	New
ECG	Checked Out	
US	New	New
XA	Checked Out	
CR	Reviewed	

1. Enable the **Preliminary Impression** list by clicking **View** in the **Study Browser** and selecting **Preliminary Impressions Tab**.
2. Once selected, the **Preliminary Impression** list appears.
3. Search for and select the study with the preliminary impression.


4. Verify that the selected study has a preliminary impression by viewing the impression on the **Preliminary Impression** list, then click the **Add icon**.



5. A new text box displays. Click in the **Review Comments** text box and type the preliminary impression review comments.




6. Select a status from the **Status** drop-down list (Agree, Disagree, or Indeterminate).



7. The preliminary impression review is now displayed on the Preliminary Impression list under the Preliminary Findings heading and below the initial preliminary impression. This has a time and date stamp for the user that created the review.

8. Click **Save** to record the review or **Cancel** to cancel out of creating a preliminary impression review.

Note: Preliminary finding reviews cannot be edited once they have been saved. All preliminary impressions and reviews for the selected study display on the Preliminary Impression list. Users can see the preliminary impression review information on the Preliminary Impression Activity list in the Study Browser.



9. Once the preliminary impression review is saved and recorded, the study has the Preliminary Impression Status of what the radiologist selected (Agree, Disagree, Indeterminate, or New)


LOCATING REVIEWED PRELIMINARY

IMPRESSIONS

Users can locate a reviewed Preliminary Impression, review the comments, and add additional review comments to the Preliminary Assessment, if applicable. This enables the user to correct their initial impression based on feedback from the radiologist.


How to locate Reviewed Preliminary Impressions in a study:

1. First, enable the Preliminary Impression list by clicking **View** in the Study Browser and selecting **Preliminary Impressions Tab**.



2. Search for and select the study with the preliminary impression and review using the Preliminary Impression Status search filter in the Study Browser.


Users can search by Agree, Disagree, or Indeterminate Preliminary Impression status to find a reviewed Preliminary Impression.



Users can also select multiple Preliminary Impression statuses to broaden their search.


2-5—IMAGES TAB

The **Images** tab displays all the series preview thumbnails in the list that are separate from the other lists/tabs.



The Image list is located as a tab next to the Series list in the Study Browser by default.


At any time, the user can double-click a series preview thumbnail in the Images list to open that series in the Image Viewer. The study/series display is dependent upon the selected modality's hanging protocol layouts.



Users can also left click and drag a series preview thumbnail image into the Image Viewer at any time.

2-6—REPORT TAB


The **Report** tab displays the selected study's report preview in a tab that is separate from the other lists/tabs. The



Report tab appears below the Filtered Studies list and to the right of the Patient History list in the Study Browser Window by default.

Print—opens a print dialog allowing the user to print the report.

Clicking **Print Report** opens a standard print dialog box.



Note: When users sign and save reports in NovaPACS, the preliminary report displays in all report locations depending on configuration. Once the report is finalized in Nova RIS, a DICOM-encapsulated PDF report is sent to NovaPACS and added to the study. This DICOM PDF report displays in any location where users access reports.


2-7—STUDY BROWSER TAB FUNCTIONALITY

Each list or section in the Study Browser has a tab at the top with its name listed. Users can click and drag tabs to rearrange them.

DRAG AND DROP LISTS/TABS


The drag-and-drop functionality gives users the option to move each list/tab (Filtered Studies, Patient History, Series, Images, Report, and Preliminary Impressions lists/tabs) to a new location in the Study Browser.

1. Click and drag the list/tab.




2. As the list is dragged, hover the cursor over the docking station icon for any section in the Study Browser.

Users can drop the list to the top, bottom, left, or right of the current docking section using the docking station icon.




Users can also drop the list in the middle of the docking station icon, which groups the list as a new tab with the current tab(s).


3. Once the destination is selected, release the mouse button and the list appears in its new location.



Note: Users may drag and drop their lists to any location and/or new window in the Study Browser. However, the ability to drag lists into independent windows and maintain an optimal display depends upon the number of monitors used. Using a single monitor to display the Study Browser's lists each as independent windows may not display all relevant details optimally.





2-8—HANGING PROTOCOL LAYOUTS



Hanging Protocol Layouts are the saved layouts associated with a specific Hanging Protocol for displaying study images in the Image Viewer.


If activated, the **Hanging Protocol Layouts** section is located at the bottom of the Study Browser window and displays thumbnail images of the layouts defined for the currently opened study. This section only displays content in the Study Browser if a study is open in the Image Viewer.

The title of this section also lists the current Hanging Protocol. This appears on all tabs in the Study Browser (i.e., the Search tab, Assigned Worklist tab, and Advanced Worklist tabs).




To activate the Hanging Protocol Layouts section in the Study Browser window, click **View** in the Study Browser and select **Hanging Protocol Layouts**. Users can also press the default hotkey **F4**.

The Hanging Protocol Layouts section collapses by default at the bottom of the Study Browser. To view the current study's hanging protocol layouts, hover the cursor over **Hanging Protocol Layouts**.



The Hanging Protocols list slides out to display as a panel.


The Hanging Protocols list displays until users click anywhere else in the Study Browser. Interacting with other Study Browser features collapses the list.



At any time, users can change the panel's slide/collapse functionality to a pinned list by clicking the **Pin icon** in the top, right-hand corner of the Hanging Protocols list. It can be unpinned and returned to a sliding panel at any time using the same pin icon.

To use Hanging Protocol Layouts, select a layout to apply it to the current study in the Image Viewer.


Note: Press the number keys on the keyboard to switch to the corresponding layout in this list.



Users can Create, Remove, Duplicate, and Edit layouts using the links on the Hanging Protocol Layouts list.

2-9—THE SEARCH WORKLIST

The default worklist on the **Study Browser's** home window is the Search worklist. The Search worklist is used with the filters on the page to locate specific information.




Note: To the right of the Search worklist is an Assigned Worklist that can be shown/hidden using the View menu. There is also an Advanced Worklist that can be used to create advanced worklists with custom settings.

CUSTOMIZING THE FILTERS GROUP


The default settings for the filters vary depending on each user's role. Users can add or remove filter fields by right clicking the filter's bar and selecting the appropriate option in the drop-down list. Filters are organized by related area, and active filters on the screen have a checkmark beside them.

Note: The Saved Filters and Data Source fields are not optional and do not appear in the optional filters list.




CUSTOMIZING AND SAVING FILTERS

The Saved Filters drop-down list allows users to create and save custom filters for settings used on a regular basis.



1. Select the preferred settings by right clicking the filters and selecting options from the drop-down lists.


2. Select **Edit Saved Filters** from the **Saved Filters** drop-down list.



3. Enter a name for the new filter in the **Save Current Filters As** field.
4. Click **OK** to save the filter. Click **Cancel** to cancel out of creating a saved filter.


APPLYING A SAVED FILTER

1. Click the **Saved Filters** drop-down list.
2. Select the saved filter from the list.





DELETING A SAVED FILTER

1. Click the **Saved Filters** drop-down list.
2. Select **Edit Saved Filters**.
3. In the **Select Saved Filters to Delete** box, click to check the box next to the unwanted filter.
4. Click **OK** to delete the filter. Click **Cancel** to cancel out of deleting the filter.




USING THE SEARCH FILTERS

The filters are a group of search fields just below the Search/Worklist tabs along the top of the Study Browser Window. These filters allow users to select the information to display in the Filtered Studies list. Patient searches can be conducted more quickly and efficiently using the filters in this tab. Users can search for studies using any combination of search filters. As information is entered into the filters, the Study Browser automatically returns search results. Users can disable the automatic search function of the Study Browser and configure the search filters to only return results after clicking **Search**. Disabling automatic search could be especially useful for improving search speed and efficiency for facilities with large amounts of customer data.




Note: The Auto Load Results checkbox is located in the search filter section under the Search button. This option is selected by default. If users deselect Auto Load Results, they must click **Search** to manually search for their studies. Otherwise, users are prompted to do so. Users can select/deselect Auto Load Results at any time.

THE DATA SOURCE FILTER



The Data Source Filter allows users to filter a search using the various sources that may contain patient data. Click the **Data Source** drop-down list and select one or more of the data source checkboxes. Users can also select **All** to select all available data sources.




Data Sources display at the bottom right of the **Study Browser** window. The indicator light to the left of the Data Source display provides a visual indicator of the status of that server.

A green icon indicates the data source connection is active and the status is available, and red signifies the data source is unavailable or disconnected. Hover the mouse cursor over the Data Source display at the bottom of the Study Browser window to see the connection status and click to see the Edit Data Source dialog box.

Users may have DICOM Pass Through enabled in NovaPACS. This allows users to securely view images and data from off-site archives and PACS servers by remotely accessing them through NovaPACS. DICOM Pass Through servers are set as a data source that interacts normally with NovaPACS. Users can search for or filter out results from a DICOM Pass Through server the same way they would any other data source.

THE TIME FRAME FILTER

The Time Frame Filter allows users to select specific time periods to filter a search. Users can search all time frames by selecting **All** in the Time Frame Filter. Users can also select from the time frames listed in the drop-down list or create a custom time frame for this filter.




Note: If users deselect the Time Frame filter from the search filter list, the time frame defaults to the last three hours and only displays studies added in that time frame.


To create a custom time frame:

- Click the **Time Frame** drop-down list.
- Select **Custom** at the bottom of the list.

- c. Click the arrow for the drop-down calendar to choose a specific date for the Start Time and use the arrows next to the time to set the exact starting time frame.
- d. Click the arrow for the drop-down calendar to choose a specific date for the End Time and use the arrows next to the time to set the exact ending time frame.
- e. Click OK to apply the new custom Time Frame Filter.



THE MODALITY FILTER




The Modality Filter allows users to search for data from a specific modality. Users can select just one modality to search or select several modalities.

Click the **Modality Filter** drop-down list and click the checkboxes next to the modality options to select or deselect them.

Note: The server selected as the current data source provides the modalities that are on the Modality Filter drop-down list. If multiple data sources are selected, the Modality Filter drop-down list contains a combination of all modalities from all data sources.

THE RADIOLOGIST FILTER

The Radiologist Filter allows users to filter studies that are associated with a Radiologist.




Enter the Radiologist's username, last name, or first name. All studies that contain those letters or that last name are returned in the search list.

Note: Users cannot search for a radiologist with both a first and last name. Search for either last name or first name.

THE REFERRING PHYSICIAN FILTER


The Referring Physician filter allows users to filter studies that are associated with a Referring Physician.



Enter the Referring Physician's username, last name, or first name. All studies that contain those letters or that last name are returned in the search list.


Note: Users cannot search for a Referring Physician with a first and last name combined. Search for either last name or first name.

THE STATUS FILTER



Users can filter results based on the study's status. Click the **Status Filter** drop-down list and choose from **All**, **Checked Out**, **Imported**, **Locked for Edit**, **New**, **Reviewed**, and/or **Verified**.


THE CHECKED OUT FILTER



The **Checked Out Filter** allows users to filter results based on who has checked out the studies.

Enter a username. All studies that contain that username are returned in the search list.

THE PATIENT GROUP FILTER



The **Patient Group filter** allows users to filter by Patient Group.


Click the **Patient Group** drop-down list and click the checkboxes next to the patient groups to select or deselect them.

THE STUDY UID FILTER

The **Study UID filter** allows users to search for studies by their study unique identification (UID).

Enter all or part of the Study UID in the search field. All studies that contain those numbers are returned in the search list.


THE PATIENT ID FILTER



The **Patient ID filter** allows users to search for studies by the patient identification number (PID).

Enter the first few characters or the complete Patient ID of the patient. All studies associated with a Patient ID beginning with the characters entered are returned in the search list.


THE LAST NAME, FIRST NAME, AND ACCESSION # FILTERS



The **Last Name, First Name, and Accession # filters** all function in the same manner. By entering the first few

characters of the patient's last name, first name, or accession number, the filters return all results that contain those letters or numbers in the Filtered Studies list.


THE ASSIGNED TO FILTER



The Assigned To filter allows users to filter studies by user for studies that have been assigned.

Click the **Assigned To** drop-down list and select an option.

PRELIMINARY IMPRESSION STATUS FILTER




The Preliminary Impression Status filter allows users to filter studies by Preliminary Impressions status.

Click the **Preliminary Impression Status** drop-down list and choose from **All**, **Agree**, **Disagree**, **Indeterminate** and/or **New**.


STUDIES PER DATA SOURCE FILTER

The Studies per Data Source filter allows users to select the number of studies to view from each source.

Click the arrow next to the drop-down list and select the quantity that best fits the search needs (50, 100, 250, 350, 500, 1000, or 2000).




TEACHING FILE SEARCH FILTERS



If **Teaching Files** have been enabled, users can add Teaching File search filters to the worklist tabs.

Note: To learn how to enable this feature and to learn more about Teaching Files, see **Chapter 5, section 5-17—Teaching File**.

Right click the search filter area, hover the cursor over **Teaching File**, and select a Teaching File search filter option from the drop-down list.




Users must set the Teaching File search filter to **Yes** before they can search for saved **Teaching Files**.

Note: When the Teaching File search filter is set to **Yes**, users can only search for Teaching Files and not original studies.

Users can also select **All** to search for both teaching files along with real studies.

STUDY OF INTEREST SEARCH FILTER



Users can also add the **Study of Interest** filter to search for studies that have been marked as *Studies of Interest*—a flag that signifies that users want to come back to it later or to create it as a *Teaching File* study at a later time.

Right click the search filter area, hover the cursor over *Study*, and select **Study of Interest** from the drop-down list.


The **Study of Interest** filter appears in the filters bar, and users must set the filter to **Yes** to search for studies of interest only.

Users can also select the **All** option to mix both studies of interest along with other studies.


AUTO-HIGHLIGHTED TEXT IN SEARCH FILTERS

When users press **Tab** to navigate between search filters, any previously entered text is automatically highlighted. This enables users to quickly type over or delete text.

1. Log into the NovaPACS Diagnostic Viewer and access the Study Browser.
2. Click a text input field in the search filters section such as the **Patient ID**, **Last Name**, or **First Name** search filter.
3. Press **Tab** to navigate to the next search filter.




Any pre-entered text is automatically highlighted.



Note: The Auto-Highlight option only works on text input search filters. Information is not highlighted on drop-down search filters such as *Data Sources* or *Time Frame*.

2-10—THE ASSIGNED WORKLIST

The *Assigned Worklist* is located to the right of the *Search* worklist when activated and is used primarily to display studies that have been assigned to the currently logged in user.



The *Assigned Worklist* is an effective tool for large groups of users that read images at the same time from a *Filtered Studies* list. The *Assigned Worklist* shows studies that have been added by the user to their own *Assigned Worklist* to view later, or that have been assigned to the user's *Assigned Worklist* by someone else. Radiologists can move studies


into their worklist for review or assign studies to another user's worklist for specific users or larger groups to review.

DISPLAY

The Assigned Worklist displays information in the same way information is displayed in the Search tab—with a Filtered Studies, Patient History, and Series list. Some menu items are not available on the Assigned Worklist, including Add to Worklist.

These options are available on both the Search tab and the Advanced Worklist tabs.

ACTIVATING THE ASSIGNED WORKLIST




Users can activate the Assigned Worklist tab by selecting **View** in the Study Browser and selecting **Worklist**. Users can also press the default hotkey **F2**. The Assigned Worklist tab appears next to the Search tab in the Study Browser.

The Assigned Worklist also appears automatically when a study is added to the worklist.

USING THE ASSIGNED WORKLIST SEARCH FILTERS

The filters are a group of search fields just below the Assigned Worklist tab along the top of the Study Browser window. These filters allow users to select the information to display in the Filtered Studies list. Patient searches can be conducted more quickly and efficiently using the filters in this tab. Users can search for studies using any combination of search filters.

THE DATA SOURCE FILTER



The Data Source filter allows users to filter a search using the various sources that may contain patient data. Click the **Data Source** drop-down list and select one or more of the data source checkboxes. Users can also select **All** to select all available data sources.

THE TIME FRAME FILTER

The Time Frame filter allows users to select specific time periods to filter a search. Users can search all time frames by selecting **All** in the Time Frame filter. Users can also select from the time frames listed in the drop-down list or create a custom time frame for this filter.

To create a custom time frame:

1. Click the **Time Frame** drop-down list.
2. Select **Custom** at the bottom of the list.

3. Click the arrow for the drop-down calendar to choose a specific date for the Start Time and use the arrows next to the time to set the exact starting time frame.
4. Click the arrow for the drop-down calendar to choose a specific date for the End Time and use the arrows next to the time to set the exact ending time frame.
5. Click **OK** to apply the new custom Time Frame filter.

THE STATUS FILTER


The Status filter allows users to filter results based on the status of the studies being searched. Click the **Status Filter** drop-down list and choose from **All**, **Checked Out**, **New**, **Reviewed**, and/or **Verified**.

THE ASSIGNED TO FILTER

The Assigned To filter allows users to filter studies by user for studies that have been assigned.

Click the **Assigned To** drop-down list and select an option.


2-11—ADVANCED WORKLISTS




Users can create sets of filters for specific conditions such as reading certain types of exams or reading at certain locations. This eliminates the need to change filters on the Search worklist constantly. These preferred filters are created using the Advanced Worklist option. The Advanced Worklist is located to the right of the Study Browser's Search worklist. The worklist is either labeled Advanced Worklist or a blank worklist.

CREATING ADVANCED WORKLISTS

1. Select the blank or empty tab to the right of the Search worklist.




2. Select **New** and the options window expands.
3. Create a name for the Worklist by typing it into the Worklist Name field.
4. Select the preferred criteria in each search field using the drop-down lists.



Note: To set a second set of search criteria, click + below the search criteria columns. An additional row of search criteria appears; set as many rows of search filters as needed. To delete a row of search criteria, click – to the left of the row.

Click **Save** at the bottom of the window. The Advanced Worklist tab appears to the right of the Search tab with the saved name. A new blank worklist tab appears next to the Advanced Worklist tab in the Study Browser window.



6. To access any Advanced Worklists that are no longer displayed, click the **Advanced Worklist** tab, click the **Advanced Worklists** heading, and select a listed Advanced Worklist.

EDITING ADVANCED WORKLISTS

1. Click **Advanced Worklist**.
2. Click the **Advanced Worklists** drop-down list. Select the Advanced Worklist to be edited.
3. Select **Edit**.
4. The Advanced Worklist window expands. Edit the search criteria as needed using the drop-down lists.
5. Click **Save** at the bottom of the window.

DELETING ADVANCED WORKLISTS

1. Click the **Advanced Worklist**.
2. Click the **Advanced Worklists** drop-down list. Select the Advanced Worklist to be deleted.
3. Select **Delete**.
4. A Delete Advanced Worklist confirmation box appears to confirm that the Advanced Worklist is being deleted. Click **OK** to delete the worklist.

2-12—CUSTOMIZING STUDY BROWSER LISTS

Users can customize options in the Study Browser lists, including column options, resizing options, and sorting or changing column orders.

RESIZING COLUMN HEADINGS

Users can resize the width of a column heading in any list to better accommodate the information within.

1. Hover the cursor over the right divider of the column heading. The cursor changes to a double-ended arrow.
2. Double-click the right divider of the column heading and the column auto-fits to its contents. Users can also click the right divider of the column heading and drag it to the right or left until the column is the desired size.

RESIZING LIST HEIGHTS

Users can resize the height of the Study Browser lists to display more or less information.

1. Hover the cursor over the divider at the edge of the list to be resized. The cursor changes to a double-ended arrow.
2. Click and drag to the preferred height.

CHOOSING COLUMN HEADINGS


1. Right-click the column heading bar to customize. A drop-down list containing the column heading options appears.
2. Click to select or deselect the column headings to appear in the list.

SORTING DATA IN A COLUMN

1. Click the column heading to automatically alphabetize the information (or sort numerically if applicable). It automatically lists in ascending order.
2. To reverse the sort order to descending, click the column heading a second time.

CHANGING THE ORDER OF COLUMNS

Users can change the order in which the columns appear in the lists. To change the column order, complete the following steps:



1. Click and drag the column heading to the preferred location and release.


Note: As the cursor passes over the other column headings, it highlights a dark blue—representing the new location of the column heading when released.

HIDE/SHOW ALL DISPLAY COLUMNS

1. In the Filtered Studies or Patient History list, right-click any display column.

A drop-down list of display column options appears.

The Show All and Hide Hideable options are at the bottom of the display column drop-down list.



2. To display all possible display columns, select **Show All**.

All display columns are added to that list.

3. To quickly remove all display columns except the Local and Open columns, select **Hide Hideable**.

All display columns except the Local and Open display columns are disabled and removed.

Note: The Local and Open display columns always display and cannot be removed.


2-13—USING THE STUDY BROWSER MENUS/CUSTOM TOOLBAR

The Study Browser's Menus and Custom Toolbar have many helpful functions. The following sections provide brief information on the Study Browser's Menu and Custom Toolbar.

WHERE TO FIND ADDITIONAL STUDY BROWSER MENU INFORMATION


The Study Browser Menu contains available actions for the Study Browser and Image Viewer.

USING THE STUDY BROWSER CUSTOM TOOLBAR



The Study Browser Toolbar is just below the menu bar at the top of the Study Browser window. The default settings for the Study Browser toolbar vary by role. Click any button on the toolbar to use its function.


CUSTOMIZING THE STUDY BROWSER'S CUSTOM TOOLBAR



Users can customize the Study Browser's Custom Toolbar to their preferences. Users can add or remove toolbar items and change the way the toolbar appears in the window.

Click **Settings** in the Study Browser and select **Customize Toolbar**. Users can also right-click any toolbar in the Study Browser window and select **Customize** from the drop-down list. The **Customize Toolbar** window opens.


2. To add an item to the toolbar, select the item to add from the Available items list on the left of the screen. Click **Add** to add it to the Current toolbar items list on the right of the screen. Users can also select items in the Current toolbar items list and click **Remove**.



3. Click **OK** to save and exit.

2-14—ACTIVATING THE STUDY BROWSER


Users can hide the **Study Browser** window when using the **Image Viewer** to allow a full-image display.




1. Move the mouse to the right or left edge of the Image Viewer window.

Note: The Study Browser auto-shows and can be displayed if Auto-show Study Browser is selected on the Options page.

2. Click the **Study Browser icon** on the Image Viewer toolbar.



3. Press the default hotkey **Spacebar** on the keyboard to close all workspaces in the Image Viewer and show the Study Browser.




2-15—CUSTOMIZING THE STUDY BROWSER COLOR AND FONT SIZE



The settings gear icon is located in the top right-hand corner of the Study Browser. Users can click this icon to customize the color theme and font size options in the Study Browser.

STUDY BROWSER COLOR THEMES




Users can quickly access Study Browser display options when needed using the gear icon. When not in use, these options only take minimal space.

Users have access to six color themes including dark gray, dark blue, pink, classic, light blue, and green. Hovering the cursor over a color displays a tooltip with the theme name.

The Classic theme enables users to change their display to a look and feel they are more accustomed to with earlier software versions.

Note: Once selected, the chosen color scheme saves to the user's preferences.

FONT SIZE ADJUSTMENT SLIDER



Users can adjust the font size of the Study Browser window using a font size slider.

By moving the font size slider to the left, the text displays smaller.

If the slider is moved to the right, the text displays larger.


Note: The Classic Theme is merely similar to the previous Study Browser look and is not an exact replica of it. The toolbar icons revert to icons used in versions prior to 8.3. Study Browser themes and displays are dependent upon monitor display and resolution settings.

INTERFACE SCALING


Click the **Gear icon** on the top right-hand side of the Study Browser.

The Viewer Scale has a drop-down box where users can select a scaling percentage.

The page adjusts to the selected scale.



CHAPTER 3—USING THE STUDY BROWSER FILE MENU



One main function of the Diagnostic Viewer’s **Study Browser** is to help users search for and view studies. The **File** menu in the Study Browser allows users to open studies to view or to exit the Diagnostic Viewer. The **File** menu is on the Study Browser’s main toolbar.

3-1—OPEN FROM DIRECTORY

Open from Directory allows users to open a group of DICOM image files from a directory that is on a computer or CD-ROM.

1. Click **File** in the Study Browser and select **Open from Directory**. The Browse for Folder dialog box opens.
2. Browse for the folder by expanding the folders in the dialog box. Click **OK**.

Note: Opening images from a file or a directory does not import the study into the NovaPACS system. Studies must be imported on the Studies tab of the NovaPACS Admin Console.

3-2—OPEN FROM FILE

Open from File allows users to open a single DICOM image file from a directory on a computer or a CD-ROM.

1. Click **File** in the Study Browser and select **Open from File**. The Open dialog box opens.
2. Browse for the image file by expanding the folders in the dialog box. Click **Open**.


3-3—EXIT THE DIAGNOSTIC VIEWER

Exit allows users to exit from the Diagnostic Viewer application completely. If users do not fully exit the NovaPACS Diagnostic Viewer correctly, they can potentially lose some of their most recent customizations that would otherwise be saved to their user preferences. A full copy of the user’s preferences is written to the server upon correctly exiting the application.

Click **File** in the Study Browser and select **Exit**.

If the option for the Diagnostic Viewer to prompt before closing the application is selected, a warning appears. Click **Yes** to exit, or **No** to return to the Diagnostic Viewer.

CHAPTER 4—USING THE WORKLIST MENU



The Study Browser’s Worklist menu contains Worklist functions that help users organize and manage workflow. The Worklist menu is located on the Study Browser’s main toolbar. The items on the Worklist menu depend on which tab is used. The Search and Advanced Worklist tabs contain fewer worklist menu items than the Assigned Worklist tab.

4-1—INTRODUCTION TO THE ASSIGNED WORKLIST

The Assigned Worklist is a tab located to the right of the Search tab when activated and is used with the filters on the page to locate specific information. The Assigned Worklist is an effective tool for large groups of users that read images at the same time from a Filtered Studies list. The Assigned Worklist shows studies that have been added by the user to view later or that have been assigned to the user by someone else. Radiologists can move studies into their worklist for review or assign studies to another user’s worklist for specific users or larger groups to review.


4-2—ADD TO WORKLIST

Users can add a study to the Worklist to organize and manage workflow. Users can also assign a study to another user’s Worklist for them to review.

ADD TO WORKLIST

Users can add studies from the Search tab and the Advanced Worklist tabs to their Assigned Worklist. This allows users to see only specific studies along with studies assigned to them by other users for convenience.

1. Select a study from the Filtered Studies or Patient History list on the Search tab or an Advanced Worklist tab. To select multiple studies to add, hold the **CTRL** key while selecting patients.



2. Click the **Add to Worklist** icon on the toolbar.

The study appears in the Assigned Worklist tab.

Note: If the study added does not meet the search criteria in the search filters, it may not appear in the Filtered Studies list on the Assigned Worklist. Review the search filters if a study is not visible in the worklist.


Users can also add selected studies to their Assigned Worklist

- Click **Worklist** in the Study Browser and select **Add to Worklist**.
- Right-click the study and select **Add Selected Studies to My Worklist** from the drop-down list.
- Users can customize a hotkey to assign a study to the Worklist.

ASSIGN

Users can assign studies to another user for review. Studies can be assigned to other users from any tab in the Study Browser (Search tab, Assigned Worklist tab, and Advanced Worklist tabs).


1. Select the study to assign from the Filtered Studies list in the Study Browser. To select multiple studies, hold **CTRL** while clicking.
2. Click **Assign** on the Filtered Studies list toolbar. The Assign Studies from _____ dialog box opens.



3. In the Assign Studies from _____ dialog box, select an option from the What role is the assignment for drop-down list.

Note: The **What role is the assignment for?** section populates with users that have that selected role.

4. Select a user from the Which user to assign to section and click **OK**.



Note: If users want the Worklist items to have a status other than **New**—so other users in the Worklist group do not see the patients in the Worklist—users can click **Checkout** for the item before or right after moving it to their Worklist.

Users can also assign selected studies to other users.

- Click **Worklist** in the Study Browser and select **Assign**.
- Right-click **Assign Selected Studies to Another User** from the drop-down list.

Users can customize a hotkey to assign a study.

4-3—OPEN WORKLIST ITEM

Users can open Worklist studies in the Image Viewer.


1. Select a Worklist study.
2. Click **Open** on one of the toolbars in the Study Browser. (The Filtered Studies, Patient History, and Series lists all have the Open toolbar button). Click **Open + Old** on the Filtered Studies or Patient History toolbar to open the selected study with prior studies. This opens in the Image Viewer.

Users can also open selected studies:

- Click **Worklist** in the Study Browser and select **Open Worklist Item**.
- Double-click a study on the Assigned Worklist tab.
- Right-click a study and select **Open Selected Studies** from the drop-down list. Select **Open Selected Studies with Older Related Studies** to open the selected study with prior studies.

Users can customize a hotkey to open a study.

4-4—NEXT WORKLIST ITEM



Next Worklist Item allows users to close all currently opened studies and open the next worklist study in the Image Viewer.


Click **Worklist** in the Study Browser and select **Next Worklist Item**. Users can also press the default hotkey **F1**.

4-5—REMOVE WORKLIST ITEM

Users can remove an item from their Assigned Worklist by selecting the study and clicking **Remove Worklist Item** from the Worklist menu.


4-6—CLEAR WORKLIST

Users can clear their Assigned Worklist using **Clear** on the Worklist menu.



Note: Users may want to ensure they don't accidentally clear their worklist by turning on **Prompt before clearing worklist** on the Options page.

CHAPTER 5—USING THE STUDY BROWSER STUDY MENU



The **Study** menu helps radiologists find additional study information sources and tools and manage caseloads. The Study menu is located on the Study Browser's main toolbar.

5-1—SELECT ALL

Select All allows users to select all of the items on the active list (Filtered Studies, Patient History, or Series List).

1. Click a list in the Study Browser to make it the active list.
2. Click **Study** in the Study Browser and click **Select All**.

The following are alternate ways to use the **Select All** option:

1. Click a list to make it the active list, then right-click the list and choose **Select All** from the drop-down list.
2. Click a list to make it the active list, then click **Select All** on the toolbar for that list.
3. Click to make a list active and press the default hotkey **CTRL + A**.
4. Click anywhere on the list to deselect the items after using the **Select All** function.


5-2—OPEN

Open allows users to open a selected study, series, or image in the Image Viewer.

Click the desired study, series, or image and click **Study** in the Study Browser. Select **Open**.

The following are alternate ways to use the **Open** option:

1. Click the desired study, series, or image from any list in the Study Browser, and click **Open** on the toolbar for that list.
2. Double-click the study, series, or image and the Image Viewer automatically opens it.



Note: In order to prevent two users from changing a study at the same time, a message alerts the second user that the study/series is currently checked out. This warning also appears when an administrator is editing the study in the Admin Console.

5-3—CLOSE WORKSPACE

Close Workspace closes the current, visible workspace in the Image Viewer while leaving any other workspaces open.

Click **Study** in the Study Browser and select **Close Workspace**.

Note: There are other ways to use **Close Workspace**. Users can set a hotkey or add a button to the Study Browser's Custom Toolbar.

5-4—CLOSE ALL WORKSPACES

Close All Workspaces closes the current, visible workspace and all other open workspaces in the Image Viewer.

1. Click **Study** in the Study Browser and select **Close All Workspaces**.
2. Press the default hotkey **CTRL + Shift + F4**.

5-5—ROUTE


Route allows users to share images with another site. Images can be routed manually from a NovaPACS data source to a specified destination.

Note: NovaPACS can be set up to automatically route prior studies to workstations based on the RIS schedule for the next day. This allows users to view necessary information without workflow interruptions or delays. Contact customer support to enable this feature.

1. Select the desired study or series to route from any list in the Study Browser.

Note: Use the **Select All** function on the Filtered Studies list to select all studies to route or hold **CTRL** and click to select multiple studies or series to route.

2. Click **Study** in the Study Browser and select **Route**. The Add Route Queue Entry dialog box opens.



3. Choose the desired study, series, and/or image(s) to route by expanding the lists and clicking the corresponding checkboxes next to the items. If the checkbox next to the study is marked without expanding the list, all of the accompanying series and images are automatically routed.

Note: Route rules can be configured in the Admin Console. If one series in a multi-modality study matches the route rule, the entire study is routed (e.g., a mammography study can include images from multiple modalities like MG and US and also non-image files such as Structured Reports).

4. Choose a destination from the **Destination** drop-down list.

Note: An Administrator must add the Destinations that appear in the drop-down list.

5. Set the routing priority for the image routing by clicking the up or down arrows next to the **Priority** field. Routing priorities range from 1-100, with one as the highest priority. The default routing Priority is three.
6. Click **OK** to route to the specified destination.

Note: Users can customize a hotkey or add a button to the Study Browser's Custom Toolbar for the Route option.

Note: Finalized reports are sent to NovaPACS as a DICOM Encapsulated PDF report—with an image type of DOC. These reports are part of the NovaPACS study and therefore can be routed to any location the study can be sent. However, the ability for a receiving system to store and view the reports is dependent on their ability to support storage and viewing of DOC DICOM images. If users route to a Novarad site using an older version of software, they may not see the DICOM PDF report as part of the study.

5-6—DOCUMENTS

Documents allows users to see any documents associated with the patient file that were scanned into the RIS.

1. Select the study from any list in the Study Browser.
2. Click **Study** in the Study Browser and select **Documents**.

The Document Viewer window opens.

5-7—EXAM REPORT


Exam Report allows users to see the report associated with the patient file if one exists. Exam Reports are only available if the system is integrated with Nova RIS.

1. Select the study from any list in the Study Browser.
2. Click **Study** in the Study Browser and select **Exam Report**.

The report for the study opens in a new window.

Note: When reports are finalized in Nova RIS, a DICOM PDF report is sent to NovaPACS and is added to the study as a new series. If no matching NovaPACS study is found, a placeholder study is created for the report. Admins should locate this study in the NovaPACS Admin Console on the Studies tab and move the report series to the correct study manually.


The following are alternate ways to use the Exam Report option:



1. Select the study from any list in the Study Browser and press the default hotkey **CTRL + R**.
2. Users can add an Exam Report button to the Study Browser's Custom Toolbar.

Note: When users sign and save reports in NovaPACS, the preliminary report displays in all report locations depending on configuration. Once the report is finalized in Nova RIS, a DICOM encapsulated PDF report is sent to NovaPACS and added to the study. This DICOM PDF report displays in any location where users access reports.

5-8—MANAGE REPORT WORKFLOW



Manage Report Workflow allows users to quickly view their reports and report statuses. **Manage Report Workflow** is particularly useful for Administrators, Transcriptionists, and Radiologists. Users can review and/or sign reports using this option.

In the Report Workflow window, users see four main boxes for Dictated Reports, Unsigned Reports, Signed Reports, and Finalized Reports. In the Other window, users can see reports that have been marked Requires Correction and can look up physicians in the system using the Physician Lookup box.

To access the Report Workflow window, click **Study** in the Study Browser and select **Manage Report Workflow**.

Note: Users can customize a hotkey or add a button to the Study Browser's Custom Toolbar for the Manage Report Workflow option.

5-9—RIS INFORMATION

RIS Information allows users to view the RIS Information dialog box for a study. The information included in the RIS Information box is submitted by the Admin, Receptionist, or Technologist in Nova RIS. RIS Information dialog boxes include patient information, Referring Physician information, order notes, and procedure information about a study.

1. Select a study from any list in the Study Browser.
2. Click **Study** in the Study Browser and select **RIS Information**.

The RIS Information dialog box opens and displays any applicable Nova RIS information for the study.

Note: Users can customize a hotkey or add a button to the Study Browser's Custom Toolbar for the RIS Information option.


5-10—MATCH STUDY TO ORDER

The Match Study to Order feature can be used to match images with RIS orders. This allows Technologists to quickly see the information they need about the patient. Below are some of the cases where this feature will be useful:

- Allows the Technologist to verify the patient's name is spelled correctly.
- Allows them to see if the patient is scheduled twice in small bowel series on the same day. This lets them organize which procedure needs to happen first.
- It helps them see if CR images and Fluoro images of the same anatomy are placed in the same order and if they need to be merged, viewed, and read on the same report.
- Matching studies to orders makes it easier to organize images in the event the patient requires follow-up images.

Users can first enable the Scheduled Orders View to see what orders have been scheduled. To enable this view:

1. Right-click the toolbar and select **Customize**.
2. In the custom toolbars window, click the **View** dropdown menu.
3. Click **Scheduled Orders**, then click **Add**.



4. Click **OK**.

Users will then be able to view the scheduled orders.

Accession	First Name	Last Name	Modality	Description	Date
5916112H16098-BF	Shannon Leigh	Luna	US	US UNHS ULTRASOUND BREAST(D)	3/27/2025 9:06:00 AM
5930677H16098-BF	Tuesesa	Lewis	MG	Diagnostic - ComboHD	3/27/2025 1:23:00 PM
690446	Harsa	Warmbier	US	US Bilateral leg screening	3/27/2025 1:00:00 PM
690450	Aljandra	Warmbier	US	US Bilateral leg screening	3/27/2025 1:30:00 PM
691184	Tawnya	Teichert	Ultrasound	US OB less than 14 weeks	3/27/2025 1:00:00 PM
691734	Jeremy	Andrus	CT	CTA Heart and calcium scoring	3/28/2025 9:00:00 AM

To access the Match Study to Order Window, right-click the study and then click the **Match Study to Order** option.

Match Study to Order
— □ ×

Study Details

First Name: PET
 Last Name: TEST 4
 DOB: 1/1/1930
 Study Description: PET/CT SKULL TO THIGH
 Patient ID: PET 4
 Patient Group: UT-UCS
 Study UID: 1.2.840.113619.2.55.3.168430348.1143!
 Study Date: 6/16/2016 11:18 AM
 Accession:

Order Details

First Name:
 Last Name:
 DOB:
 Order Description:
 Patient ID:
 Patient Group:
 Procedure ID:
 Study Date:
 Accession:

Select Order or Create RIS Order

Accession	First Name	Last Name	Modality	Description	Date
691734	Jeremy	Andrus	CT	CTA Heart and calcium scoring	3/28/2025 9:00:00 AM
691990	Paul	Thomas	CT	CT Full Body/Head	3/24/2025 3:00:00 PM
692114	Aleffio	Galletta	CT	CT Full Body/Head	3/24/2025 9:00:00 AM
692940	Launna	Weiss	CT	CT Abdomen/Pelvis w/o contrast	3/24/2025 1:00:00 PM
693150	Stephanie	Fish	CT	CT Orbit, Sella, IAC w/o contrast	3/26/2025 9:30:00 AM
693400	Glen	Garcia	CT	CT Pelvis w/o contrast	3/27/2025 1:50:00 PM
693458	Olivia	Wright	CT	CT Lower Extremity w/o contrast (Foot)	3/27/2025 12:00:00 PM

Create RIS Order and Match
Match
Cancel

To match the study to the order, click on the order that matches the study details. Then click **Match**.

If the RIS order does not exist, click **Create RIS Order and Match** and the system will create the order in RIS from the Study information.

5-11—RETRIEVE IMAGES

Retrieve Images allows users to download all image files in a selected study from an archive server to their computer. Retrieve Images operates in the background, allowing users to perform other tasks while waiting. This


option is most useful for situations where there is a slow connection between the archive server and the local Diagnostic Viewer installation.

Click **Study** in the Study Browser and select **Retrieve Images**.

Note: Users can customize a hotkey or add a button to the Study Browser's Custom Toolbar for the Retrieve Images option.

5-12—PATIENT CD/DVD/USB


Users with a CD/DVD burner can activate the **Patient CD/DVD/USB** option to burn a patient's file to a CD, DVD (from the Managed Viewer or Web Viewer), or copy to a USB device. Users can continue employing the Diagnostic Viewer while a Patient CD/DVD/USB burns in the background. Alternatively, images can be shared using Novarad's CryptoChart Image Sharing solution.



CREATING PATIENT DISK

1. Select the desired study or series to burn to disk.
2. Click **Study** in the Study Browser and select **Patient CD/DVD/USB**.


The Patient CD/DVD/USB dialog box opens.



RIS Information

Accession:
Patient:
Complaint:
Referring:
Work: Mobile:
Reason:
Order Notes:
Procedure: FOREARM LT

3. Select the desired Studies, Series, and Compression options for the disk. Users can choose to burn just the selected study, or all studies for the patient.




Studies to Include on the Medium

Selected studies only

All studies for the selected patients

Users can select to burn all series for each study, or only select series in the study/studies.

Select series for each study is helpful for sites that have scanned in documents such as screening forms and other paperwork that should not be burned to the Patient disk. Before burning the disk, users see a list of each series in the study and can select each series to be included or excluded.




Choose Series to Burn to CD/DVD

Patient, Anonymous - CR 12/21/2009 2:09:17 PM

Series ...	Description	Has Critic...	Image Co...
1	PA		1

OK Cancel

By default, diagnostic reports are burned to the Patient disk. The **Include Diagnostic Reports** checkbox is selected by default each time users burn a Patient disk.




Include Diagnostic Reports

Compression: None

Note: Not all vendors support compressed images

OK Cancel

If a site is currently merging or fixing patients and studies, they may exclude diagnostic reports from the Patient disk by deselecting this checkbox. Sites may also choose to exclude reports to not give their patients a copy before meeting with their primary physician.



Users can select compression options for the Patient disk. Novarad recommends not using any compression for patient disks whenever possible because most vendors cannot support compressed images. Unless burning multiple or very large studies, choose **None** under the Compression drop-down list.

4. After selecting the Patient disk options, click **OK** to save the change(s). Click **Cancel** to cancel out of creating a Patient disk.

In order to include the exam report on a Patient CD/DVD/USB, the user who creates the disk must have a corresponding role with permissions in the NovaPACS Admin Console.


Note: Users can add a Patient CD/DVD/USB button to the Study Browser's Custom Toolbar.

EXPORTING STUDIES TO A USB DEVICE

1. Select a study from a list in the NovaPACS Study Browser and click **Study**.
2. Select **Patient CD/DVD/USB**.

The Patient CD/DVD/USB window opens.

3. Select **USB/Drive** in the Select Storage Medium section.



After clicking the USB/Drive option, the Drive drop-down list is enabled and displays all available drives, including removable drives.

Note: The Drive drop-down list defaults to the first removable drive available.

4. Click the **Drive** drop-down list to select a different drive, if applicable.
5. Choose to export **Selected studies only** or **All studies for the selected patients**.
6. Then choose to either export **All series for each study** or **Select series for each study**. The latter is helpful for studies that have additional series such as scanned document series that shouldn't be exported.
7. Finally, decide whether to export the diagnostic report. Click **OK**.
8. The **Select Folder** window opens. Users can pinpoint where they would like to save the file.
9. After making the selection, click **Select Folder** to save the file.
10. To quickly open the study that was saved to the USB device, users can click **ViewImages.cmd**.

This loads the Novarad Patient CD Viewer application along with the study.

Name	Date modified	Type	Size
Cache	7/21/2015 11:35 A...	File folder	
DICOM	7/21/2015 11:36 A...	File folder	
ImageTemp	7/21/2015 11:38 A...	File folder	
Reports	7/21/2015 11:38 A...	File folder	
Viewer	7/21/2015 11:38 A...	File folder	
Autorun.inf	7/21/2015 11:34 A...	Setup Information	1 KB
DICOMDIR	7/21/2015 11:35 A...	File	52 KB
ViewImages.cmd	7/21/2015 11:34 A...	Windows Comma...	1 KB


There are a total of four additional folders saved to the selected USB device: Cache, DICOM, Image Temp, and Viewer.

5-13—ASSIGN PASSWORD

Assign Password allows users to assign a temporary Single-Study Access Password for a study, giving a person (i.e. a Referring Physician or patient) access to that study through the **Web Viewer**. The Single-Study Access Password gives the user access to the study's images and RIS reports.

Note: For more information on the Web Viewer, see **Chapter 1, section 1-3—Accessing the Diagnostic Viewer**.

1. Click the study on any list in the Study Browser.
2. Click **Study** in the Study Browser and select **Assign Password**.



The dialog box is titled "Assign Web Password to Selected Study". It contains the following text: "This study can be accessed through the NovaPacs Web Viewer using the following password until the date specified." Below this, there are fields for "Study:" (set to "Anonymous Patient"), "Acc:" (set to "Anonymous Patient"), "Web Password:" (with a text input field and a "new password" link), and "Expiration Date:" (with a date picker set to "Thursday, January 22, 2015"). At the bottom are "OK" and "Cancel" buttons.

3. The Assign Web Password to Selected Study dialog box appears. A randomly generated password is assigned automatically. Users can generate their own password by clicking the **new password** link.
4. Click the **Expiration Date** drop-down list to set the length of time the password is valid and click **OK**.

5-14—SCAN COMPLETE

Scan Complete allows users to change the patient's status to reflect a completed scan, which reflects on the Patient Workflow page and Scheduling page in Nova RIS.

1. Select a study from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser and select **Scan Complete**.

Note: Users can customize a hotkey for the Scan Complete option.

5-15—LAUNCH MAMMO VIEWER

Launch Mammo Viewer allows users to launch a mammography viewing tool used with NovaPACS for mammography viewing.

1. Launch the Diagnostic Viewer and open an appropriate study (i.e. an MG or breast US study).
2. Return to the Study Browser by clicking the **Study Browser** icon in the top right corner of the Image Viewer window.

3. Click **Study** in the Study Browser and select **Launch Mammo Viewer**. The configured mammography viewer launches in a new window.
4. When finished using the Mammo Viewer, users can return to the PACS by closing the window.

Note: A Mammo Viewer software must be installed and configured on the workstation to use the Launch Mammo Viewer option. If it is not configured, Launch Mammo Viewer is unavailable on the Study Menu. Users can customize a hotkey or add a custom toolbar button to the Image Viewer or Study Browser to access this option.

5-16—VIEW STUDY ON OTHER WORKSTATION

View Study on Other Workstation allows users to send studies from the Diagnostic Viewer to another computer. Contact Novarad Support for the required configuration for this option. Once configured, studies can be viewed and processed by another computer. This option is useful if users need to view studies on post-processing stations for MIP/MPR or 3D Reconstruction.

1. Select a study from a list in the Study Browser.
2. Click **Study** in the Study Browser and select **View Study on Other Workstation**.

5-17—LAUNCH ORTHOPEDIC TOOLS


Launch Orthopedic Tools allows use of the orthopedic template third party software with Novarad's Diagnostic Viewer. This optional software integration allows users to create a complete array of orthopedic pre- and post-surgical planning and mapping tools including (but not limited to) Cobb Angles, templates, overlaying, and surgical repairs.

1. Launch the Diagnostic Viewer and open an appropriate study.
2. Click the **Study Browser icon** in the top right corner of the Image Viewer window.
3. Click **Study** in the Study Browser and select **Launch Orthopedic Tools**. The configured orthopedic tool launches in a new window.
4. When finished using the orthopedic tools, users can return to the PACS by closing the window.

Note: The product launched when using the Launch Orthopedic Tools option is dependent upon which software option is configured in the PACS Admin Console. If there is not an orthopedic software tool configured, the Launch Orthopedic Tools option is unavailable on the Study menu. Users can customize a hotkey or add a custom toolbar button to the Image Viewer or Study Browser to access this option.


5-18—TEACHING FILE

Teaching File allows users to create anonymized versions of particular studies.



This option is only available under the Study menu if **Teaching Files** have been enabled. To enable Teaching Files, navigate to the Options page of the Settings menu. At the bottom of the General tab, select the box next to **Enable Teaching Files**.

CREATING TEACHING FILES



Once **Teaching Files** have been enabled, users can click a study in the Study Browser and select **Teaching File** from the Study menu to create an anonymized version of the selected study. The Teaching File Information dialog box opens, allowing users to record information about the study while maintaining the anonymity of the patient once the information has been saved.

Note: For reference, the patient's name and ID is located at the top of the dialog while the users initially fills out the dialog. Once the information has been saved, the study no longer contains the patient's personal identification, and the dialog (if reopened) does not reference the original study or patient.

Note: Users can also create Teaching Files from the Image Viewer. See section 10-31 for more information.

TEACHING FILE INFORMATION DIALOG

Users should fill out the Teaching File Information dialog box as prompted, and for convenience, users can easily move through fields by pressing the tab key. Follow these guidelines for filling out and managing the Teaching File Information dialog:

THE REQUIRED TAB

The first few fields under the Required tab allow users to select from a drop-down menu of prepopulated options. Users can click and type into any of the other available fields. Any field with a red asterisk next to it is a required field that users must fill out in order to save the new information.

The Study Description field is automatically populated, based on information from the original study. However, users can choose to edit this information by clicking in the field and typing.

Note: Study Description changes save on the user's machine but are not saved to the server.

The Findings and Notes fields allow users to add in a large amount of text or copy and paste directly from the original report.

THE ADVANCED TAB

The Advanced tab allows users to record additional information about a study, though none of the fields under this tab are required fields. Each of the open fields under the Advanced tab have a limit of 256 characters per field. =

ADDITIONAL OPTIONS AND FIELDS

On the right side of the dialog box, users can locate who created a Teaching File, as identified by their username as well as their first and last name.

After information in a Teaching File dialog is saved, if it is opened and edited later, the dialog displays the identification of the last user to modify information. This is displayed beneath the Created By field.

Also on the right side of the dialog is the Teaching File Patient Information box, which contains three fields:

- **New Patient ID:** This field displays the Teaching File Patient ID. A new ID is generated for each new Teaching File, and each ID begins with the letters TF, indicating Teaching File. The New Patient ID field is not editable and is not related to the original study.

- **Last Name:** By default, all new Teaching File patients are assigned "Patient" as their last name. Users can assign a new last name by clicking inside the field and typing. Users can then search for the saved last name to locate the file later.

first name by clicking inside the field and typing. Users can then search for the saved first name to locate the file later.

- **First Name:** By default, all new Teaching File patients are assigned "Anonymous" as their first name. Users can assign a new

If any reports are available for a study, the user can select **Get Reports from Study** beneath the Teaching File Patient Information box.

users cannot select it.

If there are no available reports for a study, this button is grayed out, and

To save information in the Teaching File Information dialog, select **OK**. To cancel changes made, select **Cancel**.

Users should note the warning on the bottom of the dialog which informs users that they are primarily responsible for ensuring that no patient identification is located anywhere on the images or within the available fields in the Teaching File Information dialog.


An anonymized copy of this study will be stored to the server. You are responsible to ensure no patient identifying information is present on the images themselves or in the entered fields.

Note: Scanned images of reports or text within a study are **not** anonymized when a **Teaching File** is created. Users should plan accordingly.

ACCESSING SAVED TEACHING FILES

To access saved **Teaching Files**, filter for **Teaching File** studies in the Study Browser. First, enable Teaching File filters from the View menu.

To reopen or edit a Teaching File Information dialog, navigate to the Study Browser, search for the file of interest, click the study, and go to **Menu** → **Study** → **Teaching File**. The dialog opens, allowing users to edit information and save changes.



Users can also view Teaching File patient information from within the Study Browser without pulling up the study dialog. Locate the study and make sure it has been selected. Next, click the **Report** tab of the Series list in the Study Browser. The information displayed under this tab is drawn directly from the Teaching File

Information dialog.


5-19—STUDY OF INTEREST

Study of Interest allows users to mark studies for future reference. This allows users to easily search for and locate interesting studies.

1. Select a study from a list in the Study Browser.
2. Click **Study** in the Study Browser and select **Study of Interest**.
3. A message appears notifying users that the study they selected has been marked as a study of interest. Select **OK** or the red **X** in the upper right corner to exit the box.

Studies are then marked as studies of interest, and users can search for studies of interest by enabling the **Is Study of Interest** filter in the filter's bar.

5-20—COLLABORATION



Collaboration allows multiple Radiologists to work simultaneously from the same worklist without having duplicate readings. The Collaboration option also allows users to participate in peer-to-peer reviews of patient studies.

Studies just entering the NovaPACS system have a status of new or imported, depending on the origin of the study. Studies with a new status are coming from the modalities into the NovaPACS system. If users have an external source for studies, such as an long-term archive (LTA) or other DICOM source, studies imported and cataloged into the NovaPACS system have a status of imported.

The Collaboration option enables users to mark various statuses such as Verified, Checked Out, and Reviewed to manage studies in collaborative environments.

1. Select a study from a list in the Study Browser.
2. Click **Study** in the Study Browser and select **Collaboration**. A submenu drops down with collaboration options.

CHECKOUT

Checkout prevents multiple users from reviewing the same study simultaneously.

1. Select a study from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser. Select **Collaboration** and then **Checkout**.

The study now displays **Checked Out** in the Status column of the **Filtered Studies** list. The Current User column reflects the username of the person who checked the study out (if enabled as a display column).

The Checked-Out status for the study updates on all collaborating users' lists. Other users may view a Checked-Out study but cannot mark it Reviewed or use the Checkout function.

Alternate ways to use the Checkout option:

1. Users can set the Study Browser to automatically checkout a study when opened. Click **Settings** in the Study Browser. Click **Options**. Click the checkmark next to Auto-checkout studies when they are opened.
2. Users can customize a hotkey for the Checkout option.
3. Users can add a Checkout button to the Study Browser's Custom Toolbar.

UNDO CHECKOUT

Undo Checkout allows users to undo the checkout process for a study. This option is helpful when the user checks out the wrong study or is unable to complete the review of a study after checkout. The study or series is restored to the status it had before checkout (i.e. New or Verified).

The **Undo Checkout** option cannot restore the status of a study if it has already been marked Reviewed. Users can only change a study with a Reviewed status using **Reset Status to New** (see section **Reset Status to New**).

1. Select the checked-out study from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser, select **Collaboration**, and then select **Undo Checkout**.

Note: There are other ways to use the Undo Checkout option. Users can set a hotkey or add a button to the Study Browser's Custom Toolbar.

VERIFIED

Verified allows users to mark studies as verified, meaning that all attached documents, images, etc., are checked and ready for review.

1. Select the study to be verified from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser, select **Collaboration**, and then select **Verified**.

Note: There are other ways to use the Verified option. Users can set a hotkey or add a button to the Study Browser's Custom Toolbar.

The study has a Verified status. The Current User column reflects the username of the person who verified the study. The Verified status updates on all collaborating users' lists.


REVIEWED

Reviewed allows users to mark a study as reviewed.

1. Select the study that has been reviewed.
2. Click **Study** in the Study Browser, select **Collaboration**, and then select **Reviewed**.

Note: There are other ways to use the Reviewed option. Users can set a hotkey or add a button to the Study Browser's Custom Toolbar.

The study now has a Reviewed status. The Current User column reflects the username of the person who reviewed the study. The Reviewed status for the study updates on all collaborating users' lists.



Note: The status of a study resets to its previous New or Verified status if the study is closed out without marking it Reviewed, and if the following options are marked by navigating to **Settings** → **Options** → **General** tab:

- Auto-checkout studies when they are opened
- Auto-review studies when they are closed
- Prompt before marking a study reviewed

RESET STATUS TO NEW

Reset Status to New allows users to restore the status of a reviewed study back to New, so it can be reviewed again.


1. Select the reviewed study from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser, select **Collaboration**, and then select **Reset Status to New**.

Note: There are other ways to use the Reset Status to New option. Users can set a hotkey or add a button to the Study Browser's Custom Toolbar.

The study's status changes from Reviewed to New in the Status column of the Filtered Studies list. The username no longer appears in the Current User column.

Note: The Reset Status to New option must be used within six hours of a study being marked Reviewed. To change the status of a study after six hours, contact an administrator. Only administrators can reset studies older than six hours.

PEER REVIEW




Peer Review allows radiologists to evaluate each other's reviewed studies, give second opinions, and provide quality assurance. Any radiologist involved in the review process can view a completed Peer Review.

1. Select a reviewed study from the **Filtered Studies** list in the Study Browser.
2. Click **Study** in the Study Browser, select **Collaboration**, and then select **Peer Review**. The Peer Review dialog box appears.
3. Select an option under *Position*.
4. Enter any comments needed in the Comment field box and click **OK**.

Note: Peer Review reports are available to NovaPACS Administrators in the Admin Console. The Administrator can choose to report an average rate of missed diagnoses over a given time period.

UTILIZATION REVIEW

Utilization Review allows radiologists to review studies and provide quality assurance.



1. Select a reviewed study from the Filtered Studies list in the Study Browser.
2. Click **Study** in the Study Browser. Hover the cursor over **Collaboration** and select **Utilization Review**.

The Utilization Information for Reviewed Study dialog box appears.

3. Select the **Position** drop-down list and select the position that applies.
4. Enter any comments needed in the Comment field box and click **OK**.

Note: Administrators can pull a Study Utilization report on the Reports tab of the NovaPACS Admin Console to see all Utilization Review and study analysis information.

Select Report: Study Utilization

Filters


Date: < all > Reviewed last name: Assessment: < all > Search

Facility: < all > Reviewed first name:

Accession: Modality: < all >

Modality	Reviewed	Assessment	Comment	Updated	Study	Facility
Positive						
MR		Positive	This is a sample Utilization Review...	8/26/2015 1:44:08 PM		
Negative - Wrong Exam						
MR		Negative - Wrong Exam		8/21/2015 1:40:51 PM		

CHAPTER 6—USING THE STUDY BROWSER VIEW MENU




The **View** menu contains options for how information is displayed in the **Study Browser** window, as well as settings for the **Hanging Protocols** in the **Image Viewer**. The **View** menu is located on the main toolbar in the **Study Browser**.

6-1—WORKLIST TAB


Worklist Tab displays the **Assigned Worklist** containing studies that have been added by the user to view later or that have been assigned to the user by another user. Large radiologist groups use this function to improve workflow.

1. Click **View** on the **Study Browser**'s main toolbar, and select **Worklist Tab**.
2. Press the default hotkey **F2**.



The **Assigned Worklist** appears as a tab next to the **Search** worklist on top of the **Filtered Studies** list. The **Assigned Worklist** looks and functions like the **Search** worklist, except there are fewer search filters to use on the **Assigned Worklist** because studies are already filtered by assignment.

6-2—MANAGE HANGING PROTOCOLS



Manage Hanging Protocols are the saved layouts associated with a specific **Hanging Protocol** for displaying study images in the **Image Viewer**.


To activate the **Manage Hanging Protocols** section in the **Study Browser** window, click **Settings** on the **Study Browser**'s main toolbar and select **Manage Hanging Protocols**.

If activated, the Hanging Protocol Layouts section is located at the bottom of the Study Browser window and displays thumbnail images of the layouts defined for the current study. Users can toggle the Hanging Protocol Layouts using the default **hotkey F4**.

The Hanging Protocol Layouts section only displays content in the Study Browser if a study is open in the Image Viewer.

The title of this section also lists the current Hanging Protocol. This section appears on all tabs in the Study Browser (i.e., the Search tab, Assigned Worklist tab, and Advanced Worklist tabs).

6-3—PATIENT HISTORY




The **Patient History** list displays all completed studies for the selected patient, including all time frames and all modalities. The Patient History list is part of the default view in the Study Browser. It is below the Filtered Studies list and to the left of the Series list in the Study Browser window.

Note: Study Browser lists can be dragged and dropped to new locations, as well as toggled on/off at any time.

Users can choose to have the Patient History list not available by clicking **View** in the Study Browser and selecting **Patient History**.

6-4—SERIES



The **Series** list is part of the default view in the Study Browser. The Series list displays all the series that are part of the selected study. It is below the Filtered Studies list and to the right of the Patient History list in the Study Browser window.

Users can choose to have the Series list unavailable by clicking **View** on the Study Browser's main toolbar and selecting **Series**.

6-5—REFRESH LISTS

The lists in the Study Browser window can be refreshed to show the most up-to-date information using **Refresh Lists**.


1. Click **View** in the Study Browser and select **Refresh Lists**.
2. Press the default hotkey **F5**.

Note: To enable the Study Browser lists to automatically update when idle, click **Settings** in the Study Browser and select **Options**. Mark the check box next to Refresh Study Browser automatically when idle. The lists automatically refresh after the user has been idle for one minute.


6-6—RESET FILTERS

Users can reset filters when searching the Filtered Studies list using **Reset Filters**.

There are multiple ways to use the Reset Filters option:




1. Click **View** in the Study Browser and select **Reset Filters**.



2. Click **Reset Filters** on the Study Browser's Custom Toolbar (if set up).

3. Press the default hotkey **CTRL + Shift + F5**.



6-7—DOCKING TABS

Users can toggle on/off the docking tabs in the NovaPACS Study Browser. This may be helpful to conserve space in the Study Browser.

1. Click **View** in the Study Browser and select **Docking Tabs**.

The tab names no longer display in the Study Browser.


2. Toggle the tab names back on by reselecting **Docking Tabs** in the **View** menu.

3. If tabs are part of a tab group, only the selected or foremost tab displays when **Docking Tabs** is deselected.

If a tab needs to be visible all the time, drag it as its own docking section (top, bottom, left, or right of a docking section) before deselecting the **Docking Tabs** option.

Note: This feature is enabled by default. Once deselected, the **Filtered Studies**, **Patient History**, **Series**, **Report**, **Images**, **Preliminary Impressions**, **Preliminary Impressions Activity**, and **Hanging Protocols** Layouts tabs are no longer visible.

6-8—RADIOLOGIST PRODUCTIVITY DISPLAY



The Radiologist Productivity Display list shows how many studies have been Assigned to the current user and how many have been Reviewed, by each modality. This helps users keep track of their productivity.

The Radiologist Productivity Display list automatically refreshes every five minutes. Users can manually refresh their list at any time by clicking the **Refresh icon** in the top, right-hand corner of the list.


1. To enable the Radiologist Productivity Display, click **View** in the Study Browser and select **Radiologist Productivity Display**.

By default, the Radiologist Productivity Display list appears to the right of the Series list. However, this list can be dragged and dropped into any location as with all lists in the Study Browser.

2. To manually refresh the list, click the **Refresh icon** in the top, right-hand corner of this list.

Note: The Radiologist Productivity Display list automatically updates every five minutes, regardless of whether users manually refresh the list.


CHAPTER 7—USING THE STUDY BROWSER SETTINGS MENU



The **Settings** menu contains options to arrange and setup the Diagnostic Viewer's features, tools, and monitors. The Settings menu is on the Study Browser's main toolbar.

7-1—OPTIONS

The **Options** page in the **Settings** menu allows users to customize options for the Study Browser and the Image Viewer.



Click **Settings** in the Study Browser and select **Options**. The **Options** page opens.

Note: The **Options** page can also be set as a button on the Study Browser's Custom Toolbar.

GENERAL TAB

The **General Tab** on the **Options** page allows users to select general options for the Study Browser and Image Viewer.

1. Smart image layout

This helps users to use the space more efficiently in the series window. When selected, this layout causes the **Image Viewer** to place images into any unused space of a series window.

2. Auto-show study browser

This hides the Study Browser window from view when the cursor is moved off of the window. The Study Browser window reappears when the cursor moves to the far left or right edge of the **Image Viewer** window.

Note: If the Auto-show Study Browser option is not selected, the Study Browser window does not automatically appear when the cursor is moved to the far left or right edge of the Image Viewer window.

3. Auto-subtract multiframe angios

Some modalities create multiframe Angio images. The subtract feature lets users take away a frame from the entire multi-frame image and only show the changes between the other frames. **Auto-subtract** performs the subtracting automatically when the multi-frame image is opened.

Note: The subtract option is available in the Image Viewer to Cardiologists only.

4. **Auto-checkout studies when they are opened**

This sets the status of a study to *Checked Out* when it opens in the Image Viewer.

5. **Auto-review studies when they are closed**

This is a sub-option available when the **Auto-checkout studies when they are opened** option is selected. This automatically sets the status of a study to *Reviewed* when it is closed.

6. **Prompt before marking a study reviewed**

This is a sub-option available when the **Auto-checkout studies when they are opened** and the **Auto-review studies when they are closed** options are both selected. This prompts users to confirm that the study can be marked as *Reviewed*.

7. **Open a new workspace for each patient**

This allows users to open multiple studies in separate workspaces.

8. **Display older related studies when opening a study**

This opens previous studies along with the selected study in the Image Viewer.

Note: The Older Related Studies that open are determined by the Anatomical Area field of the image.

9. **Prompt for related studies when needed**

This is a sub-option available when the **Display older related studies when opening a study** option is selected. This prompts users when the Diagnostic Viewer cannot determine whether an older study is related or not, based on discrepancies of the Anatomical Area field of the image.

10. **Show patient history when opening a study**

This displays a list of other studies available for a patient when their study is opened in the Image Viewer.

11. **Show patient history when opening with a barcode**

This displays a list of other studies available for the patient when a study is opened with a barcode scanner.

12. **Show study notes when opening a study**

This opens a dialog box when a study is opened in the Image Viewer.

13. **Show other users' annotations**

This shows all annotations created for a study.

14. **Reset filters when opening the study browser**

This resets the search filters back to default settings every time the Study Browser is opened.

15. **Add study split when opening a second study for a patient**

This automatically splits the **Image Viewer** window to display multiple studies. If this option is not selected, the second study opened appears as a tab at the top of the Image Viewer window.

16. **Auto-sharpen XAs**

This automatically applies edge enhancement to XA images that may need more detail.

17. Auto-local contrast enhance CRs

This applies filters to CR images to enhance details in the image(s). Selecting this option applies the Local Contrast Enhancement to all CR images when opened in the Image Viewer.

18. Use multiple windows to display RIS reports

This allows users to view more than one exam report at a time. If this option is not selected, the default setting causes a second opened report to replace the first one.

19. Show only the most recent study for each patient in filtered studies

This shows only the most recent study for each patient in the Filtered Studies list in the Study Browser window. This is the default setting for users.

20. Always sort filtered studies by priority

This sorts studies by highest priority in the Filtered Studies list in the Study Browser window.

21. Display separate window when creating a key image series

*This automatically shows the Key Image Series in the **Image Viewer** window when a Key Image is added. If this option is not selected, a Key Image Series is still created, but does not automatically display in the Image Viewer window.*

22. Reset cursor when closing a study

This resets the mouse cursor back to the default settings if a study is closed while annotating. If this option is not selected, the mouse cursor retains the same settings that it was using when the previous study was closed.

23. Auto-sync hanging protocol image layout options

There are a few options in this section:

a. Always match layout

This displays all series windows using the same image layout, regardless of plane, when one series is adjusted.

b. Only match layout if series are coplanar

This synchronizes image layouts for all coplanar series.

c. Never match layout

If this option is selected, any adjustments to the layout need to be made to each series individually.

24. Mouse control options

There are a few options in this section:

a. Mouse scroll sensitivity

Increased sensitivity makes the cursor faster, while a decreased sensitivity results in a slower cursor.

b. Minimum mouse drag distance

This determines how long of a distance the mouse must be moved to respond. A higher setting increases the distance the mouse must move to respond. Users can select a setting from 1-100.

c. Mouse wheel sensitivity

This determines how much the mouse wheel needs to be moved before it responds. Increased sensitivity makes the mouse wheel respond faster, while a decreased sensitivity results in a slower response.

25. Prompt before marking a study reviewed if it has not been completely viewed

This prompts users with a dialog box stating that all of the images in the study have not been viewed. Users must confirm whether they want the study to be marked Reviewed.

26. Prompt before saving annotations

This prompts users to confirm whether they want to save their recent annotations to the study permanently.

27. Prompt before closing the application

This prompts users to confirm whether they want to close the entire Diagnostic Viewer when exiting.

28. Prompt before clearing the worklist

This prompts users to confirm whether they want to clear the Worklist.

29. Prompt to allow cross-localization of different frames of reference

This prompts users to confirm that they want to set up a cross localization between images with different frames of reference.

30. Display spine label cross references

This displays spine label references at the top of intersecting images.

31. Refresh study browser when idle

This automatically refreshes the Study Browser lists after the user has been idle for one minute.

Note: Contact Novarad Support to change the length of the idle time before auto-refresh begins in the Study Browser.

32. Show RIS information when opening a study

This brings up a RIS Information page in the **Image Viewer** when a study is opened.

33. Fade Away after This Many Seconds

This is a sub-option available when the **Show RIS information when opening a study** option is selected. This lets the user select the duration the RIS Information page is visible when the study is opened. Users can select from 0-100 seconds. The default is 60 seconds.

34. Magnify glass size

This allows users to adjust the size of the magnifying glass tool for the **Image Viewer**. Autosize is the default setting. When using Autosize, the magnifying glass automatically adjusts according to the available series window and study space. The glass size also grows or shrinks as the magnification increases or decreases.

35. Enable teaching files

This allows users to access and create Teaching Files in the Study Browser and Image Viewer.

36. Draw arrow annotations head to tail

This allows users to change the direction that arrows are drawn when annotating images. When selected, an arrow annotation is drawn with the pointed end of the arrow positioned where the user first clicks and the tail end of the arrow positioned where the user clicks a second time. This feature is turned off by default.

37. Action when closing primary study

The options beneath this section allow users to set preferences for how NovaPACS deals with comparison studies after the user closes a primary study.

a. Always leave comparisons open

This option leaves comparison studies open after primary studies have been closed.

b. Always close comparisons

This option automatically closes comparison studies when primary studies are closed.

c. Prompt for primary from remaining unread comparisons

When users select this option, the system sorts through comparison studies and prompts users to promote a comparison study to a primary study. Users are prompted only if the comparison study has not been reviewed.

d. Auto-promote remaining unread comparisons


This option auto-promotes an unread open comparison study to a primary study without prompting the user.

38. Radius Menu

This section allows users to enable or disable the Radius Menu, which is an alternative to the standard menu that appears when users right-click in the viewer.

- a. Pallet Mode:** *This is a consistent set of options that will stay the same for every study and modality.*
- b. Context Mode:** *The context mode contains ten options that remain consistent across all studies. The outer rings of the Radius Menu in the context menu contain options specific to the study or modality the user is working with.*
- c. Disabled:** *Disables the Radius Menu.*

To enable either of the modes, select the desired mode and then click OK.



KEYBOARD SHORTCUTS TAB (HOTKEYS)

The most commonly used features in the Study Browser and **Image Viewer** come with default hotkey settings. Hotkey settings for menu functions are listed next to function names in the menu drop-down lists. However, these default settings can be changed and customized.

To customize Keyboard Shortcuts, click **Settings** in the Study Browser, and select **Options**.

In the Options dialog box, select the **Keyboard Shortcuts** tab. The page has windows for both the Study Browser menu items and the **Image Viewer** menu items.

Assign a Keyboard Shortcut to the desired menu function by expanding the menu options in the Assign shortcuts to the Study Browser menu items window.


Select a Keyboard Shortcut option from the drop-down list. Click **Clear** to remove the current setting and click **OK** to select the preferred function.

AssignShortcut to the desired menu function by expanding the menu options in the Assign shortcuts to the Image Window menu items window. Type a keystroke or keystroke combination in the Image Window Keystroke field (any key combination can be used). Click **Clear** to remove the current setting and click **OK** to confirm the assignment.

Note: The Image Window settings do not allow the same Keyboard Shortcut to be assigned for multiple functions and prompts users if a combination is already in use. However, a Keyboard Shortcut assigned to a Study Browser function can be used for Image Viewer functions also.

MOUSE SETTINGS TAB

Most functions in the **Image Viewer** can be performed with the mouse instead of the keyboard for radiologists' convenience. The most commonly used features in the Image Viewer come with default mouse settings. These default settings can be changed and customized.




1. Viewing Mouse Settings

Users can change or create mouse settings using the Options menu. To

view the mouse settings, click **Settings** in the Study Browser, select **Options**, and click the **Mouse Settings** tab.

There are three categories of mouse settings that users can customize:

- **Assign click events to actions**—lists the functions available for control by mouse click.



- **Assign drag events to actions**—lists the functions available for control by dragging the mouse.
- **Assign wheel events to actions**—lists the functions available for control by the mouse wheel. Users can set the mouse wheel to scroll through the images slice by slice by assigning the UP or DOWN wheel events to the Previous or Next Image actions.

2. Assigning Functions to Mouse Settings


To assign or change a Mouse Setting, click one of the three drop-down lists (**Mouse Click Event**, **Mouse Drag Event**, or **Mouse Wheel Event**).

Select the preferred mouse action from the corresponding drop-down list


and click **OK**.

MOUSE 3D SETTINGS


The Mouse MPR/3D Settings allow users to configure what mouse functions do in NovaPACS MPR/3D Viewer. To access these settings, click the **Mouse MPR/3D Settings** tab in the Options menu.




Select the preferred View from the **View** drop-down



menu.



Click to highlight the desired action.



Select the preferred configurations for buttons, modifiers, and mouse action.

Once the desired settings have been configured, click **OK**.


Note: The **Invert** option is available. To enable this feature, select **Invert** below Mouse Action.

IMAGE CONTEXT MENU TAB

The **Image Context Menu** allows users to add or remove options in the various modality Image Context menus.

Click **Settings** in the Study Browser and select **Options**. Click the **Image Context Menu** tab.

The **Image Context Menu** tab has two windows—an Available items window and a Current toolbar items window.



To add Image Context Functions to a specific modality's Image Context menu, click the **Modality** selection field at the top of the Image Context Menu tab page, and select a modality.

Note: If users choose all modalities, the option they modify shows up in all Image Context Menus when they right-click. However, if users choose a specific modality, the option they modify shows up in the Image Context Menu for that modality only.

Select a function to change or add from the Available items list. Click **Add** to add the function to the Current toolbar items list and click **OK**. The function is available in the Image Context Menu in the Image Viewer when the user right-clicks.

To remove Image Context Functions from a specific modality's Image Context Menu, click the **Modality** selection field at the top of the Image Context Menu tab page and select a modality.


Select the unwanted function from the Current toolbar items list. Click **Remove** to delete the function from the Image Context Menu, then click **OK**.

RE-ORDERING THE IMAGE CONTEXT FUNCTIONS

Users can customize the order of the Image Context Menu items. Items can be reset to original settings and moved to a different location or hierarchy on the list as well.

To do this, select an item on the Current toolbar items window list and click **Reset**, **Move Up**, **Move In**, **Move Out**, or **Move Down**. Click **OK**.

AUTO W/L TAB



Users can view and manage Auto Window/Level Adjustments from the Auto W/L tab of the Options page accessed from the Study Browser's Settings menu.


The Options page Adjustments box lists all Auto Window/Level Adjustments saved for the current user.

The Filters box lists all filter settings applied across a modality.

Users can add, edit, and remove Adjustments and Filters from this tab to ensure that images always display according to their preferences.

Note: There is only one set of Add, Edit, and Remove buttons. These buttons are used to manage both Adjustments and Filters.

ADDING ADJUSTMENTS



To add an adjustment, click inside the **Adjustments** box and select **Add** at the bottom of the page. The Window Level Adjustment dialog box opens.

The Window Level Adjustment dialog box allows users to manually specify the modality and image type the Auto Window/Level Adjustment applies to as well as the desired Window/Level values.


In the Modality field, specify the modality to be associated with the new adjustment. In the Low Window, High Window, Low Level, and High-Level fields, specify the range of intensity values for the images to which the new adjustment applies.

In the New Window and New Level fields, specify the values to be automatically applied to the specified images.

Note: When created directly from the Image Viewer, these values are automatically calculated and populated based on the currently selected image and applied Window/Level values.

Click **OK** when finished. The new Window/Level Adjustment is saved and appears in the Adjustments list. The new Window/Level Adjustment automatically applies to the appropriate images the next time they are opened.

In order to save a new adjustment, all fields must be completed. If fields are not completed appropriately, warning dialog boxes prompt the user to correct the error before saving.



EDITING ADJUSTMENTS

Users can edit Window/Level Adjustment by selecting the item they wish to edit from the Adjustments list and clicking **Edit**.

The Edit Parameters dialog box appears, allowing users to specify the images that the adjustment applies to as well as the Window/Level values to be applied.

Note: Users may not edit the modality of an existing Auto Window/Level Adjustment.

In the Low Window, High Window, Low Level, and High-Level fields, users may edit the range of intensity values for the images to which the new adjustment applies. In the New Window and New Level fields, users may choose the values to be automatically applied to the specified images.

Click **OK** when finished. Changes to the Window/Level Adjustment are saved and appear in the Adjustments list. The new changes automatically apply to the appropriate images the next time they are opened.

In order to save changes to an auto Window/Level Adjustment, all fields must be completed. If fields are not completed appropriately, warning dialog boxes prompt the user to correct the error before saving.

REMOVING ADJUSTMENTS

Users can delete Window/Level Adjustment by selecting the item they wish to delete from the Adjustments list and clicking **Remove**. The Window/Level Adjustment is deleted immediately and no longer automatically applies to images.

ADDING FILTERS

To add a filter, click inside the Filters box and select **Add** at the bottom of the page.

The Window Level Filter dialog box opens.

The Window Level Filter dialog box allows users to manually manage filters that apply across an entire modality. Only one filter may be created per modality.

In the Modality field, specify a modality for the new filter.

In the Low Window field, enter the minimum low window value for the filter.

In the High Window field, enter the maximum high window value for the filter.

In the Low-Level field, enter the minimum low-level value for the filter.

In the High-Level field, enter the maximum high-level value for the filter.

Enter the new window value in the New Window field.

Enter the new level value in the New Level field.


Click **OK** to confirm the changes. Click **Cancel** to exit without adding the filter.

In order to save a new filter, all fields must be completed. If fields are not completed appropriately, warning dialog boxes prompt the user to correct the error before saving.

EDITING FILTERS

Users can edit existing filters by selecting the item they wish to edit from the Filters list and clicking **Edit** or double-clicking a listed filter.

The Edit Parameters dialog box appears, allowing users to edit filter values.



Note: Users may not edit the modality of an existing filter.

Users can edit the Minimum and Maximum Pixel value fields if needed.

By default, the NovaPACS Diagnostic Viewer uses the window/level values supplied by the modality that are in the DICOM file. Users can decide whether to apply the DICOM window/level values or have the Viewer calculate optimum window/level values.

Users can edit this by clicking the **Use DICOM Values** drop-down list.

- **Always**—This always uses the window/level values supplied by the modality in the DICOM file for the selected modality type.
- **When Reasonable**—The NovaPACS Diagnostic Viewer calculates auto window/level values from the selected modality type and compares them to the values supplied in the DICOM file. Images display with the most optimal window/level values of the two.
- **Never**—This never uses the supplied DICOM window/level values for the selected modality type and the NovaPACS Diagnostic Viewer calculates optimal auto window/level values.

In addition, if users want the NovaPACS Diagnostic Viewer to always calculate auto window/level values for the selected modality type, they can select the **Always Use Calculated Values** checkbox—which takes precedence over the **Use DICOM Values** drop-down list settings.

This ignores the supplied DICOM auto window/level values and the NovaPACS Diagnostic Viewer calculates optimal auto window/level values.

If users decide to have the NovaPACS Diagnostic Viewer calculate auto window/levels instead of using the supplied DICOM auto window/level values, they can select **Calculate Each Image**.

This ensures that auto window/level values are calculated using values from each image in a series, instead of an average from the entire series. Some modality types have a high variety of values from image to image, and choosing this option ensures a more accurate calculation.

Additional settings on this page are best configured only by Novarad personnel. Contact Novarad Support for assistance in setting/adjusting auto window/level filter values.

Make any necessary changes and click **OK**. Changes to the filter are saved and appear in the Filters list. The new changes automatically apply to the appropriate Modality.

In order to save changes to a filter, all fields must be completed. If fields are not completed appropriately, warning dialog boxes prompt the user to correct the error before saving.

REMOVING FILTERS

Users can remove Filters by selecting the item they wish to delete from the Filters list and clicking **Remove**.

The Filter is removed immediately and no longer applies to images of the selected modality type.

LABELS TAB

Admins can change and customize the DICOM labels that appear on images in the Image Viewer. These labels provide useful patient, study, and series information. Font sizes and colors can be customized for the labels as well. Labels automatically adjust according to each monitor's resolution.

Note: The ability to customize or edit DICOM labels is a role-based permission and can be assigned by Admins.

VIEWING THE LABELS TAB



Click **Settings** in the Study Browser and select **Options**. Select the **Labels** tab.




The Labels page displays the default DICOM label settings. These default settings are used if a specific modality's labels are not customized. Select a modality from the **Modality** drop-down list to view its DICOM label settings. Users can create and modify their own custom labels from the Labels tab.

Note: The default color of a label is determined by the level of that label (i.e. the Patient level has a default color of light teal; the Study level has a default color of dark teal; the Series level has a default color of purple; and the Image level has a default color of gold).

Editing Label Font Size and Color




- a. Select a modality from the **Modality** drop-down list.
- b. Click **Size/Color** above the DICOM Fields section. The **Edit Default Label Properties** dialog box opens.
- c. Select a label from the **Select label level to edit** drop-down list. The **Size** and **Color** fields become available once a label is selected.
- d. Select a size for the font from the **Size** drop-down list. The default size is **Normal**.



- e. To change the font color, click the colored square next to the **Color** field. The **Color** dialog box opens allowing users to select from the **Basic Colors** options or create a custom color using the **Define Custom Colors** option.

To create a custom color, click **Define Custom Colors**. Use the cross-point tool and/or enter values in the **Hue**, **Sat**, **Lum**, **Red**, **Green**, and **Blue** custom boxes to define the desired color.




Click **Add to Custom Colors**. Select that custom color as the font color and click **OK** to save the changes.

ADDING LABELS

Admins can add a Label for a specific modality's image.

- a. Select a modality from the **Modality** drop-down list.
- b. Select the desired label from the **DICOM Fields** list.
- c. Drag the label onto the black image box on the page to add it to the labels for the modality's images.

Note: Other roles, such as Radiologist, can modify DICOM labels, when Admin enables permission to do so.




Note: Users cannot drag a label from one quadrant of the black image box to another. To move a label from one quadrant to another, users must first delete the image label by double-clicking it. Then, drag the label onto the appropriate quadrant. In addition, users cannot add duplicate labels in the same quadrant of the image box.

REMOVING LABELS

Double-click the label within the black image box to remove it.

CREATING LABELS


If users need a label that is not available in the DICOM Fields list, labels can be created and added to any specific modality's images.



Note: Users must know the specific group and element numbers located on the image to create a label. This is found using the Image Details option in the Image Viewer.

- Click **Create** located at the bottom of the DICOM Fields list. The Create DICOM Label dialog box appears.
- Enter the necessary information into the fields in the **Create DICOM Label** dialog box.

- Name**—Enter a name for the label, for example use “View” for FOV.




- Modality**—Choose a modality to make this DICOM label available for a specific modality.
- Display Label**—Enter a short, abbreviated title for the label (i.e. ACC# for Accession Number).
- DICOM Tag Fields**—Specify the DICOM group and element from which to pull information.
- Text Display Properties**—Use the **Size** drop-down list to select the text size. The default text size is Normal. Select the **Color** box to choose a basic color or create a custom color.

EDITING DICOM FIELDS

- Select a field from the **DICOM Fields** list and click **Edit**. The Edit DICOM Label dialog box opens.
- Make any needed changes to the information fields and click **OK** to save the changes. Click **Cancel** to exit without saving changes.

LEXICON TAB



A **Lexicon** value is one or more words commonly used to annotate images in the Image Viewer. Users can add phrases that are used frequently and have the option to select these from a list in the Image Viewer instead of manually typing them. An administrator can create global Lexicon values available to all users, and customized values can be created by each user available only to them.


VIEWING THE LEXICON TAB

Click **Settings** in the Study Browser and select **Options**. Click the **Lexicon** tab.

Adding Lexicon Values

1. On the Lexicon tab page, click **Add** at the bottom. The Add Lexicon Item dialog box opens.
2. Select a modality from the **Modality** drop-down list to add a Lexicon value.
3. Enter the desired text into the Lexicon Text box and click **OK** to save the new Lexicon value.

EDITING LEXICON VALUES



1. Select a Lexicon value from the Lexicon list.
2. Click **Edit** at the bottom of the page. The Edit Lexicon Item dialog box opens.
3. Make any needed changes to the Lexicon value and click **OK** to save.


DELETING LEXICON VALUES

1. Select a Lexicon value from the Lexicon list.
2. Click **Delete** at the bottom of the page to delete the value from the list.

PLUG-INS TAB

The NovaPACS Diagnostic Viewer has three **Plug-ins** options available:

- **Image Rendering Provider**—options are **DirectX**
- **Dictation Provider**—options are **Basic, File, Dragon, CMDLine, Fluency Direct, and PowerScribe360**



Note: Users can configure the PowerScribe Plugin settings by going to **Settings** → **Options** → **Plugins** → **Plugin Settings** and clicking on the **PowerScribe** tab.


- **CD Burning Provider**—options are **Primera, Legacy, and Basic**

Note: The Dictation Provider and CD Burning Provider options may require additional software.


To view and/or change the Plug-ins options, click **Settings** in the Study Browser, and select **Options**. Select the **Plug-ins** tab.

Use the drop-down lists next to the Plug-ins to select new options.

Note: In order to edit the Plug-ins options the user's credentials must be verified as the Windows Administrator.



DATA SOURCES TAB



Data Sources are the connections that allow the Diagnostic Viewer to communicate with the Database and Image Servers. Data Sources are the primary locations from which images are downloaded.

Note: Users cannot add, edit, or delete Data Sources for DICOM Pass Through servers.

VIEWING DATA SOURCES

Click **Settings** in the Study Browser and select **Options**. Select the **Data Sources** tab.

All the currently established Data Sources are listed on the Data Sources tab page. Users can change the Primary Archive and Image Cache options using the corresponding drop-down lists.

ADDING A DATA SOURCE

1. On the Data Sources page, click **Add** at the bottom. The Add Data Sources dialog box opens.
2. Enter the appropriate Friendly Name and Server Address information.
3. Select the NovaPACS or DICOM source type from the Source Type section.
4. If NovaPACS is selected as the source type, choose the appropriate Port and Compress Options. Users can select checkboxes to Stream Images, Automatically retrieve images queued for this machine (Client Requested Transfer), and Automatically login.

Note: If Automatically Login is selected, users are prompted to enter a username and password and verify a login for each data source.

5. If the DICOM source type is selected, enter the appropriate information for the Port and AE Title fields.

EDITING DATA SOURCES

1. Select the **Data Source** listing from the Data Sources list.
2. Click **Edit** at the bottom of the page. The Edit Data Source dialog box opens.
3. Make any needed changes to the Data Source information and click **OK** to confirm.


DELETING DATA SOURCES

1. Select the Data Source from the Data Sources list.
2. Click **Remove** at the bottom of the page. The Data Source is removed from the Diagnostic Viewer's Data Sources options.

TESTING DATA SOURCES

1. Select the Data Source from the Data Sources list.
2. Click **Test** at the bottom of the page.
3. The Diagnostic Viewer tests the Data Source. If the Data Source is incorrect, users receive a warning message stating that the Data Source failed, however, if it is functioning properly, users do not see any message.

7-2—W/L PRESETS




W/L Presets allows users to specify window/level preset values for each modality when using the Image Viewer. The **Image Viewer** sets default W/L values for each modality that correspond in number with the F1-F11 keys on the keyboard (i.e. F1 corresponds with the W/L setting set to one).

Note: Up to 12 W/L Presets may be specified; however, only 11 may be accessed from default hotkeys F1-F11. W/L Preset 12 can only be accessed from the menu or toolbar. The F12 keystroke is the default hotkey for Auto Window/Level Adjustments.

CUSTOMIZING THE W/L SETTINGS

Users can customize the W/L values and save them to user preferences.



1. Click **Settings** in the Study Browser and select **W/L Presets**. The Configure Window/Level Presets dialog box appears.
2. Select a modality from the **Modality** drop-down list.
3. Select the desired preset value from the **W/L Preset** drop-down list. The Preset Values fields populate with the values assigned to that W/L Preset.
4. Change the Menu/Toolbar Name by entering text in the field. The Menu/Toolbar Name appears on the W/L Toolbar in the Image Viewer.

5. The **Apply Current Window/Level** button allows users to apply the W/L settings used in the Image Viewer to the W/L Preset being customized. Adjust the window/level settings in the Image Viewer, then click **Apply Current Window/Level** to apply the settings to the W/L Preset.

CUSTOMIZING WIDTH AND CENTER FOR WINDOW/LEVEL FIELDS AND PRESETS

Users can customize the width and center values for the window/level fields and presets.

To customize the width and center window/level values, click **Settings** in the Study Browser and select **W/L Presets**. The Configure Window/Level Presets dialog box appears.

1. Click the up or down arrows next to the Window (Width) and/or Level (Center) fields to choose the values.
2. Manually enter the values in the Window (Width) and/or Level (Center) fields.
3. Click **Apply Current Window/Level** to apply the settings being used in the active study in the Image Viewer.
4. To set the W/L Preset in the Image Viewer, close the Configure Window/Level Presets dialog box and go to the Image Viewer. Make any window/level adjustments needed on the image and press **CTRL + F__** (the F key


corresponding to the changing preset, i.e. F6 for preset six). This sets the adjustments as the W/L Preset for that F key.

Note: To view the saved settings in the Configure Window/Level Presets dialog box, go back to the Study Browser, click **Settings**, and select **W/L Presets**. The Configure Window/Level Presets dialog opens. Select the changed modality and the changed preset number, and the new values populate in the Window (Width) and Level (Center) fields.

7-3—MANAGE HANGING PROTOCOLS

The **Hanging Protocol Layouts** are the saved layouts associated with a specific modality (and optionally, study description) for displaying study images and series in the Image Viewer.

Users can define multiple Hanging Protocols for each modality and define multiple Layouts for each Hanging Protocol.




Within each Layout, users can filter individual Image Boxes to display only certain types of images. This allows users to easily create, manage, and navigate between the specific views they prefer for individual study types.



Users can manage Hanging Protocols from the Study Browser menu by clicking **Settings** in the Study Browser and selecting **Manage Hanging Protocols**.

The Manage Hanging Protocols editor dialog box opens.

ADD NEW HANGING PROTOCOL



Click **Add New** at the bottom of the Study Types list. The Manage Hanging Protocol editor now allows the user to select options for the new Hanging Protocol.

Select Modality


First select a modality from the **Modality** drop-down list. This is the modality with which the new Hanging Protocol is associated. A modality must be selected in order to save a new Hanging Protocol. A Hanging Protocol can only be applied to a study

if they are both associated with the same modality.


MANAGE STUDY DESCRIPTIONS

A Hanging Protocol applies automatically to any study that matches all of its Study Descriptions. Adding Study Descriptions helps users distinguish what types of studies the Hanging Protocol was set up to display. Users may manually apply a Hanging Protocol to a study that does not match its study description, as long as both are associated with the same modality.


To add a Study Description, click **Add** to the right of the Study Descriptions list.




A small text box appears inside the Study Descriptions list. Type in an anatomical area to associate with the new Hanging Protocol. Repeat as many Study Descriptions as needed.



Note: If users edit/create a hanging protocol with a study open in the Image Viewer by clicking the gear icon, all applicable study and series descriptions pre-populate.



A unique Study Description (or combination of Study Descriptions) must be added in order to save a new Hanging Protocol for the same modality as an existing Hanging Protocol. Radiologists use different Study Descriptions to differentiate between study types. They use different Hanging Protocols to view those studies in a specific order.



To edit an existing Study Description, highlight the description and click **Edit**.


To delete an existing Study Description, highlight the description and click **Remove**.





SELECT A DISPLAY TOOL OPTION

Select an option from the **Display Tool** drop-down list to specify how to display the Hanging Protocol.

SELECT A SPINE LABELING DEFAULT



Select an option from the **Spine Labeling Defaults** drop-down list to specify the spine labeling configuration for the new Hanging Protocol.




Select **Edit Options** to customize the available spine labeling options. This allows users to add, edit, and

remove spine labeling options. It also allows users to customize font and font color for spine labeling options.


MANAGE LAYOUTS

Users can add, edit, remove, and duplicate layouts within a Hanging Protocol using the layouts section of the Manage Hanging Protocol editor.



Add Layout

Click **Create Layout** at the top of the layouts section. This adds a new, blank layout.




Edit Layout

Click **Edit Layout** at the top right of the layouts section or double click the thumbnail image of a layout. The Edit Layout window opens.

Editing Series Layout


There are three primary ways to edit the Series Layout.

1. Click the box with green gridlines, located above the image boxes.



Choose Series Layout box opens.


Choose how the series displays by hovering the cursor vertically and horizontally over the squares and clicking when the preferred layout is highlighted in light gray.



2. Click the appropriate **Split Window Tool** located in the top bar of an image box to edit the image layout within that image box or to split the series horizontally or vertically.
3. Click the red X at the top of an image box to delete that window from the series layout.

USING THE IMAGE LAYOUT SETTINGS TAB

Select the **Image Box**. This enables the **Image Layout Settings** tab and allows users to choose image layout options for the selected portion of the layout.



On the **Image Layout Settings** tab, users can edit the **Layout type**, **Image Layout**, and **Clip playback**.

The options for editing the **Layout type** are as follows:

The images in the series can be set to a **Tiled display**. This is the default view.

When the **Tiled** option is selected, users can do the following:

1. Use the **Choose Image Layout** box to select how images display within the series. Select how the images display by hovering the mouse vertically and horizontally over the squares and clicking when the preferred layout is highlighted in light gray.
2. Select **Zoom To Fit**, **Zoom To Fill**, **Custom**, or **No Zoom**.
3. Adjust the **Frame rate**, which controls the playback speed.
4. Select the **Auto-start** and **Use DICOM Frame Rate** checkboxes.

Note: Auto-start begins to play the clip as soon as the image loads. Use **DICOM Frame Rate** to play the image frame rate at the recorded rate. When not checked, the **Image Viewer** determines the image frame rate by the **Desired Frame Rate** count.

5. Select **Looping**, **Sweeping**, or **Stop** for the **Playback sequencing**.
 - a. **Looping** repeats the playback from the first image to the last and then begins again with the first image.

- b. **Sweeping** repeats the playback from the first image to the last and then from the last image to the first.
- c. **Stop** allows users to control image playback by manually advancing or reversing a single image with the forward and back buttons. After going to the final image, the series stops and is displayed in default Tiled layout.

Cine


The images in a series can be set to a Cine display. Cine is typically used for Cardiology images.

When Cine is selected, users can do as follows:

1. Set the Desired Frame Rate to control the playback speed.
2. Select the Auto-Start checkbox to begin playing the clip as soon as the image loads.
3. Select the Use DICOM Frame Rate checkbox

Note: This plays the image at the rate recorded if selected. If not selected, the Image Viewer determines the image frame rate by the Desired Frame Rate count.

4. Set Playback Sequencing to Looping, Sweeping, or Stop



The images in a series can be set to be processed by either **MPR** or **Image Fusion**. Select the **Processed** layout type, then select an option from the **Processed By** drop-down list.

If **Processed by MPR** is selected, the series automatically displays in a **Multi-Planar Reconstruction (MPR)** view. This view reconstructs multiple 2D images into a multi-planar 3D view. Novarad provides MIP/MPR functionality and support with NovaPACS.

Note: If the system has not been used with the Diagnostic Viewer's MIP/MPR support, an administrator must contact Technical Support to configure each workstation to access MIP/MPR.

If **Image Fusion** is selected in **Processed By**, two series from different modalities (typically CT and PET) fuse for a more detailed image.

When finished editing a layout, click **OK** to save the changes. Click **Cancel** to close the Edit Layout window without saving changes and return to the Manage Hanging Protocols editor.

1. Remove Layout

Select the layout's thumbnail image from the Manage Hanging Protocols editor and click **Remove layout**. The layout is deleted from the current Hanging Protocol.

2. Duplicate Layout

Select the layout's thumbnail image from the Manage Hanging Protocols editor and click **Duplicate layout**. The layout is duplicated in the current Hanging Protocol.

EDIT HANGING PROTOCOL

Click on an existing Hanging Protocol's description in the Study Types list. The Manage Hanging Protocols editor allows users to make changes to any options described in the previous sections.


REMOVE HANGING PROTOCOL



Select an existing Hanging Protocol's description and click **Remove** below the Study Types list to remove it.

When finished managing Hanging Protocols, click **Close** to exit the Manage Hanging Protocols editor and return to the Study Browser.

7-4—CONFIGURE MONITORS



Configure Monitors allows users to select monitors for viewing images. Any monitor properly connected to the user's computer appears in the **Configure Monitors** window as numbered rectangles.

Click **Settings** in the Study Browser and select **Configure Monitors**. Users can also press the default hotkey **CTRL + M**. The **Configure Monitors** dialog box opens.

Users can select a monitor or combination of monitors by clicking the rectangles. Selected monitors have a red border. Deselect monitors by clicking a selected rectangle again. Users can also select which kinds of images will be displayed on each monitor. The current options are **NovaPACS** and **NovaMG PRO**.




Note: On a three-monitor system, where the third monitor is used only for displaying the Study Browser, make sure the third monitor for the Study Browser is not selected for use with images.

Note: when configuring monitors, users cannot have a monitor between two configured monitors. The configured monitors must be contiguous.

Users can also determine whether the left monitor displays new studies or old studies by selecting the options at the bottom of the Configure Monitors dialog box. Click **OK** to save the monitor settings.

7-5—PREFERENCES




Preferences allows users to select settings for viewing and organizing images and information in the Diagnostic Viewer. Any changes the user makes to the default preferences settings are stored and automatically recalled each time that user logs into the Diagnostic Viewer.

Note: Novarad recommends that users periodically export a copy of their preferences to the desktop, a flash drive, or another data storage device.

Click **Settings** in the Study Browser and select the **Preferences** drop-down list.

RESET PREFERENCES


Reset Preferences allows users to reset user preferences to the site's default settings.



Note: The Diagnostic Viewer permanently overwrites any settings that have been changed in the Study Browser and the Image Viewer when the Reset Preferences option is selected. Users receive a warning before resetting.

Click **Settings** in the Study Browser and select the **Preferences** drop-down list. Select **Reset Preferences**.

The Reset Preferences dialog box opens.



Click the **Select the source profile for preferences** drop-down list and select the role to use for default preferences.


Once the role is selected, choose whether to reset the preferences to the **Basic, Normal, or Advanced** level. Administrators can customize and save these levels.

Click **OK** to confirm the reset.

IMPORT PREFERENCES

Import Preferences allows users to import preference settings from a saved location. This is helpful when sharing preference settings from one Diagnostic Viewer installation to another.

Click **Settings** in the Study Browser. Hover the cursor over **Preferences** and select **Import Preferences**.



Browse to find the desired preferences file and click **Open**.

Note: When exiting the Diagnostic Viewer, if imported preferences are not stored in the database, users receive a warning before the preferences are overwritten.

EXPORT PREFERENCES

Export Preferences allows users to save and export a copy of their user preferences in case they are ever lost or reset.

Click **Settings** in the Study Browser and hover the cursor over **Preferences**. Select **Export Preferences**.

The Export Preferences to File dialog box opens.

Select a destination for the file by clicking a listed location or browsing through the computer and file options.

Name the preferences file in the File Name field.

Note: Novarad recommends naming the user preferences with the username and the current date, to quickly locate them later.


Click **Save** to export a copy of the user preferences.

RESTORE FROM BACKUP

Restore from Backup allows users to restore their user preferences from preference files that have been backed up automatically on the NovaPACS server. This is helpful if users have had connection or server issues and their user preferences have been lost. Users can see up to the last ten user preference files that have been saved and can quickly restore their user preferences to a selected date's preferences.

Note: User preferences are only backed up automatically and available if the user has made any changes/customizations while logged in. If the user has not made any changes/customizations, that day's preferences are not available on the Restore from Backup list.

1. Click **Settings** in the Study Browser and hover the cursor over **Preferences**. Select **Restore from Backup**.




The Restore From Backup dialog box opens.

2. Click the **Most recent preference backups** drop-down list and select a date. The dates listed are the most recent dates when preference backups were saved to the server.

3. Click **Restore** to restore the user preferences to the selected date's backed-up user preference file.


EXPORT USER PREFERENCES

Export User Preferences enables administrators to export all of their users' preferences from their own user login. This allows admins to quickly export all users' preferences from a single location, and from a remote location if necessary.




1. Click **Settings** in the Study Browser and hover the cursor over **Preferences**. Select **Export User Preferences**.


The **Export User Preferences** window opens.



2. Select the user preference(s) to export. Users can hold down **CTRL** or **SHIFT** to select multiple user preferences.



3. Select the location for the exported preferences by clicking **Browse**.




The **Browse for Folder** window opens.

4. Browse for a location, select it, and click **OK**.
5. Click **Export** in the bottom right-hand corner of the window.

Users receive a notification when the user preferences are successfully exported.

Users can navigate to the save location to verify that the users' preferences have successfully exported.



SAVE AS DEFAULT PREFERENCES


Save as Default Preferences allows administrators to save default user preferences for all users. A new profile for CD/DVD/USB preferences can also be set for users burning Patient disks.

Note: Only an administrator can alter Default User Preferences. It is important to remember that all settings including buttons, hotkeys, and the general layout of the window are what the users can use. This includes CD/DVD/USB burning preferences.

1. Make all of the preference changes needed (i.e. hotkeys, toolbars, image context menus, mouse settings, etc.).
2. Click **Settings** in the Study Browser, and hover over **Preferences**. Select **Save as Default Preferences**.
The Set Default Preferences for Profile dialog box opens.
3. Select the user role to save preferences.
4. Select the level of preferences to save. The options are **Basic**, **Normal**, or **Advanced**.
5. Click **OK** to save the current user preferences for that role.

SAVING INDIVIDUAL PREFERENCES


Users can save specific preferences. Click **Settings** in the Study Browser and hover the cursor over **Preferences**.



Select **Import**, **Restore From Backup**, or **Reset to Default**. A separate window opens.

Click **Browse** to select a preference file.

To expand the **All Preferences** checkbox, click the arrow to its left. The window extends and reveals a list of individual preferences.




From the extended window, users can select which preferences to save by selecting or deselecting the checkboxes next to the specific preferences.

7-6— CUSTOMIZE TOOLBAR

Customize Toolbar allows users to edit the functions on the **Study Browser**'s Custom Toolbar, add or remove toolbar items, and reorder the menus.

The following sections describe what functions users can do with the Study Browser toolbars:

ADDING TOOLBAR ITEMS




1. Click **Settings** in the Study Browser and select **Customize Toolbar**. Users can also right-click on any toolbar in the Study Browser window and select **Customize**. A Custom Toolbars page opens.
2. Select the toolbar from the **Select Toolbar** drop-down list.
3. Select the desired feature from the Available Items list on the left side of the page.
4. Click **Add** to add the feature to the Current Toolbar Items list.
5. Click **Save** to save the change(s) and exit.

REMOVING TOOLBAR ITEMS

1. Click **Settings** in the Study Browser and select **Customize Toolbar**.
2. Select the toolbar from the **Select Toolbar** drop-down list.
3. Select the feature to remove from the Current Toolbar Items list.
4. Click **Remove** to remove the feature from the Current Toolbar items list.
5. Click **Save** to save the change(s) and exit.

REORDERING THE MENU



Users can customize the Study Browser toolbars by resetting the toolbar settings or reordering individual features.


1. Click **Settings** in the Study Browser and select **Customize Toolbar**.
2. Select the toolbar from the **Select Toolbar** drop-down list.
3. Select the feature from the Current toolbar items list.
4. Click **Reset**, **Move Up**, **Move In**, **Move Out**, or **Move Down** until the feature is in the desired location.
5. Click **Save** to save the change(s) and exit.

CUSTOMIZING TOOLBAR OPTIONS

1. Click **Settings** in the Study Browser and select **Customize Toolbar**.
2. Click the **Toolbar Options** tab at the top of the page.
3. Select an option from the **Button Layout** drop-down list (**Text Only**, **Image Only**, or **Image and Text**).
4. To change the toolbar layout, deselect the **Auto Toolbar Layout** checkbox.

5. In the *Toolbar Layout* field, select an option from the *Image Text Relation* drop-down list (**Overlay**, **Image Above Text**, **Text Above Image**, **Text Before Image**, or **Image Before Text**).
6. To change the image and/or text size, deselect the **Use Viewer Scale** checkbox.
7. In the *Image/Text Size* field, select an option from the *Image size* drop-down list (**Normal** or **Large**) and an option from the *Text size* drop-down list (**Normal**, **Medium**, or **Large**).
8. Click **OK** to save and exit.


CHAPTER 8—USING THE STUDY BROWSER HELP MENU



The Study Browser's **Help** menu offers users help with using the Study Browser and **Image Viewer** functions. The Help menu contains instructions, current user information, Novarad Support's contact information, the option to submit feedback, the Diagnostic Viewer User Manual, What's New documents detailing new features in the Diagnostic Viewer, and an About section.


The Help menu is located on the Study Browser's main toolbar.

8-1—INSTRUCTIONS



The **Instructions** window offers users basic information about how to open case studies with the Diagnostic Viewer. The Instructions window also offers a Shortcuts section with information on mouse functionality in the Diagnostic Viewer.


8-2—CURRENT USER INFORMATION



Current User Information shows users their role in the Diagnostic Viewer and what features and tools are available for that role.

8-3—CONTACT SUPPORT


Contact Support allows users to submit general questions to Novarad Support.



1. Click **Help** in the Study Browser and select **Contact Support**. The Contact NovaPACS Technical Support dialog box opens.
2. Enter the following information in the message window:
 - a. Name
 - b. Facility Name and Location
 - c. Contact Phone Number and/or E-mail Address
 - d. Complete description of the issue
 - e. Any troubleshooting details, including attempts to correct the issue
3. Click **Send** to email the message to Novarad Support.

8-4—SUBMIT FEEDBACK

Submit Feedback allows users to send product feedback and suggestions to Novarad Corporate.



1. Click **Help** in the Study Browser and select **Submit Feedback**. The Submit Product Feedback dialog box appears.
2. Enter the following information in the message window:
 - a. Name
 - b. Facility Name and Location
 - c. Contact Phone Number and/or E-mail Address
 - d. Complete description of the suggestion/feedback
3. Click **Send** to email the message to Novarad Support.


8-5—USER MANUAL

User Manual allows users to see the Diagnostic Viewer User Manual (in PDF). Users can search for information using the Table of Contents, Index, Find (**CTRL + F**) feature, or the PDF chapter and section bookmarks to the left of the PDF viewer window.

8-6—WHAT'S NEW

What's New allows users to see the NovaPACS What's New document detailing the newest features in that version of the NovaPACS software (as a PDF document).

8-7—ABOUT



About allows users to view Product, Version, and Copyright information about the NovaPACS Diagnostic Viewer as well as provides the user with the Novarad Support phone number.

8-8 FEATURE HIGHLIGHT DIALOG

Every 30 days, users will see the Feature Highlight Dialog pop-up. This will highlight a Feature in PACS or another Novarad product that might be of use to the user. To see the New Feature Dialog more than just once a month, users can click the Party Popper icon located next to the settings button in the Study Browser.



CHAPTER 9—USING THE IMAGE VIEWER



When users open a study from the **Study Browser**, the images display in the **Image Viewer**. The Image Viewer screen displays the images along with easy access to features on toolbars, menus, and buttons. Users can review and annotate images, customize the workspace, and complete reports from the Image Viewer.

9-1—INTRODUCTION TO THE IMAGE VIEWER

When users open a study from the Study Browser, the images display in the Image Viewer. The images, series, and study display in the **Image Viewer** window according to the default Hanging Protocols for that modality. Depending on the monitor configuration, the Image Viewer either displays in front of the Study Browser window or displays on a separate monitor, so both windows are visible.


To access the Image Viewer, double-click a study in the Study Browser window or select a study and click **Open** on that Study Browser list.

To return to the Study Browser at any time from the Image Viewer, click the **Study Browser icon** in the top right of the Image Viewer window, move the mouse cursor to the far left or right of the screen, or press the **Spacebar** (the default hotkey).

Note: Users must select the Auto-show Study Browser checkbox in the Study Browser's Options menu (under the Settings Menu) to access the Study Browser window by moving the mouse to the far left or right of the screen.

Users can access Image Viewer options from menus, toolbars, and hotkey settings. The Image Viewer's main menu is located in the top right corner of the window. Hotkey functions can be assigned to any of these menu functions. Users can access many functions from the Image Viewer's four default toolbars (Main Toolbar, Quick Reference Toolbar, Annotation Toolbar, and Window/Level Toolbar), or they can customize new toolbars.


Note: Users can also access Menu items by assigning a hotkey to that function. Many of these hotkeys are available as default hotkeys.



- **Quick Reference Toolbar**—a toolbar that has options such as **Choose Series** and **Series Layout** for quick access. This toolbar is located to the far left of the Image Viewer window, vertically aligned on the page. However, it can be moved to any other toolbar location in the Image Viewer window.
- **Switch Between Workspaces button**—drop-down list that displays a list of open studies with a checkmark by the currently open study.
- **Patient Name, Study Date, Time of Study, Anatomical Area window**—provides study information.
- **Close Current Study button**—closes the current study but keeps the Image Viewer open.
- **Main Toolbar**—located at the top of the Image Viewer window with drop-down menus (Annotate, Study, Series, Image).
- **Default Hanging Protocol drop-down menu**—contains options for the Image Viewer’s default Hanging Protocol settings.
- **Study Browser icon**—users can click this button to restore the Study Browser window.
- **Image Viewer Menu**—the main menu list in the Image Viewer that contains the Study, Series, Image, Navigation, Customize Toolbars, Tools, and Help menus.
- **Close Application button**—closes the entire Diagnostic Viewer application.
- **Annotation Toolbar**—contains annotation features. Usually located at the bottom of the Image Viewer window.
- **Window/Level Toolbar**—contains Window/Level preset buttons for images. Usually located at the bottom of the Image Viewer window.

Note: Users must activate the four Image Viewer toolbars to view them (Main Toolbar, Quick Reference Toolbar, Annotation Toolbar, and Window/Level Toolbar.) Users can add to or edit the Image Viewer toolbars according to preferences. New toolbars can also be created.

9-2—HANGING PROTOCOL DROP-DOWN MENU



The Hanging Protocol drop-down menu displays the name of the currently applied Hanging Protocol and offers five options users can choose from: Add

study description to current hanging-protocol (HP), Create new hanging-protocol for (Study Description), Save current layout, Create new layout from current layout, and HP (Available Hanging Protocol).

ADD STUDY DESCRIPTION TO CURRENT HANGING PROTOCOL

This option adds the Study Description (Anatomical Area) from the current study to the currently applied Hanging Protocol. In the future, the current Hanging Protocol is applied as the default Hanging Protocol for all studies with the added Study Description. If the user makes changes to the Layout in the Image Viewer prior to adding the Study Description, those changes are saved to the current Layout when the Study Description is added to the Hanging Protocol.

CREATE NEW HANGING PROTOCOL FOR THE OPEN STUDY'S STUDY DESCRIPTION

This option creates a new Hanging Protocol with the Study Description (Anatomical Area) from the current study added. In the future, the new Hanging Protocol is applied as the default Hanging Protocol for all studies with the added Study Description. The new Hanging Protocol is based on the Hanging Protocol that was applied when the option was selected and contains all of the same Layouts as the original Hanging Protocol.

If the user makes changes to the Layout in the Image Viewer prior to selecting this option, those changes are saved to the current Layout in the new Hanging Protocol. The original Hanging Protocol is not modified.

SAVE CURRENT LAYOUT

Make any changes to the Layout in the Image Viewer and select **Save Current Layout** from the drop-down menu. All changes are saved to the current Layout in the current Hanging Protocol. If the user makes changes using the drag-and-drop images or the resize image box features, those changes are saved to the layout. However, filters will remain in their original placement.

CREATE NEW LAYOUT FROM CURRENT LAYOUT

This option allows users to add a new Layout to the current Hanging Protocol directly from the Image Viewer. The new Layout is based on the Layout that was applied when the option was selected. Select **Create New Layout from Current Layout**. With the new Layout applied, make any changes to the new Layout and select **Save Current Layout**. If users make changes to the original Layout in the Image Viewer prior to selecting the Create New Layout from Current Layout option, those changes are saved to the original Layout as well as the new Layout.

APPLY AVAILABLE HANGING PROTOCOLS

After the previous options, all Hanging Protocols available for the open study are listed on the drop-down menu. Selecting a Hanging Protocol from the list applies it to the current study. Although Hanging Protocols associated with a study description are generally set up with a particular type of study in mind, users can apply any Hanging Protocol to the current study that is associated with the same modality, regardless of study description.

9-3—CUSTOMIZING IMAGE VIEWER TOOLBARS

Users can edit the functions on each of the Image Viewer's toolbars, add new functions, and create new toolbars.

1. Right-click any toolbar in the Image Viewer window and select **Customize** from the drop-down list. The Custom Toolbars page opens.
2. Select a toolbar from the **Select Toolbar** drop-down list.
3. Choose an option to the right of the **Select Toolbar** drop-down list:
 - **New Toolbar**—allows users to add a new toolbar
 - **Edit Name**—allows users to edit the toolbar name
 - **Remove Toolbar**—allows users to delete the toolbar

ADDING TOOLBAR ITEMS

1. Select the toolbar from the **Select Toolbar** drop-down list in the Custom Toolbars page.
2. Select the desired feature from the Available Items list on the left side of the page.
3. Click **Add** to add the feature to the Current Toolbar Items list.
4. Click **Save** to save the change(s) and exit.

REMOVING TOOLBAR ITEMS

1. Select the toolbar from the **Select Toolbar** drop-down list in the Custom Toolbars page.
2. Select the unwanted feature from the Current Toolbar Items list.
3. Click **Remove** to remove the feature from the Current Toolbar items list.
4. Click **Save** to save the change(s) and exit.

REORDERING THE MENU

Users can customize the Image Viewer toolbars by resetting the toolbar settings or reordering features.

1. Select the toolbar from the **Select Toolbar** drop-down list in the Customize Toolbars page.
2. Select the feature from the Current toolbar items list.
3. Click **Reset**, **Move Up**, **Move In**, **Move Out**, or **Move Down** until the feature is in the desired location.
4. Click **Save** to save the change(s) and exit.

CUSTOMIZING TOOLBAR OPTIONS

1. Click the **Toolbar Options** tab on the Custom Toolbars page.
2. Select an option from the **Button Layout** drop-down list (**Text Only**, **Image Only**, and **Image and Text**).
3. To change the toolbar layout, deselect the **Auto Toolbar Layout** checkbox.

4. In the *Toolbar Layout* field, select an option from the **Image text relation** drop-down list.
5. To change the image and/or text size, deselect the **Auto Size** checkbox.
6. In the *Image/Text Size* field, select an option from the **Image size** drop-down list (**Normal** or **Large**) and an option from the **Text size** drop-down list (**Normal**, **Medium**, or **Large**).
7. Click **OK** to save and exit when finished.

9-4—MOVING TOOLBARS

Users can move each of the four default toolbars or any created toolbar to the top, bottom, or sides of the Image Viewer window.

1. Right-click any toolbar in the Image Viewer window, and make sure the **Lock Toolbars** option is not checked in the drop-down list.
2. Drag the toolbar by clicking the gray dots to the left of the buttons.
3. Drop the toolbar in the new location.


Toolbars can also be hidden from the Image Viewer window by right clicking any toolbar and deselecting the toolbar name in the drop-down list.

9-5—LOCKING TOOLBARS

Users can select **Lock Toolbars** to lock the toolbars so they cannot be dragged and dropped to a new location. The toolbars remain in their locations unless the user deselects the **Lock Toolbars** option.

Right-click any toolbar in the Image Viewer window and select **Lock Toolbars** from the drop-down list.

9-6—ACCESSING MENU OPTIONS



The *Menu* list in the Image Viewer contains the *Study*, *Series*, *Image*, *Navigation*, *Customize Toolbars*, *Tools*, and *Help* menus. Chapters 10 through 13 offer more details on the Image Viewer's menu functions.

9-7—IMAGE VIEWER ICON BUTTONS



In the Image Viewer window, each image in a series or study displays with icons in the corner that users can access for different functions.


- A. Return to Study Browser

- B. **Display Image in 3D icon**—displays applicable images (CT, MR, PT, XA) in 3D, using Nova3D+
- C. **Scroll to the previous image with an annotation.**
- D. **Annotation icon** – launches the Annotation menu, where users select which type of annotation to use.
- E. **Scroll to the next image with an annotation.**
- F. **Choose Image Layout icon**—opens the Choose Image Layout dialog box to allow users to choose how images are laid out in a series.
- G. **Split Window Horizontally icon**—splits the layout window horizontally into separate windows.
- H. **Split Window Vertically icon**—splits the layout window vertically into separate windows.
- I. **Edit HP icon**—opens the Edit Layout dialog box.
- J. **Close Series icon**—closes only the study or series currently open.
Note: If there is only one series displayed, the Close Series icon is disabled.
- K. **Maximize/Minimize Window Layout icon**—maximizes or minimizes only the study or series currently open.
- L. **Close Window Layout icon**—closes the window layout from view in the Image Viewer window.

9-8—DRAGGING AND DROPPING IMAGES

In the Image Viewer window, users can drag and drop images into any series space on any monitor. This feature helps to facilitate comparison between images in different studies. Users can also drag the side of an Image Box to resize it.


DRAGGING FROM IMAGE BOX



1. Click and hold the top bar of the image display box of the image to be moved (where the Series ID number and Split Window Tools are).
2. Drag it to any series space on any monitor.
3. The image then displays in the selected image box.

DRAGGING FROM CHOOSE SERIES THUMBNAIL

1. Access the Choose Series thumbnail box by clicking **Menu** → **Study** → **Choose Series** or by launching it from the quick access toolbar.



2. Hover the cursor over the thumbnail of the image to be moved.

3. Click and drag it to any series space on any monitor.
4. The image then displays in the selected image box.


DRAGGING IMAGES WITHIN A STUDY

Images that are dragged and dropped to an empty series space within the same study moves to the new display location without any additional labeling or conditions. Only the visible image moves, not the entire series.


Note: If an image is dragged and dropped to an image box that already contains a series, the entire series swaps display boxes with the other.

PATIENT HISTORY SERIES CHOOSER

The Patient History Series Chooser can be used to drag and drop series thumbnails into the image box without having to use the Study Browser. The Patient History Series Chooser contains all of the series thumbnails for both the primary and prior studies and is located at the bottom of the Image Viewer. The Series Thumbnails are grouped by Study with the Study groupings displayed in order based on the study date.



Users can toggle the Patient History Series Chooser On and Off by clicking the gear icon in the left-hand corner of the image viewer and selecting **Patient History Mode**.




Study groups can be collapsed or expanded by clicking the [+] or [-] button in the top left corner of each group.

There is a settings menu in the Patient History Series Chooser that users can access by clicking the blue gear icon on the right-hand side.

Then users will have the option to filter the thumbnails they see using the **Filters** drop down.


The **Max Studies** drop down lets users select the max number of studies they see and select the minimum visible thumbnails.



The **Minimum Visible Thumbnails** drop-down lets users configure how many thumbnails are visible when the study group is collapsed.

DRAGGING IMAGES BETWEEN STUDIES

To drag images to a display box in another study, both studies must be from the same patient. Images that are dragged and dropped to a display box in a different study continues to display in their original study. Only the visible image displays in the new study, not the entire series. Images moved from a comparison study to display in a new study retain all comparison labeling. The moved image continues to display in the new image box even if its original study is closed. Images moved to display in a new study are not saved in the new study or deleted from their original study.



RESIZING IMAGES

To resize an Image Box, hover the cursor over any side of the Image Box that is adjacent to another Image Box. A white, double-pointed arrow appears. Click and drag the side of the Image Box to the desired position. Both of the adjacent Image Boxes are resized.

SAVING CHANGES

Users can save changes made in the Image Viewer to the Layout with the **Hanging Protocol** drop-down list option **Save Current Layout**. All changes that alter image box size or shape are saved. Although content moves with an image box as it is dragged and dropped during viewing, those changes to content placement cannot be permanently saved as part of the permanent Layout. Content placement resets when the layout is reapplied. Content Filters cannot be dragged and dropped and remain in their original placement. To manage content placement in a permanent Layout, use the **Filters** options in the Manage Layout window.

9-9—COPY TO CLIPBOARD

Users can copy an image as it appears in the Image Viewer to the Window Clipboard to paste it onto an email, another document, etc.


FROM THE RADIUS MENU

Select a study that has the image to copy from the study browser.

Click the image to copy. The selected image has a red border.

Right-click to populate the Radius Menu.


Search for **Copy To Clipboard** on the search menu.



Open the desired program and Paste the image by pressing **CTRL + V**. by right-clicking and selecting the **Paste icon**.

The image is pasted without any PHI.

FROM THE IMAGE MENU



Select the study containing the desired image from the study browser.

Click the image to be copied. The highlighted image has a red border.

Select **Menu**.


Hover the cursor over **Image**.

Select **Copy to clipboard**.

The image is copied to the clipboard and can be pasted into another document, application, etc. without any PHI.

CHAPTER 10—USING THE IMAGE VIEWER STUDY MENU


The Image Viewer's **Study** menu contains many functions to help users as they work with images. The Study menu is located on the Image Viewer's Menu list.



Note: The Image Viewer's toolbars can be customized to contain frequently used functions from any menu to increase efficiency. Users can also add any options or functions to Image Context Menus for quick access in the Image Viewer.

10-1—PATIENT HISTORY

Patient History in the Image Viewer allows users to see all of the studies administered to that patient from all time periods and modalities.



In addition, the Patient History window displays the Series, Images, and Report lists as independent lists, and users can drag/drop their lists to any location as they would in the Study Browser.

Users must either manually open the Patient History window in the Image Viewer or choose to have the Patient History window automatically open when opening studies on the Options page.

Manually

1. Open a study in the Image Viewer.
2. Click **Menu**, hover over **Study**, and select **Patient History**.

The Patient History window opens.

Automatically

1. Click **Settings** in the Study Browser and select **Options**.
2. Click the checkbox next to Show patient history when opening a study and click **OK**.
3. Open any study and the Patient History window displays.

Editing the layout of the Patient History window.

1. Once the **Patient History** window is open in the Image Viewer, users can make it any size by clicking and dragging a corner of the window.

The NovaPACS Diagnostic Viewer automatically remembers the size and location of the user's windows in the Image Viewer.


2. As with the Study Browser, if the user clicks a list's tab, the list can be dragged and dropped into any docking panel section.

Lists can be dragged to the top, bottom, left, and right of any section.

Lists can also be dropped into the center of a docking panel section, which adds the list as a tab to the right of an existing list.




Optionally, users can drag lists into their own, independent lists.



Note: If the user has configured their list display on one machine by dragging and dropping and then logs into another machine, they may need to separately configure that machine's lists as well. Users should only need to do this once for each machine from which they access Novarad products.

10-2—SELECT ALL IMAGES IN A STUDY




Select All Series allows users to select all of the series in an active study. This is helpful when wanting to apply changes to multiple series at the same time.

Click the **Menu** in the Image Viewer and select **Study**. Select **Select All Images in a Study**.

Each series is now highlighted in red, showing they are active.

10-3—CLOSE STUDY

Close Study allows users to close the active study while leaving the Image Viewer open. Click **Menu** in the Image Viewer and select **Study**. Select **Close Study**.



Note: Users can also press the red X button next to the study and patient information in the Image Viewer to close the active study while leaving the Image Viewer open.

10-4—CLOSE STUDY AND SHOW STUDY BROWSER

Close Study and Show Study Browser allows users to close the active study and view the Study Browser window in front of the Image Viewer window. Click **Menu** in the Image Viewer and select **Study**. Select **Close Study and Show Study Browser**.

Note: If other studies are open in other workspaces, they remain open. The Study Browser window still shows in front of the Image Viewer window.

10-5—CLOSE WORKSPACE

Close Workspace allows users to close the active workspace. Click **Menu** in the Image Viewer and select **Study**. Select **Close Workspace**.

Note: If other workspaces are open in the Image Viewer, they remain open.

10-6—CLOSE ALL WORKSPACES

Close All Workspaces allows users to close all of the workspaces open in the Image Viewer. Click **Menu** in the Image Viewer and select **Study**. Select **Close All Workspaces**.

Users can also press the default hotkey **CTRL + Shift + F4**.

10-7—CLOSE ALL WORKSPACES AND SHOW STUDY BROWSER

Close All Workspaces and Show Study Browser allows users to close all open workspaces in the Image Viewer and view the Study Browser window. Click **Menu** in the Image Viewer and select **Study**. Select **Close All Workspaces and Show Study Browser**.

Users can also press the default hotkey the **Spacebar**.

10-8—MARK REVIEWED, CLOSE WORKSPACE, AND SHOW STUDY BROWSER

Mark Reviewed, Close Workspace, and Show Study Browser allows users to mark a study Reviewed when closing it and view the Study Browser window.

Note: If the Prompt before marking a study reviewed option or Prompt before marking a study reviewed if it has not been completely reviewed option is selected, a dialog box appears asking the user to confirm they want to mark the study Reviewed. If other studies are open in other workspaces, they remain open. The Study Browser window appears in front of the Image Viewer window.

Click **Menu** in the Image Viewer and select a **Study**. Select **Mark Reviewed, Close Workspace, and Show Study Browser**.

10-9—MANAGE HANGING PROTOCOLS

Navigate to **Menu** → **Study** → **Presentation** → **Manage Hanging Protocols**. The Manage Hanging Protocols editor dialog box opens.

ADD NEW HANGING PROTOCOL

Click **Add New** at the bottom of the Study Types list. The Manage Hanging Protocol editor now allows users to select options for the new Hanging Protocol.

1. Select Modality

First, select a modality from the **Modality** drop-down list. This is the modality with which the new Hanging Protocol is associated. A modality must be selected in order to save a new Hanging Protocol.


2. Manage Study Descriptions

Next, manage any Study Descriptions to be associated with the new Hanging Protocol. A Hanging Protocol is applied automatically to any study that matches all of its Study Descriptions. Adding Study Descriptions also helps users distinguish what types of studies the Hanging Protocol was set up to display.

To add a Study Description click **Add** to the right of the Study Descriptions list. A small text box appears inside the Study Descriptions list. Type in an anatomical area to be associated with the new Hanging Protocol. Repeat as many Study Descriptions as needed. A unique Study Description (or combination of Study Descriptions) must be added in order to save a new Hanging Protocol associated with the same modality as an existing Hanging Protocol.

To edit an existing Study Description, highlight the description and click **Edit** to the right of the Study Descriptions list.

To delete an existing Study Description, highlight the description and click **Remove** to the right of the Study Descriptions list.




3. Select a Display Tool Option

Select an option from the **Display Tool** drop-down list to specify how to display the Hanging Protocol.

4. Select a Spine Labeling Default

Select an option from the **Spine Labeling Defaults** drop-down list to specify how to configure spine labeling for the new Hanging Protocol.



Select **Edit Options** to the right of the **Spine Labeling Defaults** drop-down list to customize the available spine labeling options. The **Spine Labeling Options** editor allows users to add, edit, and remove spine labeling options. Users can also customize font and font color for spine labeling options.

MANAGE LAYOUTS

Users can add, edit, remove, and duplicate layouts within a Hanging Protocol using the **Presentation Groups (layouts)** section of the **Manage Hanging Protocol** editor.

To add a new **Layout** to a **Hanging Protocol**, click **Add New** at the bottom of the **Study Types** list. The **Manage Hanging Protocol** editor allows users to select options for the new **Hanging Protocol Layout**.

1. Add Layout


Click **Create layout** to the top right of the layouts section. The system adds a new, blank **Layout**.

2. Edit Layout


Click **Edit Layout** to the top right of the layouts section or double click the thumbnail image of a layout. The **Edit Layout** window opens.

Editing Series Layout

There are three primary ways to edit the **Series Layout**.




1. Click the box with green grid lines located above the image boxes.




The **Choose Series Layout** box opens.

Select how to display the series by hovering the cursor vertically and horizontally over the squares and clicking when the preferred layout is highlighted light gray.




2. Or click **Split Window Tools** located in the top bar of an image box to edit the image layout within that image box or to split the series horizontally or vertically.
3. Or click the red X at the top of an image box to remove that window from the series layout.



Using the Imagebox (Series) Tab

First, select the Image Box to be edited. This enables the Imagebox (Series) tab and allows users to choose image layout options for the selected portion of the layout.



On the Imagebox (Series) tab, users can edit the View, Image Layout, and Clip Playback.

Tiled

The images in the series can be set to be in a Tiled display. This is the default view.

When the Tiled option is selected, users can do the following:

1. Use the Choose Image Layout box to select how images display within the series. Hover the cursor vertically and horizontally over the squares and click when the preferred layout is highlighted light gray.
2. Zoom To Fit, To Fill, Custom, or No Zoom.
3. Adjust the Frame rate, which controls the playback speed.
4. Select the Auto-start and Use DICOM Frame Rate checkboxes.
5. Select Looping, Sweeping, or Stop for Playback Sequencing.

Cine


The images in a series can be set to a Cine display. Cine is typically used for Cardiology images.

When the Cine option is selected, users can do the following:

1. Set the Desired Frame Rate, which controls the playback speed.
2. Select the Auto-Start checkbox, which begins to play the clip as soon as the image loads.
3. Select the Use DICOM Frame Rate checkbox.

4. Set the Playback Sequencing to Looping, Sweeping, or Stop


Processed By



The images in a series can be set to be processed by either **MPR** or by **Image Fusion**. Select **Processed by** from the drop-down list after selecting the **Processed View**.

If users select the **Processed by MPR** option, the series automatically displays in a **Multi-Planar Reconstruction (MPR)** view. This view reconstructs multiple 2D images into a multi-planar 3D view. Novarad provides MIP/MPR functionality and support with NovaPACS.

Note: If the system has not been used with the Diagnostic Viewer's MIP/MPR support, an administrator must contact Novarad Support to configure each workstation to access MIP/MPR.



If **Processed by Image Fusion** is selected, two series from different modalities (typically CT and PET) fuse for a more detailed image.

Using the Filters Tab

Select the **Image Box** to edit. This enables the **Filters** tab and allows users to choose **Content Filter** options for the selected portion of the layout.


Users can edit options for the **Filters** tab in the **Edit Layout** window. The **Filters** tab allows users to **Add**, **Edit**, and **Remove** **Content Filters** that manage which files display in individual **Image Boxes**.

Adding Content Filters

Click **Add** to the right of the **Content Filters** list. This enables the **Edit Filter** box.

Editing Content Filters

Add a new filter or select an existing filter from the **Content Filters** list and click **Edit**. The **Edit Filter** options allow users to **Filter by Attribute Presence** or **Filter by Image Plane**.



1. To Filter by Attribute Presence, select the checkbox next to the option's title. This enables the rest of the options associated with this filter. Select whether to include content if the attribute is present or absent. Select the Attribute to filter by from the drop-down list. The VR value is entered automatically. Click **Apply**. New filters now appear in the Content Filters list.

Note: The Image Details window displays all available attributes and their values for a selected image. Users can refer to this window if they have difficulty setting up filters.


2. To Filter by Attribute Value select the Attribute to filter by from the drop-down list. This enables the rest of the options associated with this filter. The VR value is entered automatically.

Note: Do not select the Filter by Attribute Presence checkbox when filtering by Attribute Value.


Select the desired Operation from the drop-down list. The operation specifies the relationship between the Attribute and the Value. Select the desired Value from the drop-down list (or type it). Click **Apply**. The new filter now appears in the Content Filters list.

3. To Filter by Image Plane, select the checkbox next to the option's title. This enables the rest of the options associated with this filter. Select whether to include content that is a Member of or Not a Member of an image plane from the Operation drop-down list. Select the Image Plane to filter by from the Value drop-down list. Click **Apply**. New filters now appear in the Content Filters list.


Note: Not all images are associated with an Image Plane. If an image plane is associated with an image the data is stored in the image's Series Description. Users can view an image's Series Description from the Series List on the Study Browser.



Users can also view an image's Series Description from the Image Details window.



4. Users may also select whether they would like the rest of a series or other content that is associated with a filtered image but does not meet filter requirements to be displayed in the same Image Box.



Deleting Content Filters

Select an existing Content Filter from the Content Filters list and click **Remove**.

Display Multiframe as Series

To view multiframe images as their own series for a specific modality, open *Manage Hanging Protocols* inside the *Image Viewer*. Either select a pre-existing hanging protocol or create a new hanging protocol. Check the checkbox marked *Display Multiframe as Series*.



When finished editing a layout, click **OK** or **Cancel** to close the *Edit Layout* window and return to the *Manage Hanging Protocols* editor.

3. Remove Layout

Select the layout's thumbnail image from the *Manage Hanging Protocols* editor and click **Remove layout**. The selected layout is deleted from the current *Hanging Protocol*.

4. Duplicate Layout

Select the layout's thumbnail image from the *Manage Hanging Protocols* editor and click **Duplicate layout**. The selected layout is duplicated in the current *Hanging Protocol*.

EDIT HANGING PROTOCOL

Click an existing *Hanging Protocol*'s description in the *Study Types* list. The *Manage Hanging Protocols* editor allows users to make changes to any options described in the previous sections.

REMOVE HANGING PROTOCOL

Select an existing *Hanging Protocol*'s description and click **Remove** below the *Study Types* list.

When finished managing *Hanging Protocols*, click **Close** to exit the *Manage Hanging Protocols* editor and return to the *Image Viewer*.

10-10—ADD SPLIT

Add Split allows users to split the *Image Viewer*'s viewing area in half. The new viewing space opens to the right of the current study. If the user has more than one study open from a patient (open in tabs at the top of the *Image Viewer* window), the new viewing space automatically populates the second study. Users can also click and drag an open image into the new window.

Click **Menu** in the *Image Viewer* and select **Study**. Select **Add Split**.

10-11—CHOOSE SERIES

Choose Series allows users to see all the series in the current study in a thumbnail view list. Users can select other series to view from this list.

Click **Menu** in the Image Viewer and select **Study**. Select **Choose Series**.

The Choose Series list disappears from view after interacting with the Image Viewer unless the **Stay Visible** checkbox is selected at the top of the list. The Choose Series list can be moved anywhere in the viewing area by clicking and dragging it to a new location. Users can resize the Choose Series list box by hovering the cursor in the corner of the box until the cursor becomes a double-ended arrow and then clicking and dragging to the desired size.


1. Drag the thumbnail image from the Choose Series list into the series window.
2. Select an image window in the viewing area. It highlights red. Double-click the thumbnail image from the Choose Series list and it opens in that selected image window.

10-12—CHOOSE SERIES BOX DOCKING


The Choose Series box in the Image Viewer enables users to quickly open a study's individual series—by clicking the preview thumbnails or dragging and dropping. Users can dock their Choose Series box so it is in the most convenient location while viewing a study's images. Users can dock the Choose Series box to the top, bottom, left, or right of the Image Viewer window.

1. Click **Menu** → **Study** → **Choose Series** in the Image Viewer.

Note: If users dock the Choose Series box or select the Stay Visible option, they do not need to select the Choose Series box option again. It opens and displays automatically for all studies.




2. Click the **gear icon** in the top, right-hand corner of the Choose Series box. Select the option to **Dock Left**, **Dock Right**, **Dock Top**, or **Dock Bottom** to dock the Choose Series box.



Once docked, users can choose to **Undock** or select a new docking option at any time.


All studies, series, images and toolbars display in the Image Viewer along with the docked Choose Series box.



3. If users open a current and comparison study together, each study's Choose Series box displays with its respective study. Each Choose Series box can be independently docked/undocked.

Note: If users undock the Choose Series box, the option to Stay Visible is available when they click the **gear icon**. If the Choose Series box is docked, the Undock option is available when users click the **gear icon**.

10-13—AUTO SYNC



Auto Sync allows users to set images to synchronize as they are scrolled (in coplanar series). If the Auto Sync option is selected, it remains on for all subsequent studies until deselected.

Click **Menu** in the Image Viewer and select **Study**. Select **Auto Sync**. Users can also press the default hotkey **S**.

Scroll through synchronized images by either:

1. Left-clicking an image and dragging the cursor up or down.
2. Pressing the left or right arrow keys on the keyboard.

The synchronized images in the other windows scroll in sync with the image.

10-14—CROSS LOCALIZE

Cross Localize allows users to view a line of intersection between two orthogonal images. A linear reference marker appears in all open series windows that correlate to the selected image. Images will still be cross-localized as the user scrolls through them.

1. Select an image to Cross Localize. More than one series needs to be viewable in the Image Viewer window.
2. Click **Menu** in the Image Viewer and select **Study**. Select **Cross Localize**. Users can also press the default hotkey **X**.

The linear reference marker appears in all visible series that correlate to the selected image.

Users can adjust the angle of obliqueness when using Cross Localize lines in MPR or Fusion view. After enabling Cross Localize lines, there are rotational red dots on the ends of the lines. Click these dots and drag the line in any direction. This changes the oblique angle of the corresponding selected image. A small dot near the middle of the cross-localize line indicates the center of rotation in the oblique view.

Note: When Radiologists activate the Cross Localize feature for a specific modality, the system remembers and saves whether the feature is turned on or off for that specific modality as a user preference. The state of Cross Localize is independently saved as a user preference for each modality.

10-15—CROSS LOCALIZE

Point Cross Localize allows users to select a point on an image, and the Image Viewer finds the intersecting images in all orthogonal series. The Auto Sync option must be activated to use Point Cross Localize.

1. Select an image with which to use Point Cross Localize.
2. Click **Menu** in the Image Viewer and select **Study**. Select **Point Cross Localize**.
3. Right-click the specific point in the selected image. A bullseye symbol appears on that point; the same point is located in the other images in the Image Viewer window.

10-16—SET SYNC

Set Sync allows users to link series from different studies to a frame reference of choice. If the user wants to compare two separate studies side-by-side, they can choose the specific image in a series to start synchronization using Set Sync.

The Set Sync option enables Auto Sync. After Set Sync is set, the Image Viewer synchronizes the image slices as the user scrolls.

Note: The Set Sync option affects how the two series synchronize with each other; the Auto Sync within each series remains the same.

1. Open the comparative studies so the Image Viewer displays both.
2. Scroll through the images in each series until the images you want to synchronize are both visible and selected (highlighted red).
3. Click **Menu** in the Image Viewer and select **Study**. Select **Set Sync**.

10-17—AGGREGATE

Aggregate allows users to select images from various studies or series and display them in a single series window. The Aggregate series created is temporary and does not exist after the window closes.

1. Open images from one or more series in the Image Viewer window.
2. Click **Menu** in the Image Viewer and select **Study**. Select **Aggregate**. An empty Aggregate window opens.

3. Click the top bar of the series (where the Series ID number and Split Window Tools are). Hold **CTRL** while dragging and dropping the series into the Aggregate window. Users can select any number of series to drop into the **Aggregate** window.

10-18—PRINT

Print allows users to print images. Users can choose to print to a default printer or select another printer. This allows users to print to paper or film.


SELECTING A PRINTER

1. Click **Menu** in the Image Viewer and select **Study**. Select **Print**. The Print Wizard dialog box opens.

Note: If the user is printing specific images, select the images before selecting the Print option (images highlight in red when selected). Otherwise, all images in the active study or series will print.

2. Click **Print using a DICOM printer I will select**.
3. Select an option under the **Which images will be printed?** section. Users can choose between printing all open images in the active series or study or printing selected images in the active series or study.
4. Select an option for printing clips. Users can choose to print only the active frame or all frames.
5. To print each series on a new page, select the **Start each series on a new page** checkbox.
6. Click **Next** when finished. The Printer Options dialog box appears.
7. Select a printer from the **Printer** drop-down list.

Note: If the needed printer is not on the **Printer** drop-down list, contact an Administrator. Printers and their printing formats are configured by the Administrator, Novarad Support, and the printer's vendor.



8. Click the black **Image Format** window to select the format.
9. Select format options from the drop-down lists and checkboxes.


Note: The printer's abilities and the vendor settings determine the printing formats.

10. The options to print image and overlays are selected by default. Deselect the checkboxes if needed.
11. Click **Next**. A final page opens. Users can click **Back** to make any changes or click **Finish** to print.

PRINTING TO THE DEFAULT PRINTER


1. Click **Menu** in the Image Viewer and select **Study**. Select **Print**. The Print Wizard dialog box opens.
2. Select **Print using my computer's default printer**.

3. Select an option under the **Which images will be printed?** section. Users can choose between printing all open images in the active series or study or printing selected images in the active series or study.
4. Select an option for printing clips. Users can choose to print only the active frame or all frames.
5. To print each series on a new page, select the **Start each series on a new page** checkbox.



6. Click **Next** when finished. The Printer Options dialog box appears.
7. Click the black **Page Layout** window to select the layout.
8. Complete the **Company Information** section to have it printed at the top of each page.
9. The options to print image labels and overlays are selected by default. Deselect the checkboxes if needed.
10. Select a Zoom option (**Fit** or **Fill**).
11. Click **Next**. A final page opens. Users can click **Back** to make any changes or click **Finish** to print.

10-19—ROUTE



Route allows users to share images with other sites or locations. Images can be routed manually from a NovaPACS data source to a specified destination.

1. To route a study series, click **Menu** in the Image Viewer and select **Study**. Select **Route**. The **Add Route Queue for** dialog box opens.
2. Choose the study, series, and/or image(s) to route by expanding the lists and clicking the corresponding

checkboxes next to the items. If the checkbox next to the study is marked without expanding the list, all of the accompanying series and images route automatically.

Note: Route rules can be configured in the Admin Console. If one series in a multi-modality study matches the route rule, the entire study routes (e.g., a mammography study can include images from multiple modalities like MG and US, and also non-image files such as Structured Reports).

3. Choose a destination from the **Destination** drop-down list.

Note: An Administrator must add the Destinations that appear in the drop-down list. If the needed destination isn't on the list, talk to an Administrator.

4. Set the routing priority for the image routing by clicking the up or down arrows next to the **Priority** field. Routing priorities range from 1-100, with 1 as the highest priority. The default routing Priority is 3.
5. Click **OK** to route to the specified destination.

10-20—DOCUMENTS

Documents allows users to see any documents associated with the patient file that were scanned into the RIS.

Click **Menu** in the Image Viewer and select **Study**. Select **Documents**.

10-21—EXAM REPORT

Exam Report allows users to see the report associated with the patient file if one exists. Exam Reports are only available if the system is integrated with Nova RIS.

Note: When reports are finalized in Nova RIS, a DICOM PDF report is sent to NovaPACS and is added to the study as a new series. If no matching NovaPACS study is found, a placeholder study is created for the report. Admins must locate this study in the NovaPACS Admin Console on the Studies tab and move the report series to the correct study manually.

Click **Menu** in the Image Viewer and select **Study**. Select **Exam Report**.

Note: PDF reports can be viewed but not edited. If authorized to do so, users can edit HTML reports.

Note: When users sign and save reports in NovaPACS, the preliminary report displays in all report locations depending on configuration. Once the report is finalized in Nova RIS, a DICOM-encapsulated PDF report is sent to NovaPACS and added to the study. This DICOM PDF report displays in any location where users access reports.


10-22—MANAGE REPORTS WORKFLOW

Manage Reports Workflow allows users to see a list of reports and their statuses. This option is particularly useful for Radiologists, Residents, Administrators, and Transcriptionists.

In the Report Workflow window, users see four main boxes for **Dictated Reports**, **Unsigned Reports**, **Signed Reports**, and **Finalized Reports**. In the Other section, users can see reports that have been marked Requires Correction and can look up physicians in the system using the Physician Lookup box.

Click **Menu** in the Image Viewer and select **Study**. Select **Manage Reports Workflow**.

10-23—RIS INFORMATION



RIS Information allows users to view the RIS Information dialog box for a study. The information included in the RIS Information box is submitted by the Admin, Receptionist, or Technologist in Nova RIS. RIS Information dialog boxes include patient information, referring physician information, order notes, and procedure information about a study.

Open a study and click **Menu** → **Study** → **RIS Information**.

The RIS Information dialog box opens.

Note: Information contained on the RIS Information dialog box only displays if it has been added to a Nova RIS order form when scheduling. Referring physician contact information must be added to the referring physician record in Nova RIS before it displays on the RIS Information dialog box.


10-24—DICTATE

Dictate allows users to dictate a report for the selected study. The following sections describe dictation options for the NovaPACS Diagnostic Viewer.

DICTATION TYPE

The Dictate settings are set by Novarad in user preferences. The options available are: **Basic**, **File**, **Dragon**, and **CMDLine** (found on the Plug-Ins tab of the Options page, under the Study Browser's Settings Menu).

LAUNCH DICTATION



Click **Menu** in the Image Viewer and select **Study**. Select **Dictate**.


Users can also press the default hotkey **D**.

Note: If a prior study is made active and dictation is launched, users receive a warning dialog box asking if they are sure they want to dictate the study.

DICTATING AN ADDENDUM (REPORT EXISTS)

A dialog box alerts users if there are existing dictations for the study. If the user dictates for a study that already has a report, it is added as a report addendum.

DICTIONATION MICROPHONE SETTINGS CHECK



The Diagnostic Viewer checks the default microphone settings before starting dictation and notifies the user if the microphone is muted or the recording volume is too low. Users should make any needed adjustments before beginning dictation.

Note: The Diagnostic Viewer only checks the microphone settings once per application session, not each time dictation launches.

DICTIONATIONS WAITING TO BE UPLOADED

The dictation dialog box also lists the number of dictations waiting to be uploaded, in the bottom right-hand corner of the dictation window.

This helps notify radiologists that dictation files have been saved and are being automatically uploaded to Nova RIS through the Dictation Launcher service in the event that they are experiencing network or connection issues.

FLUENCY DIRECT DICTIONATION INTEGRATION

Fluency Direct allows users to insert patient and exam data into the diagnostic report using voice commands. These commands can be used in the Instant Dictation window or within a Standard Report.

The commands insert the specified data at the current cursor position.

Voice Command Syntax Details

Fluency Direct recognizes specific voice commands enclosed in caret symbols (^). These commands allow radiologists to automatically insert patient and study data into the report.

Syntax Rules

- All commands must start with ^ and end with ^.
- Format: ^Novarad[DataType]^
- Say the full command clearly (e.g., “Novarad patient name”).
- The system replaces the command with the actual data when recognized.

Available Commands

Voice Command	Replaced With	Example Output
^NovaradModality^	Study Modality	CT

Voice Command	Replaced With	Example Output
<i>^NovaradStudyDescription^</i>	<i>Study Description</i>	<i>CT Chest with Contrast</i>
<i>^NovaradPatientAge^</i>	<i>Patient Age</i>	<i>54</i>
<i>^NovaradPatientName^</i>	<i>Patient Name</i>	<i>John Smith</i>
<i>^NovaradPatientGender^</i>	<i>Patient Gender</i>	<i>Male</i>
<i>^NovaradPhysicianReason^</i>	<i>Physician Reason</i>	<i>Shortness of breath</i>

Usage Notes

- *Commands can be inserted anywhere in the dictation at the cursor location.*
- *If the data does not exist for the current study, the command is typically left unchanged.*
- *Pause briefly before and after saying a command for best recognition.*
- *Commands work in both the Instant Dictation window and Standard Reports.*

INSERTING NOVA RIS STANDARD REPORTS

When using Dragon dictation, the user can insert a RIS Standard Report at any point in the current dictation.


TRANSCRIPTION BUTTONS

*The M*Modal Fluency Direct Function is mapped to actions in the Instant Transcription Window. The functions are as follows:*

- *Launch Dictation “D” from Novarad Diagnostic Viewer*
- *Record (hold microphone button or click the MIC/VU meter control in the Instant Transcription window)*
- *Play*
- *Pause*
- *Skip Forward*
- *Skip Backward*
- *Save*
- *Sign & Save*

The Instant Transcription Window’s size and position can be saved as a user preference. When dictating, move or resize the form, and when users close and reopen the form, the size and location will be saved.

CHANGING DICTATION FONT DISPLAY SIZE



There is a font display slider bar in the Instant Transcription window when using Dragon or M*Modal dictation in the NovaPACS Diagnostic Viewer.

1. Log into the PACS as a Radiologist, Resident, or Cardiologist, open a new study, and launch Dictation using the default hotkey **D**, or using a configured Speech mic/Nuance mic button. The Instant Transcription window opens.
2. The font display slider bar appears just above the text display section.
3. By moving the slider to the right, the font display size in the Instant Transcription text window increases; moving the slider to the left decreases the font display text size.

Note: Altering the font display size in the Instant Transcription window does not affect the text display on any diagnostic report for the patient—diagnostic report text display is dependent upon the formatting used in the reporting system.

10-25—LAUNCH MAMMO VIEWER

Launch Mammo Viewer allows users to launch a mammo viewing tool used with NovaPACS for mammography viewing.

1. Launch the Diagnostic Viewer and open an appropriate study (i.e., an MG or breast US study).
2. Click **Menu** in the Image Viewer, select **Study**, then select **Launch Mammo Viewer**. The configured mammo viewer launches in a new window.
3. When finished using the mammo viewer, users can return to the PACS by closing the window.

Note: A mammo viewer software must be installed and configured on the workstation to use the Launch Mammo Viewer option. If it is not configured, the Launch Mammo Viewer option is not available on the Study Menu. Users can customize a hotkey or add a custom toolbar button to the Image Viewer or Study Browser to access this option.

10-26—VIEW STUDY ON OTHER WORKSTATION

View Study on Other Workstation allows users to send studies from the Diagnostic Viewer to another computer. Contact Technical Support for the required configuration for this option. Once this is configured, studies can be viewed and processed by another computer. This option is useful if users need to view studies on post-processing stations for MIP/MPR or 3D Reconstruction.

Click **Menu** in the Image Viewer and select **Study**. Select **View Study on Other Workstation**.

10-27—LAUNCH ORTHOPEDIC TOOLS

Launch Orthopedic Tools allows use of the orthopedic template third-party software with Novarad's Diagnostic Viewer. This optional software integration allows users to create a complete array of orthopedic pre- and post-surgical planning and mapping tools including (but not limited to) Cobb Angles, templates, overlaying, and surgical repairs.

1. Launch the Diagnostic Viewer and open an appropriate study.
2. Click **Menu** in the Image Viewer, select **Study**, then select **Launch Orthopedic Tools**. The configured orthopedic tool launches in a new window.
3. When finished using the orthopedic tools, users can return to the PACS by closing the window.

Note: The product launched when using the Launch Orthopedic Tools option is dependent upon which software option is configured in the PACS Admin Console. If an orthopedic software tool is not configured, the Launch Orthopedic Tools option is unavailable on the Study Menu. Users can customize a hotkey or add a custom toolbar button to the Image Viewer or Study Browser to access this option.

10-28—ASSIGN PASSWORD

Assign Password allows users to assign a temporary **Single-Study Access Password** for a study, allowing a person (i.e. a Referring Physician or patient) to access that study through the Internet.

1. Click the Study to be assigned Single-Study Access.
2. Click **Menu** in the Image Viewer and select **Study**. Select **Assign Password**.
3. The Assign Web Password to Selected Study dialog box appears. A randomly generated password is assigned automatically. Users can generate their own password by clicking the **New Password** link.
4. Click the drop-down list next to Expiration Date to set the length of time the password is valid and click **OK**.

10-29—NOTES

Notes allows users to view notes that have been made for the selected study. The Notes option allows various Diagnostic Viewer users (technologists, radiologists, etc.) to communicate with each other about the study. The Notes option is not part of the patient's medical record, does not print, and does not transfer to a patient disk.

Click **Menu** in the Image Viewer and select **Study**. Select **Notes**. Users can also press the default hotkey **N**.

10-30—COLLABORATION

Collaboration allows multiple radiologists to work simultaneously from the same worklist without having duplicate readings. The Collaboration option also allows users to participate in peer-to-peer reviews of patient studies.

Studies just coming into the NovaPACS system have a status of **New** or **Imported**, depending on the origin of the study. Studies with a **New** status are coming from modalities into the NovaPACS system. If users have an external

source for their studies, such as a long-term archive (LTA) or other DICOM source, studies imported and cataloged into the NovaPACS system have a status of *Imported*.

The *Collaboration* options enable users to mark various statuses such as *Verified*, *Checked Out*, and *Reviewed* to manage studies in collaborative environments.

CHECKOUT

Checkout prevents multiple users from reviewing the same study simultaneously.

Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and select **Checkout**.

The study displays *Checked Out* in the *Status* column of the *Filtered Studies* list in the *Study Browser*. The *Current User* column reflects the username of the person who checked the study out. The *Checked-Out* status for the study updates on all collaborating users' lists. Other users may view a *Checked-Out* study but cannot mark it reviewed or use the *Checkout* function on it.

UNDO CHECKOUT

Undo Checkout allows users to undo the checkout process for a study. This function is helpful when the user checks out the wrong study or the user is unable to complete the review of a study after checkout. The *Undo Checkout* function cannot restore the status of a study if it has already been marked *Reviewed*. Users can only change a study with a *Reviewed* status using the *Reset Status to New* function (see section **Reset Status to New** below).

Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Undo**.

The study or series is restored to the status it had before checkout (i.e. *New* or *Verified*).

VERIFIED

Verified allows users to mark studies as verified, which means that all attached documents, images, etc., are checked and ready for review.

Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Verified**.

The study now has *Verified* in the *Status* column of the *Filtered Studies* list. The *Current User* column reflects the username of the person who verified the study. The *Verified* status for the study updates on all collaborating users' lists.

REVIEWED

Reviewed allows users to mark a study as reviewed.

Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Reviewed**.

The study now displays *Reviewed* in the *Status* column of the *Filtered Studies* list in the *Study Browser*. The *Current User* column reflects the user name of the person who reviewed the study. The *Reviewed* status for the study updates on all collaborating users' lists.

Note: The status of a study resets to its previous *New* or *Verified* status if the study is closed out without marking it *Reviewed* and the following options are marked in the Study Browser's Settings menu under the *General* tab:

- Auto-checkout studies when they are opened
- Auto-review studies when they are closed
- Prompt before marking a study reviewed

RESET STATUS TO NEW

Reset Status to New allows users to restore the status of a reviewed study back to *New*, so it can be reviewed again.

Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Reset Status to New**.

The study's status changes from *Reviewed* to *New* in the Status column of the Filtered Studies list in the Study Browser. The username no longer appears in the Current User column.

Note: The *Reset Status to New* function must be used within four hours of a study being marked *Reviewed*. To change the status of a study after four hours, contact Novarad Support.

PEER REVIEW

Peer Review allows Radiologists to evaluate each other's reviewed studies, give second opinions, and provide quality assurance. Any Radiologist involved in the review process can view a completed Peer Review.

1. Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Peer Review**. The Peer Review dialog box appears.
2. Click an option for the **Position** field.
3. Enter any comments needed in the **Comment** field box and click **OK**.

Note: Peer Review reports are available to NovaPACS Administrators in the Admin Console. The Administrator can choose to report an average rate of missed diagnoses over a given time period.

UTILIZATION REVIEW

Utilization Review allows Radiologists to review film to provide quality assurance.

1. Click **Menu** in the Image Viewer and select **Study**. Select **Collaboration** and click **Utilization Review**. The Utilization Information for Reviewed Study dialog box appears.
2. Select the applicable position from the **Utilization Heading Position** field drop-down list.
3. Enter any comments needed in the **Comment** field box and click **OK**.

10-31 TEACHING FILE

Teaching File allows users to create anonymized **Teaching Files** based on existing studies in NovaPACS. These anonymized files do not contain personal identifiers but maintain the layout and structure of the original studies,

and all previously made annotations are preserved when the new file is created. Once a Teaching File has been created, the new file is completely independent of the original study, and unique patient IDs are generated to prevent Teaching Files from being associated with original studies.

Note: Users can only create Teaching Files by enabling Teaching Files in the Options tab.



For this option to function correctly, users must first select individual images to add to the new Teaching File. Select the **Teaching File icon** to add individual images from a study to a new Teaching File.

When users select this icon, images are added to a Teaching File series window where users can review the images added to the file.

Note: Individual images in comparison studies cannot be added to Teaching Files. To create Teaching Files of images in a comparison study, close the study and reopen it as the primary study.


After reviewing the images in the viewing window, users can decide to remove images from the Teaching File. To do this, navigate back to the Image Viewer and select the **Teaching File icon**. Images are then removed from the Teaching File Images window.

After selecting images, navigate to the **Study** drop-down menu and select **Teaching File**. The Teaching File Information dialog box appears, allowing users to enter study and patient information.

Users can also access the Teaching File Information dialog directly after selecting individual images by hitting their spacebar and exiting the study. The dialog appears automatically.

When users create **Teaching Files** by selecting images individually, the dialog box allows them to click **Include Selection From Study** or **Include Entire Study**. **Include Selection From Study** will create the file from only the files that were selected. **Include Entire Study** creates the file from all of the images in the study.

Mark as in Progress allows users to select some images without creating the teaching file. They can come back to that study at any time, finish selecting images, and then create the teaching file.



Note: **Teaching Files** originally created from an entire study do not provide users with the option to only include selected images in a Teaching File.

10-32—RESEARCH

Research allows users to look in the GlobalRad database for cases that have similar anatomical/pathological findings. This can help with more accurate diagnoses and more efficient research.

Once activated, the Research option launches GlobalRad—a database of studies from around the world that have been made anonymously, complete with clinical findings, diagnoses, etc. Novarad users can sign up for a GlobalRad account and submit studies using the Upload to GlobalRad option. GlobalRad preserves anonymity of patient information as it is submitted. For any assistance in this process, please contact Novarad Support.

Click **Menu** in the Image Viewer and select **Study**. Select **Research**.

10-33—UPLOAD TO GLOBALRAD

Upload to GlobalRad allows users to upload all selected images in the active study to GlobalRad—a database of studies that have been made anonymous for research purposes.

Click **Menu** in the Image Viewer and select **Study**. Select **Upload to GlobalRad**.

10-34—SAVE VIEW SETTINGS

Save View Settings allows users to save the view-related changes (**W/L, Zoom, Pan**, etc.) that were made while viewing the active study/series. If the user does not select **Save View Settings**, the Diagnostic Viewer does not save any view-related changes.

Click **Menu** in the Image Viewer and select **Study**. Select **Save View Settings**.

10-35—SHOW ALL SERIES

Show All Series allows users to view all of the series for the current patient in the Image Viewer window.

Click **Menu** in the Image Viewer and select **Study**. Select **Show All Series**.

10-36—SERIES LAYOUT

Series Layout allows users to alter the series layout for the active series using a Choose Series Layout window.

1. Click **Study** in the Image Viewer. Select **Series Layout**.
2. Choose how to display the series by hovering the cursor vertically and horizontally over the squares and clicking when the preferred layout is highlighted in light gray.

10-37—HANGING PROTOCOL LAYOUTS

Hanging Protocol Layouts can be applied to a study from the Image Viewer in a number of ways.

- **Apply Next Layout** allows users to iterate to the next available Layout for the current Hanging Protocol. Navigate to **Menu** → **Study** → **Apply Next Layout**. Users may also set up a hotkey for this option. Using **Apply Next Layout** allows users to navigate sequentially through Layouts until they reach the last available Layout for the current Hanging Protocol. This option does not allow users to cycle back to the beginning.

- **Apply Previous Layout** allows users to navigate to the previous Layout for the current Hanging Protocol. Navigate to **Menu** → **Study** → **Apply Previous Layout**. Users may also set up a hotkey for this option. Using **Apply Previous Layout** allows users to navigate backwards through Layouts until they reach the first available Layout for the current Hanging Protocol. This option does not allow users to cycle back to the end.
- **Hanging Protocol Layouts** allows users to apply a preset Hanging Protocol Layouts layout to the active study/series. It also allows users to open the Hanging Protocol Layouts window to preview the set Hanging Protocol Layouts for that modality before choosing to apply one.
 - a) Click **Study** in the Image Viewer. Select **Hanging Protocol Layouts**. A drop-down list opens.
 - b) Select **View Hanging Protocol Layouts** to preview the existing Hanging Protocol Layouts for the modality and select the thumbnail image to apply to the active study.

OR

- a) Select a **Hanging Protocol Layout** from the drop-down list to apply it to the active study. These numbered Hanging Protocols correspond with the available Layouts in the current Hanging Protocol as well as the number keys on the keyboard (i.e. Hanging Protocol 1 corresponds with the first available Layout in the current Hanging Protocol as well as the number one on the keyboard). Press a number key on the keyboard to quickly apply its associated layout.

Note: These numbers are configurable in the Study Browser's Options, under the Keyboard Shortcuts tab.

10-39—MINIMIZE

Minimize allows users to minimize the Image Viewer window.

Click **Study** in the Image Viewer. Select **Minimize**.

10-40—EXIT

Exit allows users to exit and close the entire Diagnostic Viewer. If users have selected the Prompt before closing the application option from the Study Browser's Options list, a prompt appears when selecting **Exit**.

Note: Closing the Diagnostic Viewer using the Exit option closes studies properly and saves all changes to patient images. Use this exit option or the Study Browser's exit application to ensure proper study closing and saving.

Click **Study** in the Image Viewer. Select **Exit**.

CHAPTER 11—USING THE IMAGE VIEWER SERIES MENU

The Image Viewer's **Series Menu** contains many functions to help users as they examine and mark findings in a series. The Series menu is located on the Image Viewer's Menu.

Note: The Image Viewer toolbars can be customized to contain frequently used functions from any menu to increase efficiency. Users can also add any options or functions to Image Context Menus for quick access in the Image Viewer.

11-1—SELECT ALL IMAGES

Select All Images allows users to select all images in the active series. This option is helpful if users want to make changes to all the images in a series at the same time.

Click **Menu** in the Image Viewer and select **Series**. Click **Select All Images**. Users can also press the default hotkey **CTRL + A**.

11-2—CLOSE

Close allows users to close the active series. All other open series remain open.

Click **Menu** in the Image Viewer and select **Series**. Select **Close**.

11-3—VIEWS

Views allows users to change views on the active study.

TILED

The Tiled view is the default view that presents the series with the images and clips in a tiled format. The clips in the Tiled view each have their own controllers and play independently of each other.


Note: Users can pause and/or scroll through the frames of a Cine clip in a Tiled Layout.

Click **Menu** in the Image Viewer and select **Series**. Select the **Views** drop-down list and select **Tiled**.

CINE

The cine view displays all of the images in a series in a played format, including the first frame of each image.

Note: If clips are opened using a cine view, users see a message stating that the cine view is used to cine a series of images and that the Tiled view can play clips. Users are encouraged to switch to a Tiled view.



Click **Menu** in the Image Viewer and select **Series**. Select the **Views** drop-down list and select **Cine**.

To adjust the Cine or Clip Frame Rate:

1. Click the frame rate number found on the cine or clip controller bar. A slider bar appears.
2. Use the slider bar to adjust the frame rate. Moving the slider down decreases the frame rate. Moving the slider up increases the frame rate.

MPR VIEW

The MPR (Multi-Planar Reconstruction) view reconstructs an axial plane into sagittal and coronal planes to provide more comprehensive views of images.

Note: Novarad provides MIP/MPR integration and support. If the system is not configured to use the Diagnostic Viewer with MIP/MPR, an administrator should contact Novarad Support for assistance configuring the workstation to access MIP/MPR. MIP is Maximum Intensity Projection.

Click **Menu** in the Image Viewer and select **Series**. Select the **Views** drop-down list and click **MPR**.

There are four MPR views. Three of the views display a specific orientation (axial, coronal, or sagittal). The fourth view has several buttons that, when used, modify the view or the functionality of the left mouse button.

Users can quickly change/reset the fourth MPR view by clicking an orientation button (axial, sagittal, or coronal).

The following list contains a few of the MIP/MPR options:

- Rapid scroll through the series by left clicking and dragging the mouse.
- Change the slice (slab) thickness by scrolling the mouse wheel.
- Change Window/Level by right-clicking and dragging the mouse.
- Zoom by pressing and holding both the left-mouse and right-mouse buttons and dragging.
- Right-click to bring up an MPR Image Context Menu.
- Capture Series allows users to set the parameters for the capture. The capture starts at the green line and ends at the red. Drag those lines around the image for a more precise capture. Users can also change the desired slice spacing.

Users also have the ability to create an Angle Measurement in the MPR View. To create this annotation, follow the steps below:

1. Click **Menu**.
2. Under the Annotations option, click **Measure Angle**.
3. Draw the lines start to end by clicking the right-mouse button and dragging.
4. Users can drag the end points to adjust the angles.


The system will compute by matching the two start points and measuring the angle between the endpoints.

IMAGE FUSION

The Image Fusion view displays CT and PET images together, along with a row of images with fused (overlaid) CT and PET images.

Note: Image Fusion requires a separate download that an administrator should run before launching the Diagnostic Viewer. This is a one-time download per computer.

1. Open a study with more than one series
2. Click **Menu** in the Image Viewer and select **Series**. Select the **Views** drop-down list and select **Image Fusion**. The **Select Series for Fusion** dialog box appears.
3. Select the two series to fuse. Both series selected need to have the same number of images.



In the Image Viewer, the fused PET/CT series shows in three rows. The first two rows display the PET and CT series, and the third row displays the fused series.

Right-click and the Image Context menu appears with options including **Colorize**, **Measure Line**, and **Capture Image**.

Note: Users can individually customize the W/L values of the images in each row. If users modify the W/L values of the fused images, the Diagnostic Viewer revises both the CT version in the first row and the PET


version in the middle row. If users alter W/L values of just the second row of the PET version of the images, the fused images on the bottom row only adjust the PET portion.

To configure the mouse controls to navigate through the PET study more efficiently, see section **7-1 Options** for more information.

SUV MEASUREMENTS WITH IMAGE FUSION

Users can now make SUV (standardized uptake value) measurements with Image Fusion in NovaPACS.

1. Launch the Diagnostic Viewer and open a study to measure in the Image Viewer.
2. Click **Menu** in the Image Viewer, select **Series**, and select **Views**.
3. Select **Image Fusion** from the drop-down list. The **Select Series for Fusion** dialog box appears.
4. Select two series to fuse.
5. Once the Image Fusion series is opened, right-click an image to measure.
6. Select **Sphere VOI** (volume-of-interest) from the drop-down list.



7. Right-click the image and complete the measurement. When the mouse is released, the SUV measurements display on the image.
8. To save the measurements, right-click the image and select **Capture Image** from the drop-down list. This creates a Fusion Captures series.

Note: The Viewer will display SUV values when users Window/Level a PET image. The lower image bound will not go below zero when window/leveling PET images to avoid washing out the image.

11-4—SPLIT HORIZONTAL

Split Horizontal allows users to split the active window in the Image Viewer horizontally. The next series in succession appears in the new window when it is split.

Click **Menu** in the Image Viewer and select **Series**. Select **Split Horizontal**. Users can also click the **Split Horizontal icon** in the upper right corner of each active window.

11-5—SPLIT VERTICAL

Split Vertical allows users to split the active window in the Image Viewer vertically. The next series in succession appears in the new window when it is split.

Click **Menu** in the Image Viewer and select **Series**. Select **Split Vertical**. Users can also click the **Split Vertical icon** in the upper right corner of each active window.

11-6—CHOOSE IMAGE LAYOUT


Choose Image Layout allows users to choose how the images in a series are arranged.

Click **Menu** in the Image Viewer and select **Series**. Select **Choose Image Layout**. Users can also click the **Choose Image Layout icon** in the upper right corner of each active window.

The Choose Image Layout dialog box opens. Select how to display images by moving the cursor vertically and horizontally over the squares and clicking when the preferred layout is highlighted light gray.

11-7—CALCIUM SCORING

Calcium Scoring allows users to examine CT studies for calcium buildup in the coronary arteries. Normally, the coronary arteries do not contain calcium. Calcium in the coronary arteries is a sign of Coronary Artery Disease (CAD). Cardiac Calcium Scoring, or coronary artery calcium scoring, uses CT studies to check for the buildup of calcium in the coronary arteries to detect early stages and severity of heart disease.



1. Open a CT study in the Image Viewer.
2. Click **Menu** in the Image Viewer and select **Series**. Select **Calcium Scoring**. The Calcium Scoring dialog box appears.
3. Select an **Artery Button** in the Calcium Scoring dialog box:

LM: Left Main

LCX: Left Circumflex

LAD: Left Anterior Descending

DIAG: Diagonal

RCA: Right Coronary

PDA: Posterior Descending

Other: Other

4. Complete the Calcium Scoring annotation on the image for each artery to measure.
5. When the Calcium Scoring annotation is complete, click **Create Report** to create a web page report.

Note: This web page report does not save in the Diagnostic Viewer or Nova RIS.

11-8—SORT IMAGES

Sort Images allows users to change the sort order of the images in a series. The default sort order is sort By Image Number Ascending.


Click **Menu** in the Image Viewer and select **Series**. Select the **Sort Images** drop-down list, and choose a sort option:

- **By Image Number Ascending**—sorts the images numerically starting from the first image.
- **By Image Number Descending**—sorts the images numerically starting from the last image.
- **By Slice Location Ascending**—sorts the images by slice location starting from the last slice taken.
- **By Slice Location Descending**—sorts the images by slice location starting from the first slice taken.

11-9—COLORIZE

Colorize allows users to identify areas of interest in radiology images using colorization. This helps radiologists to spot or validate areas of interest in images by converting grey scale intensity to colors.

1. Open a study in the Image Viewer. Colorization works best on studies with a lot of detail.
2. Click **Menu** in the Image Viewer and select **Series**. Select the **Colorize** drop-down list and choose a color.



11-10—CREATE NEW VIEW

Create New View allows users to open a second view of the active series in a new series window. This option is useful when comparing different images within a series or comparing the same image with different window/level values.

Click **Menu** in the Image Viewer and select **Series**. Select **Create New View**.

To see the new view, scroll to the next image or use the **Split Window** tools to split the series window horizontally or vertically to see the image and the new view together.

11-11—SUBTRACT

Subtract allows users to subtract the pixels of the selected image from all other images in the series. When subtracting is complete, the series shows the pixels from all other images in the series.


Cardiologists can subtract an angio study (XA) and see the contrast through the veins once the image with bones and soft tissue is subtracted from the series.

Note: The Subtract option is available to Cardiologists only.

Click **Menu** in the Image Viewer and select **Series**. Select **Subtract**.

11-12—ADVANCED ANALYSIS

Using the Advanced Analysis option, users can generate the maximum and minimum slope for each image slice over time. This is an important feature for Radiologists because it allows them to see the areas with the steepest rate of change (maximum) or the slowest rate of change (minimum). Being able to quickly see the rate of change is helpful in a lot of cases, including but not limited to:



Maximum Slope:

- Differentiating between malignant and benign tumors.
- Mapping brain functions during cognitive or motor tasks.
- Detecting areas of impaired blood flow in the heart.
- Identifying areas of acute ischemia

Minimum Slope:


- Identifying stable, non-reactive tissue such as scar tissue or chronic lesions.
- Indicates regions with poor blood flow or low perfusion, which is common in areas suffering from chronic ischemia or vascular disease.
- Diagnosing neurodegenerative diseases.
- Assessing stable vascular abnormalities such as an aneurysm or arterial stenosis.

To calculate the Maximum or Minimum slope:

1. Open the desired study.
2. Click Menu in the Image Viewer.
3. Navigate to Series.
4. Click Advanced Analysis from the side menu that appears.

The Advanced Analysis window will appear.

5. Select the desired operation.
6. Click OK.



Users can also subtract Images from the series with this feature.

Subtract Options:

In the Advanced Analysis window, users can also subtract images from the image series. Subtracting an image will create a transparent image overlay. This is helpful because it allows users to isolate and highlight the differences between two images, making it easier to detect subtle changes or abnormalities over time. This process can help identify tumors, inflammation, or vascular changes by removing static elements and emphasizing areas of contrast enhancement or other dynamic tissue alterations.

A. Subtract selected image from the others in the series:

Choosing this option will subtract the currently selected image from the series and create a transparent overlay that can be used to compare images.

B. Subtract images of the selected time from the others in the series by slice location:

This option will allow users to select a group of images by their time distinction and subject from the series, which will create a transparent image overlay that can be used for comparison.

Results: Users can decide how they want the subtraction results displayed. They can choose to display them as a new series with a new name or in the existing series.

Original Series: This displays the number of images, the number of slice locations, and the number of images ignored in the original series.

Result Series: Displays the number of images in the Result Series.


11-13—REMOVE IMAGES

Remove Images allows users to remove selected images from a **Key Image** series or **Aggregate** series. This option does not remove images from a regular study.

Click **Menu** in the Image Viewer and select **Series**. Select **Remove Images**.

11-14—EJECTION FRACTION

Ejection Fraction allows users to measure the capacity at which the heart is pumping. Ejection Fraction is the percentage of blood that pumps out of a filled ventricle with each heartbeat. The Ejection Fraction option is usually measured in the left ventricle (LV)—the heart's main pumping chamber. A normal Ejection Fraction percentage is 55 to 70 percent.




Note: This feature is only for Cardiologists.

1. Open an image in the Image Viewer (typically MR) that shows the ventricles of the heart.
2. Navigate through the images in the active series to find an image of the left ventricle in the Systolic Phase.
3. Click **Menu** in the Image Viewer and select **Series**. Select **Ejection Fraction**. The Calculate Ejection Fraction dialog box appears.
4. Use the **Next Frame** and/or **Previous Frame** buttons to find an image of the left ventricle in the Systolic Phase.
5. To measure the ventricle at its smallest size, right-click and drag the mouse on the image. Outline or draw a circle around the left ventricle in the Systolic Phase.
6. Click **Next** when done. Use the **Next Frame** and/or **Previous Frame** buttons to find an image of the left ventricle in the Diastolic Phase.
7. To measure the ventricle at its largest size, right-click and drag the mouse on the image. Outline or draw a circle around the left ventricle in the Diastolic Phase.
8. Click **Next** to see an Ejection Fraction Report.
9. Click **Save** to save the Ejection Fraction Report.

11-15—PERCENT STENOSIS

Percent Stenosis allows users to check for the presence of a stenosis (abnormal narrowing) in a blood vessel. Calculating the Percent Stenosis tells a Radiologist how narrow the blood vessel is compared to a normal section of the blood vessel.



Note: This feature is only for Cardiologists.

1. Open an Angio study (XA).

Note: An image without any pixel spacing must be calibrated before an accurate Percent Stenosis can be calculated.

2. Navigate through the images in the active series to find an image showing the normal width of the blood vessel.

3. Click **Menu** in the Image Viewer and select **Series**. Select **Percent Stenosis**. The Calculate Percent Stenosis dialog box appears.
4. To measure the vessel at its widest point, right-click and drag the mouse on the image. The cursor begins in the Measure Length mode.
5. Draw a line extending from one border to the next, from lumen wall to lumen wall, on the fully patent or open lumen. Click **Next** when finished.
6. To measure the vessel at its most narrow point, right-click and drag the mouse on the image. This marks the distance of the obstruction.
7. To measure the lesion length, right-click and drag the mouse on the lesion in the image.
8. Click **Next** to see the Percent Stenosis report.
9. Click **Save** to save the Percent Stenosis report.

CHAPTER 12—USING THE IMAGE VIEWER IMAGE MENU

The Image Viewer's **Image** menu contains many functions to help users as they view images. Users can select from many annotation options and identify critical findings. The Image menu is located on the Image Viewer's Menu. When configured, the Image menu can also be located on the Image Viewer's Main Toolbar.

Note: The Image Viewer toolbars can be customized to contain frequently used functions from any menu to increase efficiency. Users can also add any options or functions to Image Context Menus for quick access in the Image Viewer.

12-1—ZOOM

Zoom allows users to change the zoom or magnification of an image.

Click **Menu** in the Image Viewer and select **Image**. Select the **Zoom** drop-down list.

Users have the choice of four different Zoom options:

- **1:1** - Displays the image in its original size from the modality. Each pixel of the image is represented by one pixel on the monitor. The size and resolution of the monitor affects how large the image appears.
- **To Fit** - The image scales until the entire image is visible and touching either the horizontal or vertical edges of the view box. This is the default zoom setting.
- **To Fill** - The image scales until the image touches both horizontal and vertical edges of the view box. This results in parts of the image being cut off. Monitor orientation, image orientation, and hanging protocol settings affect how much of the image is visible.
- **Manual Zoom** - Users can hold the right and left mouse buttons simultaneously while dragging the mouse up to zoom in and down to zoom out. The default hotkey is **Z**.


Note: Undo the Zoom changes by manually zooming out with the mouse. The entire image then shows.

12-2—ZOOM MAP

When users zoom into an image in the NovaPACS Diagnostic Viewer, they see an image thumbnail displayed in the bottom right-hand corner of the image called a **Zoom Map**. The Zoom Map shows users the area they are viewing in respect to the entire image. The Zoom Map is clearly marked with a yellow box pinpointing the area where users have zoomed in.

1. Log into the NovaPACS Diagnostic Viewer and open an image.
2. With an image open in the Image Viewer, press and hold the left and right mouse buttons simultaneously and drag the mouse up to zoom in.

The Zoom Map automatically displays in the bottom, right-hand corner of the image.



Note: To zoom out, click both mouse buttons and drag down.

3. The Zoom Map feature is enabled by default for all users. To disable the Zoom Map, click **Menu** in the Image Viewer.
4. Hover the cursor over **Image** → **Zoom** and deselect **Zoom Map**.

Note: The Zoom Map feature is enabled by default. If the user disables/enables this feature, it saves to their user's preferences.

12-3—RELOAD

Reload allows users to reload the image file, removing any changes in orientation, magnification, filtering, or pixel inversion.

Click **Menu** in the Image Viewer and select **Image**. Select **Reload** from the list. Users can also press the default hotkey, the **Backspace** key.

12-4—INVERT

Invert allows users to invert the grey scale shades in an image to the colors of a negative image where the black and white colors switch. The Invert option aids in visualization of angiography and other vessel/lumen exams and runoffs.

Click **Menu** in the Image Viewer and select **Image**. Select **Invert**. Users can also press the default hotkey **I**.

12-5—SUBTRACT CLIP FRAMES


Subtract Clip Frames allows users to subtract the active frame of a clip from the other frames in the clip. The Diagnostic Viewer creates a new clip containing the subtracted frames and adds it to the series.

Click **Menu** in the Image Viewer and select **Image**. Select **Subtract Clip Frames**.


12-6—CALIBRATE MEASUREMENT

Calibrate Measurement allows users to calibrate and measure studies that do not come with pixel spacing as part of their DICOM value.


1. Click **Menu** in the Image Viewer and select **Image**. Select **Calibrate Measurement**.



2. A dialog box appears, warning users of the possible inaccuracies of calibrating measurements by hand. This warning dialog box appears each time users open the Calibrate Measurement option unless the **Don't show this again** checkbox is marked. Click **OK**.




3. The cursor changes to an annotation measurement cursor. Right-click on the image and drag the mouse to draw a line of measurement.



4. In the Measurement Calibration dialog box, choose the desired measurement value—users can choose from **Use Millimeters** or **Use French (1/3 mm)** value. Select or enter the measurement value.
5. To apply the image calibration to the entire series, click the **Apply to Series** checkbox. This calibrates all of the images in the currently selected series.
6. Click **OK** to save the calibration setting.

Note: This calibration is only valid for the current image unless the Apply to Series checkbox has been selected.

12-7—DETAILS



Details allow users to view the list of DICOM fields that the modality transmitted when sending the study to NovaPACS.

Click **Menu** in the Image Viewer and select **Image**. Select **Details**. The DICOM Viewer dialog box opens.

12-8—FLIP/ROTATE

Flip/Rotate allows users to change flip or rotate an image or images.

Click **Menu** in the Image Viewer and select **Image**. Select the **Flip/Rotate** drop-down list. Users have the choice of four different Flip/Rotate options:

- **Flip Horizontal**
- **Flip Vertical**
- **Rotate Right**
- **Rotate Left**

Note: Users can flip/rotate an entire series by choosing **Select All Images** from the Image Viewer's Series menu before using the Flip/Rotate option.

Note: The Reload, Hide, Invert, and Rotate toolbar buttons are available in the web viewer as image specific features that only affect the feature to which they were applied.

12-9—ANNOTATE

Annotate allows users to perform functions like annotating, marking, highlighting, and marking regions of interest on images. There are many annotation options in the Annotate list.

Note: Admin users that do not also have a Radiologist, Resident, or Cardiologist role can annotate, but the annotations are not saved.

ACCESSING THE ANNOTATE OPTIONS

1. Click **Menu** in the Image Viewer and select **Image**. Select **Annotate** to open the drop-down list.
2. If the Image Viewer's Main Toolbar is activated, click **Image**, and select **Annotate** to open the drop-down list.
3. Right-click an image in the Image Viewer. The Image Context Menu appears containing many annotation options by default.
4. Activate the Default Annotation Toolbar by right-clicking any toolbar in the Image Viewer and selecting **Default Annotation Toolbar**.

USING ANNOTATIONS

1. Select an annotation tool from the **Annotate** drop-down list or toolbar.
2. The mouse cursor changes to reflect the selected annotation. Activate the annotation by right-clicking the image and moving the mouse.
3. Perform the action specific to the annotation option.
4. Reset the cursor when finished by selecting **Reset Cursor** from the **Annotate** drop-down list or toolbar. Users can also reset the cursor by pressing the default hotkey **Esc**.

ANNOTATION OPTIONS

Users have many annotation options to choose from when working with images. Annotations can be added to any image, even in Cine view.

COLOR

Color allows users to select a color for annotation marks from a color palette.

FONT

Font allows users to select the font, size, and style for any text annotations.

RESET CURSOR

Reset Cursor allows users to turn off any selected annotation mode and restores the mouse cursor to the default setting. **Double-Right-Clicking** with the mouse or using the default hotkey **Esc** also resets the cursor.

ROI

ROI (Region of Interest) measures Hounsfield values (or the density of a particular area on the image) of the average and deviation values to aid in diagnosis. The ROI annotation is helpful in tumor diagnosis and identifying abnormalities in the anatomy.

Select the **Annotate** drop-down list, select the **ROI** drop-down list, and choose an option:

- **ROI Polygon**
- **ROI Rectangle**
- **ROI Circle** (default hotkey is **R**)
- **ROI Ellipse**
- **ROI Point**

Users can also add 2D ROI annotations to PET/Fusion Images by right-clicking on a PET/Fusion Image and selecting Annotation. To add a 2D circle, select 'Circle.'

PEN

Pen allows users to make annotations using a freehand drawing tool. The default hotkey is **W**.

POLYGON

Polygon allows users to create a freehand polygon shape. Any side left undrawn is auto completed. This allows users to highlight an area on an image.

RECTANGLE

Rectangle allows users to draw a rectangle shape to highlight an area on an image.

CIRCLE

Circle allows users to draw a circle shape to highlight an area on an image.

On an MPR image, the circle annotation will only cover the circular area on that specific slice; unlike the sphere annotation, it will not cover the area on any other slice. This will also display the minimum, maximum, medium, and average values for the area.

ELLIPSE

Ellipse allows users to draw an oval shape to highlight an area on an image.

POINT

Point allows users to draw small dot shapes to highlight areas on an image.

ARROW

Arrow allows users to create arrows for pointing out specific features on images. The default hotkey is **A**. With the arrow annotation selected, right-click and drag to draw the arrow. By default, arrows are drawn tail to head, meaning that the tail end of the arrow is positioned where the user first clicks and the pointed end of the arrow is dragged to the desired position before releasing the mouse. Users can reverse the way arrow annotations are drawn in the Study Browser's Option menu.

COORDINATES

Coordinates allows users to determine and annotate the coordinates of a specific point on an image. This tool is useful when counting and comparing vertebrae from one plane to another.

Click the **Annotate** drop-down list and select **Coordinates**. Right-click the precise point to establish coordinates. To verify that the coordinates or anatomy are of equal position on comparable images, repeat these steps on those images.

TEXT

Text allows users to add text annotations to images.

1. Select the **Annotate** drop-down list and select **Text**. Users can also press the default hotkey **T**.
2. Right-click the image to mark the starting point. The Insert Text dialog box opens.
3. Enter new text for the annotation or select pre-existing text from the lexicon.
4. To add the new text to the lexicon, select the checkbox for **Add this annotation to the lexicon**.
5. Users can edit the existing lexicon text and/or add new lexicon items by clicking **Edit Lexicon**.
6. Click **OK** to add the text annotation to the image.

SPINE LABELING

Spine Labeling allows users to annotate vertebral body and disc mapping.

1. Select the **Annotate** drop-down list and select **Spine Labeling**. The Spine Labeling dialog box appears.
2. Select the **Vertebral Body** or **Disc** tab to identify the type of label.
3. Select the starting point for the annotation (i.e. C3).
4. Select the **Up** or **Down** option for the direction of labeling. This marks the spine labels from either superior to inferior or inferior to superior.
5. Right-click the image to mark the starting point or right-click and drag to draw a line on the border of the area. If the user draws a line on the border, then draws a second line, the Diagnostic Viewer automatically connects the ends of the lines to form a four-sided polygon; the label text is automatically centered in the polygon.
6. Continue marking the vertebrae in the image. Do not close the menu until the labeling is complete.

To use the **Spine Label Cross Reference**:

1. Click **Settings** in the Study Browser and select **Options**.
2. On the **General** tab, select **Display spine label cross references**.

Spine label cross-references show at the top of intersecting images. An intersecting image is any orthogonal in the study (with the same frame of reference) that intersects a spine label in another image.

If it is a simple label, it intersects the label text. If it is a polygon label, it intersects any part of the polygon (i.e. if a sagittal image of the spine is labeled, any axial images that intersect a spine label displays the label at the top of the axial image for reference).

MEASURE LENGTH

Measure Length allows users to measure lengths on images.

1. Select the **Annotate** drop-down list and select **Measure Length**. Users can also press the default hotkey **M**.
2. Right-click the image to mark the starting point.
3. Drag the cursor the length of the area to measure, and release. A line appears with the measurement in millimeters.

Note: The Image Viewer can only calculate measurements for images with pixel spacing information. If no pixel information exists, users must Calibrate Measurements before measuring.

MEASURE PATH

Measure Path allows users to measure a freehand line. This tool gives a more accurate measurement of non-linear objects or anatomy.

Click the **Annotate** drop-down list and select **Measure Path**. On the image, right-click and drag the cursor along the path to measure. The measurement displays in millimeters when the path is complete.

MEASURE AREA

Measure Area allows users to trace around anatomical areas on an image to calculate the area and circumference. The area and circumference values display on the image in millimeters.

MEASURE ANGLE

Measure Angle allows users to measure the angles of structures and other elements in an image. This is a useful tool for scoliosis and fracture angles.

1. Select the **Annotate** drop-down list and select **Measure Angle**. Users can also press the default hotkey **CTRL + M**.
2. Right-click and drag to draw the first line.
3. Move the mouse cursor to the fulcrum or hinge point of the angle and right-click and drag to draw another line. An arc fills in between the two lines established and the degree of the angle is displayed at the fulcrum point.

Users can also complete the following Measure Angle options:

- Cross two lines in the middle
- Draw two lines that do not touch but would intersect if they were to extend farther
- Draw two parallel lines

Note: Angle measurements can be undone one line at a time using the Undo annotation tool.

COBB ANGLE

Cobb Angle allows users to perform complex angle measurements for spinal curvature and scoliosis studies along with other angle measurements for Orthopedics diagnosis and reading.

1. Select the **Annotate** drop-down list and select **Cobb Angle**.
2. Right-click the image and drag the mouse to begin using the Cobb Angle annotation tool.
3. Right-click and drag the mouse to draw as many lines as needed. The angles between the lines are measured and displayed in degrees.

VOICE ANNOTATIONS

The Voice Annotation Dialog is a powerful feature that allows users to record or dictate a voice annotation in conjunction with other annotation marking tools. The Image Viewer records mouse movements over the image along with the voice recording, allowing users to point out and mark relevant pathology on the screen while describing it.

1. Select the **Annotate** drop-down list and select **Record**. The Voice Annotation Player dialog box appears, and the recording begins immediately. Users need a microphone for this tool.
2. Dictate the information for the image. The mouse movements on the screen are also recorded.
3. Click **Save** to end the recording session and save the recording. The recorded annotation is saved to the image. A Voice label appears in the corner of the image once recorded.

The dialog will present Patient and Study information. It will also display a position indicator (to inform users how long the recording is). There are also Pause, Play, Record, Save and Cancel.

PLAY

Play allows users to play back information that was recorded with the Record annotation tool. The recording plays back in the Image Viewer when opened.

SHOW SCALE

Show Scale allows users to display a scale on images to help judge sizes. The scale displays to the bottom right of the image. The scale is in cm units and changes to mm units if the image is zoomed in a certain amount.

Note: The Show Scale option does not save to the image, but when displayed, it prints with the image.

HIDE ALL

Hide All allows users to hide all computer-generated and user-generated annotations and labels on the selected image. This option remains on or off unless selected.

ERASE

Erase allows users to erase the user-created annotations from the selected image.

Note: This option only erases annotations created by the current user during the current viewing session in the Image Viewer. To erase saved annotations, see **History** below.

ERASE ALL

Erase All allows users to delete all the user-created annotations from the selected image.

Note: This option only erases annotations created by the current user during the current viewing session in the Image Viewer. To erase saved annotations, see **History** below.

UNDO

Undo allows users to remove the last annotation made to the selected image.

REDO

Redo allows users to reinstate the last annotation removed.

HISTORY

History allows users to open a history of all annotations created for the selected image. The Annotation History dialog box opens. The four columns in the dialog box are as follows: **Annotation Type**, **Creation Time**, **Created By**, and **Can Delete**. If the Can Delete column says Yes for an annotation, users can select the annotation and click **Delete** to remove it.

EDITING ANNOTATIONS

Users may quickly and easily edit annotations for images without having to delete or undo the original annotation. Annotations may be edited only during the current viewing session. Once an annotation has been saved, it cannot be edited. A small padlock symbol appears next to the cursor if an annotation cannot be edited. Users may edit annotations by clicking and dragging or using the **Edit Annotations** drop-down list.

CLICK AND DRAG

Users can click and drag annotations in the Image Viewer to move an existing annotation, move text independently from an existing annotation, or to resize or reshape an existing annotation.

MOVING ANNOTATIONS

To move an existing annotation, hover the cursor over the annotation until a hand-shaped cursor appears. Click and hold the left mouse button. Drag the annotation to the desired location and release the mouse. If an annotation makes calculations based on location, the calculations automatically update once the annotation is moved.

MOVING TEXT

If an annotation includes associated text, the text box can be moved independently from the main annotation. Hover the cursor over the text until a hand-shaped cursor appears. Click and hold the left mouse button. Drag the text to the desired location and release the mouse. The main annotation remains in place while the associated text is moved. If the main annotation is moved, the associated text moves with it. Users cannot edit text in an existing annotation; however, text based on calculations update to reflect changes to the annotation.

RESIZING OR RESHAPING ANNOTATIONS

Some annotations can be resized or reshaped. Typically, if an annotation has multiple lines or figures, each item can be resized/reshaped individually. To resize or reshape an annotation, hover the cursor over an existing annotation until a hand-shaped cursor appears. Follow the line of the annotation with the mouse until a small, white circle appears next to the cursor (usually at the ends or corners of an annotation). Click and hold the left mouse button. Drag the cursor to resize or reshape the annotation. Release the mouse. The edited annotation remains at the selected size or shape. Any calculations based on the size, shape, or positioning of the annotation update to reflect changes.

Note: Annotations that are drawn freehand or that only represent one point cannot be resized or reshaped.

DROP-DOWN MENU

Users can right click an existing annotation to access a drop-down list that allows them to delete the annotation, change the color of the annotation, or change the font of text used in an annotation.

CHANGING COLOR

To change the color of an annotation, hover the cursor over the annotation until a hand-shaped cursor appears. Right-click the annotation. A drop-down menu appears. Select **Color**. The Color dialog box

appears. Select the desired color and click **OK**. The color of the annotation automatically updates in the Viewer.

CHANGING FONT

To change the font of a text annotation, hover the cursor over the text until a hand-shaped cursor appears. Right-click the annotation. A drop-down menu appears. Select **Font**. The Font dialog box appears. Select the desired font and click **OK**. The annotation font automatically updates in the Viewer.

DELETE ANNOTATIONS


Users can delete annotations by right-clicking the annotation and selecting **Delete Annotation** from the drop-down menu. The selected annotation deletes from the image in the Viewer. Only annotations from the current viewing session can be deleted using this option.

12-10—NAVIGATE THROUGH ANNOTATED IMAGES

Users can quickly navigate through annotated images in the Image Viewer in the NovaPACS Diagnostic Viewer.

In the header bar for the series window in the Image Viewer, there are annotation navigation icons. Users can click the left arrow icon to navigate to the previous annotated image or click the right arrow icon to navigate to the next annotated image.

Each series box displayed in the Image Viewer window has these arrow icons.



The Image Viewer navigates to the next or previous image in the user's current series first. If there are no more available annotated images in the current series, the next available series containing annotated images appears.

12-11—MOVE IMAGES

Move Images allows users to transfer images in and out of Aggregate and/or Key Image Series.

Click **Menu** in the Image Viewer and select **Image**. Select **Move Images**. Right-click and drag to move selected images in or out of Aggregate and Key Image series.

12-12—EXPORT TO FILE

Export to File allows users to save an image to a file for export or other uses (i.e. exporting images to PowerPoint or sending fetal ultrasound images with a patient). The Export to File option supports DICOM, BMP, GIF, TIFF, and JPEG file formats.

1. Select an image from an open study to export.
2. Click **Menu** in the Image Viewer and select **Image**. Select **Export to File**. The Save As dialog box opens.
3. Browse for the export location by expanding the folders in the **Save in** drop-down list.
4. When the correct location is selected, click **Save** to export the image.

Note: The Diagnostic Viewer only exports image files.

12-13—MAGNIFICATION

Magnification allows users to activate the Magnify Glass feature. Click **Menu** in the Image Viewer and select the **Magnification** drop-down list. Users have the choice of four different Magnification options:

- **Magnify Glass**—turns the magnifying glass on or off. The default hotkey is **G**.
- **Magnify Glass Size**—increases or decreases the size of the magnifying glass window. The default hotkey is **Shift + G**.
- **Magnify Up**—increases the size of the magnifying glass window.
- **Magnify Down**—decreases the size of the magnifying glass window.

12-14—UNDO ALL SELECTIONS

Undo All Selections allows users to deselect any selected images in the Image Viewer.

Click the Image Viewer's **Menu** and select **Image**. Select **Undo All Selections**.

12-15—KEY IMAGE SERIES

A **Key Image Series** is a series that contains certain key images that the user wants to group together. When a Key Image Series is created, the Diagnostic Viewer saves this series for future viewing.

- **Key Images**—shows or hides the Key Image Series in the Image Viewer window.
- **Add to Key Images**—adds the selected image to the Key Image Series.
- **Remove from Key Images**—removes the selected Key Image from the Key Image Series.

To view the Key Image Series, click **Menu** in the Image Viewer and select **Image**. Select the **Key Image Series** drop-down list and select **Key Images**. Users can also press the default hotkey **K**.

To add to the Key Image Series, click **Menu** in the Image Viewer and select **Image**. Select the **Key Image Series** drop-down list and select **Add to Key Images**. Users can also click the **key icon** at the top of each image in the Image Viewer.

Note: Users cannot add files to an existing Key Image Series if the Key Image Series isn't already open.

To remove an image from a Key Image Series, select the key image in the series. Click **Menu** in the Image Viewer and select **Image**. Select the **Key Image Series** drop-down list and select **Remove from Key Images**.

12-16—CRITICAL FINDINGS

Critical Findings allows radiologists to mark images, series, and/or studies that contain a critical. Critical findings can contain up to ten diagnoses.

ADD CRITICAL FINDINGS

1. Select an image to mark.
2. Click **Menu** in the Image Viewer and select **Image**. Select the **Critical Findings** drop-down list and select **Add Critical Findings**. Users can also click the **critical findings triangle icon** at the top of each image in the Image Viewer. The Add Critical Findings dialog box opens.
3. Select from the list of available Critical Findings.
4. Click **Save** and yellow text labeling the Critical Finding appears in the top-right corner of the image in the Image Viewer.

Note: If the user is both a radiologist and an admin, an Add a New Critical Finding Value link appears. After selecting this option, the Add New Critical Finding dialog box appears. Enter text for the Critical Finding and click **OK**. The new value appears in the available critical findings list and can now be selected. Click **Save** to add the finding to the image.

In the Add Critical Findings dialog box, users can also choose to do the following:

1. **Add image to the Key Image Series**—adds the selected image to a Key Image Series.
2. **Record Findings**—opens the Findings Record Editor dialog box where users can record the critical findings and mark any physician notification attempts.

RECORD FINDINGS

1. Click **Menu** in the Image Viewer and select **Image**. Select the **Critical Findings** drop-down list and select **Record Findings**. The Findings Record Editor dialog box opens.
2. Click **Add** to add a finding. The Physician Notification Editor dialog box opens.
3. Mark any notification attempts and click **Add**.
4. Click **Save** to save the recorded finding.

The following are alternate ways to use the Record Findings option:

1. Click **Image** in the Image Viewer and select **Critical Findings**. Click **Record Findings** at the bottom of the Add Critical Findings dialog box.
2. Click the triangle **Critical Findings** symbol on an image in the Image Viewer. Click **Record Findings** at the bottom of the Add Critical Findings dialog box.

EDIT PHYSICIAN INFO

When using the Record Findings option, if the Physician field is blank or the needed Physician is not on the list, users can search for and edit Physician info.

1. Click **Menu** in the Image Viewer and select **Image**. Select the **Critical Findings** drop-down list and select **Record Findings**. The Findings Record Editor dialog box opens.
2. Click **Add** to add a finding. The Physician Notification Editor dialog box opens.
3. In the Physician Notification Editor dialog box, click the **magnifying glass icon** next to the **Physician** drop-down list. The Physician Search dialog box opens.
4. Enter the Physician's information in the search fields and click **Search**.
5. When the specific Physician appears in the Physicians list, select it and click **OK**. The Physician's information populates in the Physician Notification Editor dialog box.

To revise physician notification information:

1. Click **Menu** in the Image Viewer and select **Image**. Select the **Critical Findings** drop-down list and select **Record Findings**. The Findings Record Editor dialog box opens.
2. Click **Add** to add a finding. The Physician Notification Editor dialog box opens.
3. In the Findings Record Editor dialog box, select the name of the physician in the **Physicians Notifications** field by highlighting it with the mouse, and then select **Edit**. The Physicians Notification Editor dialog box opens.
4. Make any needed changes and click **Add**.
5. Select **Store**. A dialog box appears warning that the critical findings and physician notification will become a permanent record. Click **Yes** to save.

To delete the physician notification information:

1. Click **Menu** in the Image Viewer and select **Image**. Select the **Critical Findings** drop-down list and select **Record Findings**. The Findings Record Editor dialog box opens.
2. Click **Add** to add a finding. The Physician Notification Editor dialog box opens.
3. In the **Findings Record Editor** dialog box, select the name of the physician in the **Physicians Notifications** field by highlighting it with the mouse. Click **Remove**. The name of the Physician is deleted.
4. Select **Store**. A dialog box appears warning that the critical findings and physician notification will become a permanent record. Click **Yes** to save.

12-17—EQUALIZE HISTOGRAM

Equalize Histogram allows users to enhance the contrast of an image by equalizing the image's histogram. This option is helpful when viewing images of the lungs. Use the Equalize Histogram option on a PA. The lung apices and angles are easily visible, the lung markings are clear, and the mediastinum is not over- or underexposed.

Note: It is important to note when a bucky takes an image. A bucky has a grid that improves the image quality. Portable chest x-ray studies only have one view; therefore, a PA and those images do not have a grid.

Click **Menu** in the Image Viewer and select **Image**. Select **Equalize Histogram**.

12-18—SET WIN/LEVEL

Set Win/Level allows users to manually enter a window/level value for a selected image.

Click **Menu** in the Image Viewer and select **Image**. Select **Set Win/Level**. The Set Window/Level dialog box opens. Manually enter the level numbers in the fields for the window/level settings.

12-19—SMOOTH

Smooth Error! Bookmark not defined. allows users to soften the edge enhancement applied to the selected image(s). The Smooth option can be applied multiple times until the desired effect is created.

Click **Menu** in the Image Viewer and select **Image**. Select **Smooth**.

Note: The Smooth and Sharpen options work in direct opposition to each other. Users can apply these two options to correct each other.

12-20—SHARPEN

Sharpen allows users to sharpen the detail or edge enhancement applied to the selected image(s). The Sharpen option can be applied multiple times until the desired effect is created.

Click **Menu** in the Image Viewer and select **Image**. Select **Sharpen**.

Note: The Smooth and Sharpen options work in direct opposition to each other. Users can apply these two options to correct each other.

12-21—LOCAL CONTRAST ENHANCEMENT

Local Contrast Enhancement allows users to enhance details in CR images, specifically the bone in the images.

Click **Menu** in the Image Viewer and select **Image**. Select **Local Contrast Enhancement**. Users can also press the default hotkey **Shift + E**.

Note: Users can also set the option to apply the Local Contrast Enhancement automatically when CR images are opened.

CHAPTER 13—MISCELLANEOUS INFORMATION

The Image Viewer's Navigation, Customize Toolbars, and Help menus contain many functions to help users as they view images. Users can navigate through images, series, and studies, customize the various **Image Viewer** toolbars, and explore help topics. The Navigation, Customize Toolbars, and Help menus are located on the Image Viewer's Menu list. If set up, these menus can also be found on the Image Viewer's Main Toolbar.

Note: The Image Viewer's toolbars can be customized to contain frequently used functions from any menu to increase efficiency. Users can also add any options or functions to Image Context Menus for quick access in the Image Viewer.

13-1—NAVIGATION

The **Navigation Menu** contains a convenient list of options to help users navigate to the Study Browser and among images in the Image Viewer. The Navigation menu is located on the Image Viewer's Menu. If set up, this menu can also be found on the Image Viewer's Main Toolbar.


Click **Menu** in the Image Viewer and select **Navigation**. Navigate between worklist items, workspaces, studies, series, and images using the Navigation menu options.

Note: Any of the Navigation Menu options can be assigned a hotkey.

- A. **Study Browser**—launches the Study Browser window.
- B. **Next Worklist Item**—closes the current workspace and opens the next worklist study. The name of the worklist being used appears in parentheses next to this menu item.
- C. **Mark Reviewed, Close, and Open Next Worklist Item**—marks the current study as Reviewed and closes it. It then opens the next worklist item from the worklist being used. The name of the worklist being used appears in parentheses next to this menu item.
- D. **Previous Workspace**—changes the view to the study located on the previous workspace.
- E. **Next Workspace**—changes the view to the study located on the next workspace.
- F. **Previous Study**—switches to the previous study within the same viewer. When multiple comparisons are open together, users see each study as a tab in the comparison viewer.
- G. **Next Study**—switches to the next study
Note: If the first or last study tab is selected, and the user tries to navigate forward or backward past that first or last study, the navigation does not cycle around the study tabs.
- H. **Previous Series**—changes the view to the previous series in the study. The default hotkey is the **up arrow**.
- I. **Next Series**—changes the view to the next series in the study. The default hotkey is the **down arrow**.
- J. **Previous Page**—changes the view to the previous page of images in the current series. The default hotkey is **Page Up**.
- K. **Next Page**—changes the view to the next page of images in the current series. The default hotkey is **Page Down**.
- L. **Previous Image**—changes the view to the previous image in the current series. The default hotkey is the **left arrow**.
- M. **Next Image**—changes the view to the next image in the current series. The default hotkey is the **right arrow**.
- N. **First Image**—changes the view to the first image in the current series. The default hotkey is **Home**.
- O. **Last Image**—changes the view to the last image in the current series. The default hotkey is **End**.

13-2— THE RADIUS MENU

The Radius Menu is an alternative to the standard right-click menu. Instead of the options appearing in the side menu, they will appear in concentric circles around the user's mouse. Users may find this helpful because many of the features will appear closer to the mouse than they do with the standard right-click menu.



To enable the Radius Menu:

1. Click Settings.
2. Click Options.
3. Make sure the General Tab is selected.
4. Select either Palette Mode or Context mode underneath the Radius Menu section of the General Tab.

WHEN PALETTE MODE IS SELECTED, A SET OF OPTIONS WILL APPEAR. THE SELECTION OF OPTIONS WILL REMAIN CONSISTENT ACROSS ALL MODALITIES AND STUDIES.

Note: When using the Radius Menu in the 3D Viewer, only Palette mode will be available.

When Context mode is selected, a smaller set of consistent options will appear, and the rest will vary depending on the specific images being viewed.


5. Click OK when the desired mode has been selected.

Users can customize the Radius by selecting the gray arrow directly below the menu.



A selection menu will appear. Here, Users can search for the feature they want to add to the Radius Menu by dragging the feature to an empty spot on the menu.

Note: If there is no empty spot to drag a feature to in the Radius Menu, users will have to right-click the feature and select the category they want to add it to from the menu that appears.



Users can also right-click any of the features and a menu containing all of the menu categories will appear. To add the feature to the menu, click on the category to house the desired feature.

To remove a feature from the menu, right-click the feature and deselect the category in the menu that appears.

RADIUS MENU IN THE 3D VIEWER

When users enable the Radius Menu it will automatically enable for the 3D Viewer as well. It is important to note that only Palette Mode will work in the 3D Viewer, not the Context Mode.

13-3—ADDITIONAL SOFTWARE PRODUCTS

There are additional software products that can be used in conjunction with the NovaPACS Diagnostic Viewer. Users can contact Novarad Support for more information.

NOVAORTHO® ORTHOPEDIC PACS

Web-based Orthopedic PACS for image and information management that includes the following key features:

- Orthopedic trauma tools, templates, and pre-op planning tools
- Fast and accurate desktop CR with cart workstation and three-year warranty and maintenance
- On-site archive and redundant archive for image storage
- Images and pre-op plans archived in NovaPACS for fast and easy retrieval
- Admin Console
- Installation and on-site training
- 24/7 technical support service

NOVASCAN® DIGITIZE TO DICOM

This software tool allows users to create a DICOM file from images produced by non-DICOM equipment. NovaScan software accepts image data from intermediary imaging equipment like film digitizers and video capture cards. NovaScan generates files that contain image pixel data and patient/study information that users send to a DICOM server for review by DICOM viewers.

Users can run NovaScan from the NovaPACS Study Browser. This feature was implemented to allow users to scan documents directly from NovaPACS without having to open any separate applications.

A Scan Documents link has been added within NovaPACS to launch NovaScan. The Scan Documents link appears in the Study menu, the right-click menu, in the toolbar items under the Documents section, and in the keyboard shortcuts section.

Selecting any of the Scan Documents options launches NovaScan in the same manner as NovaPACS is running.

- If NovaPACS is a managed client, the system attempts to run an installed NovaScan.
- If NovaPACS is click-once, run the deploy for NovaScan. If accessing NovaScan through the web viewer, the NovaScan click-once version runs.

If NovaScan is already running when a new session is launched, the current one closes and the new session launches. This prevents multiple sessions from opening simultaneously.

NOVAPRO® PACS FOR OFFSITE RADIOLOGIST OR RAD GROUP

NovaPro provides any Radiologist or Radiology group the ability to receive DICOM images from unlimited, unaffiliated imaging producing practices with no fear of intermingling of images and/or information. Distributes the resulting reports automatically via email and/or fax or through direct HL7 messaging into the customer's own HIS, RIS, and or EMR/EHR. The system automatically scrubs images from the archive when reaching the pre-determined image storage limit. NovaPro provides this functionality without needing to store or maintain the image files. The provider can configure the system to be a complete PACS solution allowing them to store and maintain images and reports.

NOVAPRO® SINGLE USER

This product comes fully functional to tie into a group's workflow with no added expense.

TRAUMACAD® ORTHO PRE-OP PLANNING PLUG-IN

TraumaCad allows doctors to use digital images on-screen. Physicians can perform measurements, fix prostheses, simulate osteotomies, and visualize fracture reductions. Procedure-specific measurement tools and semi-automatic wizards guide the measurement process step-by-step. TraumaCad allows doctors to replicate every action currently performed by tracing on x-ray films, only more conveniently and much faster. TraumaCad has the largest digital template library with over 90,000 template images available. There are ten planning modules including Auto-Hip, Knee, and Upper Limb. Integrated orthopedic hardware lists facilitate pre-op planning.


NOVA RIS (RADIOLOGY INFORMATION SYSTEM)

Streamlines workflow, provides timestamps, tracks, and reports for every aspect of the process. Completely web-based and runs on existing workstation hardware. Provides the radiology department with the following features:

- Drag-and-drop, web-based scheduling, workflow tracking, order entry, interfaces, reports, dictation, and voice recognition
- Referring physician exam requests and orders entered from their offices go directly into Nova RIS
- Alert pop-ups for billing and insurance preauthorization requirements, contraindications by patient, prep instructions, and exam orders
- Transcription module includes foot pedal control, playback, dictation, and editing within Nova RIS
- Patient reminder letters, e-sign, auto-refresh of technology screen for easy patient tracking, and many other user-friendly features
- ONC-ATCH ambulatory certification and integration with certified EMR/HER modules to provide complete Meaningful Use compliance.
- Integration with the new Nova RIS Billing module.

NOVAVAULT™ DATA PROTECTION

Archive, backup, and storage solution service designed to provide on-site and off-site storage and off-line disaster recovery of data. There are three versions of NovaVault. These products help sites with long-term access to critical



clinical information, managed services, and disaster recovery services. This product may also be installed at sites not currently utilizing NovaPACS.


NOVACARDIO® IMAGING PLATFORM

NovaCardio is accessible anywhere and offers an efficient workflow and powerful features. NovaCardio contains four modules:

- 1) The ECG Viewer, which launches when users select a study with the ECG modality that produces EKG strips
- 2) Motion Cardiac Imaging Viewer to view ultrasound studies containing stress echo images
- 3) NovaPACS with a NovaCardio license to view cardiac catheterization, vascular ultrasounds, and nuclear medicine imaging
- 4) EncaptureMD Structured Reports.



Additional features of the NovaCardio Imaging Platform include:



- **EncaptureMD SR**—With full support for both invasive and non-invasive cardiology procedures, the NovaCardio with EncaptureMD Echocardiography Structured Reporting Application is the fastest way to create all types of accredited cardiology-related reports. EncaptureMD simplifies the documentation process and eliminates the need for transcription and dictation services.

- **ECG Management System**—The new NovaCardio ECG Viewer has been added to NovaPACS and allows users to read ECG studies and create reports. Cardiologists launch the ECG Viewer

from within the NovaPACS Study Browser, which allows users to view waveforms and data, compare prior studies, and save, create, edit, and distribute reports.

- **Syntermed's Emory Cardiac Toolbox™**—NovaPACS integrates with Syntermed™ nuclear studies viewer and Emory Cardiac Toolbox. This allows Cardiologists to streamline their workflow by launching nuclear medicine studies directly from the NovaPACS Study Browser.

NOVADOSE® IONIZING RADIATION DOSE INFORMATION SYSTEM

Novarad provides integrated Radiation Dose Monitoring that runs on existing infrastructure and is modality vendor neutral. NovaDose calculates a patient's ionizing dose history for CT, XA, MG, and DR modalities and provides alerts based on a customer-defined dose threshold when a test is ordered. When used with both NovaPACS and Nova RIS, NovaDose automatically adds a complete study dose summary for all series to the RIS Report. In addition to patient monitoring, department directors, CIO's, and physicists use NovaDose to quickly identify dosing performance issues by tech, procedure, and modality. Simple queries produce comparison charts instantly to provide data to assist with decision making. For example,

- Tech A consistently doses higher than other techs when performing Head CT.
- CT Machine #4 is set higher than other CT machines.

NOVA RIS BILLING MODULE

Novarad's Billing Client product fully integrates with Nova RIS and offers radiology practices a comprehensive billing service that emphasizes verify, code, process, and collect workflows. With features like automatic coding, complete transparent/anytime access to billing status in detail, and verification of patient benefits at check-in, Nova RIS Billing improves efficiency and cashflow. The billing module allows for the following:

- Manage patient accounts
- View unbilled procedures
- Seamlessly bill insurance and patients
- Add, edit, or delete claims
- Explanation of Benefits (EOBs)
- Batch processing

- *Electronic Data Interchange (EDI)*

NOVA RIS EHR MODULE

Novarad's ambulatory EHR solution allows Nova RIS users to maintain an efficient and effective workflow within the RIS while providing full compliance for Meaningful Use stages. This product allows users to record, modify, and retrieve patient medical, billing, insurance, and demographic data. A compliance dashboard allows physicians to quickly review their compliance progress visually.

NOVAMG® PRO MAMMOGRAPHY

NovaMG PRO is a multi-modality workstation with an integrated mammography viewer that allows users to read a wide variety of mammography studies, as well as other studies, all in one place.

CRYPTOCHART

CryptoChart is a mobile app that provides the following features:

- *Securely share medical images without passwords or logins*
- *Full DICOM viewing system*
- *Mobile compatible with PACS integration*

CRYPTOSHARE

CryptoShare is a mobile app that provides the following features:

- *No username or password needed*
- *Integrated with MobileRAD and WhatsApp: Fully mobile environment*
- *Fast, trackable image sharing that is compatible with PACS*

MOBILERAD

MobileRAD is a mobile app that provides the following features:

- *Efficient and secure*
- *Easy access to schedules and patient information*
- *Secure communication and documentation*

ALERTVIEW

AlertView is a mobile app that provides the following features:

- *Automatically send reports/images*
- *Easily connect with referring physicians*
- *Customizable settings*

SNAPVIEW

SnapView is a mobile app that provides the following features:

- *Securely capture photos/videos/documents*
- *Organized by department*
- *Rapid and secure documentation*
- *No Personal Health Information (PHI) left behind*

INDEX

Accessing Menu Options.....	99	Choose Series Box Docking	112
Accessing the Diagnostic Viewer	2	Choosing Column Headings.....	31
Add Critical Findings.....	148	Cine	83, 109, 128
Add Split.....	112	Circle	141
Add Study Split when Opening a Second Study for a Patient.....	62	Clear Worklist.....	38
Add to Worklist.....	36	Close.....	128
Adding a Data Source	76	Close All Workspaces.....	40
Adding Labels.....	72	Close All Workspaces and Show Study Browser	106
Adding Lexicon Values.....	74	Close Application.....	95
Adding Toolbar Items	89, 97	Close Current Study.....	95
Additional Software Products.....	153	Close Study.....	105
Aggregate	115	Close Workspace.....	40, 105
Always Match Layout.....	62	Cobb Angle.....	143
Always Sort Filtered Studies by Priority	62	Collaboration	51, 123
Annotate.....	139	Color	140
Annotation Toolbar.....	95	Colorize	132
Apply to Series.....	138	Compression	45
Applying a Saved Filter	22	Configure Default HPs.....	78
Arrow.....	141	Configure Monitors.....	84
Assign Password	46, 122	Contact Support	92
Auto Sync.....	114	Coordinates	141
Auto W/L Tab	68	Create DICOM Label.....	73
Auto-checkout Studies when they are Opened	61	Create New View	133
Auto-Local Contrast Enhance CRs	62	Creating Advanced Worklists.....	30
Automatically login with this account.....	4	Creating Labels	73
Auto-review Studies when they are Closed	61	Critical Findings.....	148
Auto-sharpen XAs.....	62	Cross Localize.....	114
Auto-show Study Browser.....	60	Current toolbar items.....	33
Auto-subtract	60	Current User Information.....	91
Auto-sync Hanging Protocol Image Layout Options	62	Customize Toolbar	32, 89
Calcium Scoring	131	Customizing and Saving Filters.....	21
Calibrate Measurement.....	137	Customizing Image Viewer Toolbars.....	97
Changing Dictation Font Display Size.....	121	Customizing Study Browser Lists.....	31
Changing the Order of Columns	31, 32	Customizing the Filters Group	21, 28
Checkout.....	51, 123	Customizing the Study Browser's Custom Toolbar	32
Choose Image Layout.....	131	Customizing Toolbar Options	98
Choose Series.....	112	Data Source Filter	23

Data Sources Tab	75	Hide All	144
Deleting a Saved Filter	22	History	1, 7, 39, 56, 61, 144
Deleting Advanced Worklists	30	Hotkeys	64
Details	138	HTML5 Viewer	3
Diagnostic Viewer	4	Image Context Menu	66
Dictate	119	Image Context Menu Tab	66
Display Older Related Studies when Opening a Study	61	Image Fusion	109, 130
Display Separate Window when Creating a Key Image Series	62	Image Menu	136
Display Spine Label Cross References	63	Image Viewer Icon Buttons	99
Docking Tabs	19, 58	Images Tab	17
Documents	41, 118	Include Diagnostic Reports	44
Drag and Drop Lists/Tabs	18	Instructions	91
Draw Arrow Annotations Head to Tail	64	Introduction to the Worklist	36
Edit Physician Info	149	Invert	137
Editing Advanced Worklists	30	Key Image Series	147
Editing Data Sources	76	Keyboard Shortcuts Tab	64
Editing the Series Filter	110	Keyboard Shortcuts Tab (Hotkeys)	64
Ejection Fraction	134	Labels Tab	71
Ellipse	141	Launch Mammo Viewer	121
Equalize Histogram	150	Launch Mammo Viewer	46
Erase	144	Launch Orthopedic Tools	47
Erase All	144	Launching the Web Viewer	4
Exam Report	41, 118	Lexicon Tab	74
Exit	127	Local Contrast Enhancement	151
Exit the Diagnostic Viewer	35	Locking Toolbars	98
Export Preferences	86	Log into Datasources	5
Export to File	147	Logging in Automatically	4
Export User Preferences	87	Logging in to a Data Source	4, 6
Fade Away after This Many Seconds	63	Login	3
File Menu	35	Looping	83
Filter	23, 29	Magnification	147
Filtered Studies	7, 51, 53, 62, 123, 124	Magnify Glass Size	64
Filters	21, 28, 36	Manage Hanging Protocols	55
Flip/Rotate	139	Manage Hanging Protocols	55
Font	140	Manage Report Workflow	42
General Tab	60	Manage Reports Workflow	118
Hanging Protocol	21	Managed Viewer	3, 43
Hanging Protocol Layouts	19, 78	Measure Angle	143
Help Menu	91	Measure Area	143
		Measure Length	142
		Measure Path	142

Minimize	127	Percent Stenosis	134
Minimum Mouse Drag Distance	63	Play	144
Most recent preference backups	87	Plug-ins Tab	75
Mouse Control Options	62	Point	141
Mouse Scroll Sensitivity	63	Point Cross Localize	115
Mouse Settings Tab	65	Polygon	140
Mouse Wheel Sensitivity	63	Preferences	85
Move Images	147	Print	115, 116
Moving Toolbars	98	Prompt Before Clearing the Worklist	63
MPR	4, 47, 83, 109, 129	Prompt Before Closing the Application	63
MPR View	129	Prompt before Marking a Study Reviewed	61
Navigate through Annotated Images	146	Prompt Before Marking a Study Reviewed if it has not Been Completely Reviewed	63
Navigation Menu	152	Prompt Before Saving Annotations	63
Never Match Layout	62	Prompt for Related Studies when Needed	61
Next Worklist Item	38	Prompt to Allow Cross Localization of Different Frames of Reference	63
Notes	122	Quick Reference Toolbar	95
Nova RIS	154	Record Findings	148
Nova RIS Billing Module	156	Rectangle	140
Nova RIS EHR Module	157	Redo	144
NovaCardio	155	Refresh Lists	57
NovaDose	156	Refresh Study Browser when Idle	63
NovaMG PRO	157	Reload	137
NovaOrtho	153	Remember Me	4
NovaPACS Diagnostic Viewer System	2	Remove Images	133
NovaPro	154	Remove Worklist Item	38
NovaPro Single User	154	Removing Toolbar Items	97
NovaScan	153	Re-ordering the Menu	90, 97
NovaVault	155	Report tab	18
Only Match Layout if Series are Coplanar	62	Research	126
Open	39	Reset Cursor	140
Open a New Workspace for Each Patient	61	Reset Cursor when Closing a Study	62
Open from Directory	35	Reset Filters	57
Open from File	35	Reset Filters when Opening the Study Browser	61
Open Worklist Item	37	Reset Status to New	53, 124
Options	60	Resizing Column Headings	31
Patient CD/DVD	43	Resizing List Heights	31, 34
Patient History	11, 103	Restore from Backup	87
Patient History list	56	Retrieve Images	43
Peer Review	124	Reviewed	52, 124
Peer Reviews	53		
Pen	140		

RIS Information.....	42, 119	Subtract Clip Frames.....	137
ROI (Region of Interest)	140	SUV Measurements with Image Fusion	130
Route	40, 117	Sweeping	83
Routing Priority	41	Switch Between Workspaces	95
Save View Settings	126	Teaching File Filters.....	27
Scale	144	Teaching Files.....	27, 48, 50, 64, 125
Scan Complete.....	46	Text.....	141
Select All.....	39, 40	The Advanced Worklist Tab	29
Select All Images.....	128	The Assigned Worklist Tab	28
Select All Series.....	105	The Data Source Filter	29
Series list.....	13	The Modality Filter	24
Series Menu.....	128	Tiled	82, 109, 128
Set Sync	115	Time Frame Filter.....	23, 29
Set Win/Level	150	Undo.....	144
Settings Menu	60	Undo All Selections.....	147
Sharpen	150	Undo Checkout.....	52, 123
Show All Series	126	Upload to GlobalRad.....	126
Show Only the Most Recent Study for Each Patient in Filtered Studies.....	62	Use Multiple Windows to Display RIS Reports.....	62
Show Other Users' Annotations.....	61	User Manual	92
Show Patient History when Opening a Study...	61	Using the Image Viewer	94
Show RIS Information when Opening a Study .	63	Using the Patient ID, Last Name, and Accession # Filters	26
Show Scale.....	144	Using the Study Browser Custom Toolbar.....	32
Show Study Notes when Opening a Study.....	61	Using the Study Browser File Menu	35
Single Study Access.....	3	Using the Study Browser Menus/Custom Toolbar	32
Single-Study Access	46, 122	Utilization Review.....	54, 125
Smart Image Layout	60	Verified	52, 123
Smart Worklist	10	View Menu	55
Smooth.....	150	View Study on Other Workstation.....	47, 121
Sort Images.....	132	Views	128
Sorting Data in a Column	31	W/L Presets.....	77
Spine Labeling.....	141	Web Viewer.....	3, 4, 46
Split Horizontal.....	131	What's New	92
Split Vertical.....	131	Where to Find Additional Study Browser Menu Information.....	32
Stop.....	83	Window/Level Toolbar	95
Study Browser Color Themes.....	34	Worklist.....	28, 36, 55
Study Browser icon	95	Worklist Menu	36
Study Menu	39, 40, 41, 42, 43, 46, 103	Zoom	136
Study of Interest	27, 50		
Submit Feedback	92		
Subtract.....	133		

DOCUMENT REVISION HISTORY

<i>Version</i>	<i>Revision</i>	<i>Date</i>	<i>Description of Revision</i>
9.1.5	1.0	11/2025	<ul style="list-style-type: none">• Improvements were made to the authentication process.• Enhancements were made to the Search filter by including Include/Exclude controls.
9.1.6	1.0	05/2026	<ul style="list-style-type: none">• Enhancements were made to Fluency Direct with the addition of a new section on the integration